



VERSION 2
SEP 07, 2023

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.rm7vz3pb2gx1/v2

Protocol Citation: Elena Essel, Matthias Meyer, Petra Korlevic 2023. A method for the temperature-controlled extraction of DNA from ancient bones . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vz3pb2gx1/v2> Version created by [Elena Essel](#)

MANUSCRIPT CITATION:
A method for the temperature-controlled extraction of DNA from ancient bones
Elena Essel, Petra Korlević, and Matthias Meyer
BioTechniques 2021 71:1, 382-386
<https://doi.org/10.2144/btn-2021-0025>

🌐 A method for the temperature-controlled extraction of DNA from ancient bones V.2

🔗 Version 1 is forked from [A method for the temperature-controlled extraction of DNA from ancient bones](#)

Elena

Essel¹, Matthias Meyer¹, Petra Korlevic^{1,2}

¹Department of Evolutionary Genetics, Max Planck Institute for Evolutionary Anthropology, Deutscher Platz 6, D-04103 Leipzig, Germany;

²EMBL-EBI, Wellcome Genome Campus, Hinxton, Cambridgeshire, CB10 1SD, UK

MPI EVA Ancient DNA Core Unit



Elena Essel

Department of Evolutionary Genetics, Max Planck Institute fo...

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Sep 07, 2023

Last Modified: Sep 07, 2023

PROTOCOL integer ID: 87493

Keywords: Ancient DNA, sequential DNA extraction, contamination removal, endogenous DNA, archaeological material

ABSTRACT

We here provide a protocol for the decontamination of ancient bones and teeth that is based on a temperature-controlled, sequential release of DNA. DNA can be extracted from all fractions generated with this method and the fraction with the highest proportion of endogenous DNA identified for further analysis. The protocol proceeds through repeated incubation of the sample powder in phosphate buffer at 37, 60 and 90 °C, followed by the complete lysis of the residual sample powder. As DNA is denatured at high temperature, subsequent DNA extraction and library preparation has to be performed using methods optimized for single-stranded DNA.

MATERIALS

Reagents

✕ Sodium phosphate, 0.5M buffer soln., pH 7.0 Thermo Scientific Catalog #AAJ63791AP

✕ Water for HPLC Sigma Aldrich Catalog #270733

✕ EDTA solution pH 8.0 (0.5 M) for molecular biology AppliChem Catalog #A4892,1000

✕ Tris buffer pH 8.0 (1 M) for molecular biology AppliChem Catalog #A4577,0500

✕ Proteinase K 100 mg Sigma-aldrich Catalog #3115879001

✕ TWEEN® 20 Sigma Aldrich Catalog #T2700-100ML

Consumables and equipment

✕ DNA LoBind Tubes 2.0 mL Eppendorf Catalog #0030108078

✕ DNA LoBind Tubes 2.0 mL Eppendorf Catalog #0030108078

✕ Ceramic beads 2.8 mm VWR International Catalog #432-0292

✕ 50 ml CELLSTAR® Polypropylene Tube 30/115 MM Conical Bottom Blue screw cap sterile skirt greiner bio-one Catalog #210261

✕ Parafilm M 10 cm wid neoLab Catalog #3-1012

Equipment		
Thermomixer		NAME
HLC		BRAND
52 82 00133		SKU

Equipment		
Incubator		NAME
Memmert		BRAND
Incubator IN55		SKU


Equipment		
Tube rotator		NAME
VWR		BRAND
444-0500		SKU

Equipment		
UV cross-linker		NAME
Vilber		BRAND
Bio-Link BLX 254		SKU


Equipment		
Vortex mixer		NAME
Scientific Industries		BRAND
SI-0236		SKU

Equipment		
Centrifuge		NAME
Bench centrifuge		TYPE
Eppendorf		BRAND
5424		SKU


PROTOCOL MATERIALS

 TWEEN® 20 Merck MilliporeSigma (Sigma-Aldrich) Catalog #T2700-100ML

In 2 steps

 EDTA solution pH 8.0 (0.5 M) for molecular biology AppliChem Catalog #A4892,1000


Step 4

 Proteinase K 100 mg Merck MilliporeSigma (Sigma-Aldrich) Catalog #3115879001


Step 4

 Ceramic beads 2.8 mm VWR International Catalog #432-0292


Step 6

 Sodium phosphate, 0.5M buffer soln., pH 7.0 Thermo Scientific Catalog #AAJ63791AP

Materials, Step 2

 Water for HPLC Merck MilliporeSigma (Sigma-Aldrich) Catalog #270733

In 2 steps

 Tris buffer pH 8.0 (1 M) for molecular biology AppliChem Catalog #A4577,0500

Step 3

Buffer preparation

1


Note


All buffers are irradiated with UV-C light at a dose of 7 kJ/cm² using a cross-linker.

2

Sodium-phosphate buffer (0.5 M sodium phosphate, pH 7.0, 0.1 % Tween 20) is prepared by combining the following reagents:



 49.5 mL


 Sodium phosphate, 0.5M buffer soln., pH 7.0 Thermo Scientific Catalog #AAJ63791AP


 50 µL



 TWEEN® 20 Sigma Aldrich Catalog #T2700-100ML

- 3 Tris-Tween wash buffer (10 mM Tris-HCl, pH 8.0, 0.1% Tween-20) is prepared by combining the following reagents:



 49.5 mL  Water for HPLC Sigma
Aldrich Catalog #270733


 0.5 mL


 Tris buffer pH 8.0 (1 M) for molecular biology Sigma
Aldrich Catalog #A4577,0500



 50 µL  TWEEN® 20 Sigma
Aldrich Catalog #T2700-100ML


- 4 Lysis buffer (0.45 M EDTA, pH 8.0, 0.05% Tween-20 and 0.25 mg/ml proteinase K) is prepared by combining the following reagents:


 3.725 mL  Water for HPLC Sigma
Aldrich Catalog #270733

 45 mL

 EDTA solution pH 8.0 (0.5 M) for molecular biology Sigma
Aldrich Catalog #A4892,1000

 25 µL  TWEEN® 20 Sigma
Aldrich Catalog #T2700-100ML


 1.25 mL 10 mg/ml proteinase K solution in water (prepared from

 Proteinase K 100 mg Sigma
Aldrich Catalog #3115879001)

Note



Proteinase K is added after UV irradiation

Sample preparation

- 5 In an ancient DNA cleanroom, remove approximately  50 mg of sample powder from each specimen using a sterile dentist drill and transfer the powder to a 2.0 ml DNA LoBind tube.
- 6 To facilitate resuspension of the bone powder during the subsequent incubation and wash steps,




add 3-4  Ceramic beads 2.8 mm Sigma Aldrich Catalog #432-0292 to the sample material.

Temperature-controlled phosphate treatment

- 7 Add  0.5 mL sodium phosphate buffer to the sample powder, completely resuspend the powder by thorough vortexing, and incubate the tube in a thermo block adjusted to the desired temperature  900 rpm, 00:15:00

Note

Temperature-controlled phosphate treatment steps

-  37 °C 2 times
-  60 °C 2 times
-  90 °C 2 times

Note

At least one negative control (tube without sample material) should be included in each experiment and carried through all subsequent steps).

- 8 Transfer tubes to a tabletop centrifuge and spin for 2 min at maximum speed (e.g., 16,400g/13,200 rpm).
- 9 Transfer supernatant to a 1.5 mL LoBind tube and store at -20 °C until the day of DNA extraction.


Note

Beads facilitate the resuspension of the sample powder after centrifugation steps, but make it harder to remove supernatant.
Pipette slowly and carefully.



- 10 Repeat steps 7-9 once at each temperature (for a total of 2 wash steps).

Note

For the 90 °C incubation, make sure the liquid in the tube reaches 90 °C by the end of the 15 min incubation time. If necessary, set the thermo block to a higher temperature.

- 11 The temperature-controlled phosphate treatment is followed by a room-temperature wash step with  1 mL Tris-Tween buffer at the end of the last temperature cycle. Completely resuspend the powder by thorough vortexing.
- 12 Transfer tubes to a tabletop centrifuge and spin for 2 min at maximum speed (e.g., 16,400g/13,200 rpm)
- 13 Transfer supernatant to a 1.5 mL LoBind tube and store at -20 °C until the day of DNA extraction.

Final digestion of sample material

- 14 Add  1 mL of lysis buffer to the sample powder, completely resuspended the powder by vortexing, and incubate overnight (8 – 16 h) with rotation at  37 °C


Note

Wrap the tube with parafilm to prevent leaking.

- 15 Transfer tubes to a tabletop centrifuge and spin for 2 min at maximum speed (commonly at 16,400 g/13,200 rpm).

- 16 Transfer supernatant to a 1.5 mL LoBind tube and proceed to DNA extraction or store the tube at -20 °C until the day of DNA extraction.

DNA purification of phosphate fractions and final lysate

- 17 Thaw the sodium phosphate fractions (and lysates if necessary) at  37 °C in a thermo block with gentle shaking.

Note

Make sure the liquid is fully thawed and any crystals have completely dissolved.

Note

If desired, DNA extraction can also be performed on the Tris-Tween buffer, but DNA yields are expected to be extremely low.

- 18 For the sodium phosphate fractions, purify 100 µl of the supernatant, and for the final lysate, purify 500 µl using binding buffer 'G' of the DNA extraction method described in Glocke and Meyer (2017). Final volume of all DNA extracts is 50 µl.

CITATION

Glocke I, Meyer M (2017). Extending the spectrum of DNA sequences retrieved from ancient bones and teeth.. Genome research.

LINK

<https://doi.org/10.1101/gr.219675.116>

Library preparation, sequencing, and data processing

- 19 Prepare DNA libraries using 20% of the DNA extract as input, following the protocol for library preparation, quantification and indexing by Gansauge et al. (2020).

CITATION

Gansauge MT, Aximu-Petri A, Nagel S, Meyer M (2020). Manual and automated preparation of single-stranded DNA libraries for the sequencing of DNA from ancient biological remains and other sources of highly degraded DNA.. Nature protocols.

LINK

<https://doi.org/10.1038/s41596-020-0338-0>

- 20** Perform shallow shotgun sequencing on Illumina's MiSeq or HiSeq2500 platforms (or other Illumina platforms) using a paired-end double-index configuration (2x 76 + 2x 7 cycles).

CITATION

Kircher M, Sawyer S, Meyer M (2012). Double indexing overcomes inaccuracies in multiplex sequencing on the Illumina platform.. Nucleic acids research.

LINK

<https://doi.org/10.1093/nar/gkr771>

Sequence analysis

- 21** Trim adapters and merge overlapping paired-end reads into single-molecule sequences using leeHom.

CITATION

Renaud G, Stenzel U, Kelso J (2014). leeHom: adaptor trimming and merging for Illumina sequencing reads.. Nucleic acids research.

LINK

<https://doi.org/10.1093/nar/gku699>

- 22** Use the Burrows-Wheeler Aligner (BWA, <https://github.com/mpieva/network-aware-bwa>) to align merged sequences to a suitable reference genome (e.g. turTru1.75, bosTauUMD3.1, loxAfr4) using ancient parameters (" -n 0.01 -o 2 -l 16500") allowing more mismatches and

indels.

CITATION

Li H, Durbin R (2010). Fast and accurate long-read alignment with Burrows-Wheeler transform.. Bioinformatics (Oxford, England).

LINK

<https://doi.org/10.1093/bioinformatics/btp698>

CITATION

Meyer M, Kircher M, Gansauge MT, Li H, Racimo F, Mallick S, Schraiber JG, Jay F, Prüfer K, de Filippo C, Sudmant PH, Alkan C, Fu Q, Do R, Rohland N, Tandon A, Siebauer M, Green RE, Bryc K, Briggs AW, Stenzel U, Dabney J, Shendure J, Kitzman J, Hammer MF, Shunkov MV, Derevianko AP, Patterson N, Andrés AM, Eichler EE, Slatkin M, Reich D, Kelso J, Pääbo S (2012). A high-coverage genome sequence from an archaic Denisovan individual.. Science (New York, N.Y.).

LINK

<https://doi.org/10.1126/science.1224344>

- 23 Restrict further analyses to sequences of length 35 bp and above to avoid spurious alignments of short sequences with random similarity to the reference genome.
- 24 Merge sequences with the same start- and end-coordinate into one consensus sequence using bam-rmdup (<https://github.com/mpieva/biohazard-tools>).
- 25 Generate summary statistics using samtools and choose the library with the highest proportion of endogenous DNA for further sequencing. Prepare additional libraries from remaining DNA extract if necessary.

CITATION

Li H, Handsaker B, Wysoker A, Fennell T, Ruan J, Homer N, Marth G, Abecasis G, Durbin R, 1000 Genome Project Data Processing Subgroup. (2009). The Sequence Alignment/Map format and SAMtools.. Bioinformatics (Oxford, England).

LINK

<https://doi.org/10.1093/bioinformatics/btp352>