



Apr 07, 2022

# Tissue Sample Preparation for LC-MS Analysis

Joanne Chan<sup>1</sup>, Ruiqi Jian<sup>1</sup>, Lihua Jiang<sup>1</sup><sup>1</sup>Stanford University

1

[dx.doi.org/10.17504/protocols.io.bgf5jqtq6](https://dx.doi.org/10.17504/protocols.io.bgf5jqtq6)

Human BioMolecular Atlas Program (HuBMAP) Method Development Community  
Tech. support email: [Jeff.spraggins@vanderbilt.edu](mailto:Jeff.spraggins@vanderbilt.edu)

Joanne Chan

This protocol describes the procedure for protein extraction, enzymatic digestion, sample cleanup and dTMT labeling for LC-MS analysis. Proper sample preparation and clean-up are extremely important to ensuring quality LC-MS data and reproducibility. Tissues tend to vary more than cells and extra care should be taken as MS is sensitive to any sample variations. In addition, proper clean-up is a preventative measure to less damage to the instruments and superior data. Here, we described each step in detail on how to process tissue samples prior to analysis.

DOI

[dx.doi.org/10.17504/protocols.io.bgf5jqtq6](https://dx.doi.org/10.17504/protocols.io.bgf5jqtq6)

Joanne Chan, Ruiqi Jian, Lihua Jiang 2022. Tissue Sample Preparation for LC-MS Analysis.

**protocols.io**<https://dx.doi.org/10.17504/protocols.io.bgf5jqtq6>

HuBMAP, BIOMIC, MSRC, Stanford, Proteomics, Peptides, Protein Precipitation, Desalting, Tissue Lysis, Protein Digestion, TMT labeling

protocol ,

May 15, 2020

Apr 07, 2022

37085

Please make sure at any point in the experiment, do not leave the samples unattended. If not actively performing a step, put your samples on ice.

## MATERIALS

[Lysing matrix D tubes](#) **MP**

**Biomedicals Catalog #6913100**

[2-Chloroacetamide \(CAA\)](#) **Sigma**

**Aldrich Catalog #22790**

[Trizma® base](#) **Sigma**

**Aldrich Catalog #93350**

[Acetone](#) **Sigma**

**Aldrich Catalog #179124**

[Guanidine-HCl](#) **Thermo**

**Fisher Catalog #24110**

[Pierce™ Trypsin Protease, MS Grade](#) **Thermo**

**Fisher Catalog #90059**

[TMT10plex™ Isobaric Label Reagent Set, 1 x 0.8 mg](#) **Thermo**

**Fisher Catalog #90110**

[Tris\(2-carboxyethyl\)phosphine hydrochloride \(TCEP\)](#) **Sigma**

**Aldrich Catalog #C4706**

[Lysyl Endopeptidase Mass Spectrometry Grade \(Lys-C\)](#) **Fujifilm Wako Pure**

**Chemical Catalog #12505061**

[Oasis HLB 1 cc Vac Cartridge 10 mg Sorbent per Cartridge 30 µm](#)

**100/pk Waters Catalog #186000383**

All the waste needs to go to biohazard.

Ensure you have enough tissue before preparation. We recommend using at least 30 mg of tissue to start with.

## Protein Lysate Extraction

4h

- 1 Place the tissue chunk onto a cell culture plate on ice and mince with disposable scalpels to the best of your abilities.
- 2 Add **500 µL** lysis buffer (**6 Molarity (M)** GdmCl, **10 Milimolar (mM)** TCEP, **40 Milimolar (mM)** CAA, **100 Milimolar (mM)** Tris **pH 8.5**) to the minced tissue sample and transfer it to a Matrix D tube with ceramic beads.
- 3 Beat the tissue using FastPrep-24™ at 4.5 M/S, **00:00:40** at **4 °C** 2 times .
- 4 Centrifuge the sample for **00:05:00** at **7000 x g** , and transfer the sample to a new tube, and sonicating at 40% and 60% for **00:00:20** each, respectively on ice.

- 5 Heat for **00:05:00** at **95 °C** and vortex every **00:01:00**.
- 6 Centrifuge for **00:05:00**, **00:10:00**, and **00:10:00** minutes respectively at **12000 x g**, and collect the supernatant at each step.

#### Enzymatic Digestion 16h

- 7 Take **5 µL** of the sample out and dilute to measure its protein concentration using BCA kit.
- 8 Take **100 µg - 300 µg** of the sample to perform a two-step enzymatic digestion using **25 Milimolar (mM) Tris pH 8.5** as the dilution buffer.
- 9 Endoproteinase Lys-C digestion: Add 3x volume of dilution buffer to the original lysate. Add Lys-C to lysate at 1:100 (w:w) and incubate at **37 °C** for **02:00:00**.
- 10 Trypsin digestion: Add an additional 6x volume of dilution buffer to the original lysate volume for an overall <sup>14h</sup>9x dilution of the original sample volume. Add Trypsin to lysate at 1:50 (w:w) and incubate at **37 °C** **Overnight**.
- 11 Digested samples can be stored at **-80 °C**.

#### Desalting the Peptide 8h

- 12 Desalt the digest with Waters Oasis® HLB Cartridge.
  - 12.1 Condition the column with **1 mL** each of **100 % (v/v)** acetonitrile, **80 % (v/v)** acetonitrile with **0.1 % (v/v)** acetic acid, and **0.1 % (v/v)** acetic acid sequentially.
  - 12.2 Adjust the sample to **1 % (v/v)** TFA, and load the sample onto the desalting column.
  - 12.3 Wash the sample with **200 µL** of **0.1 % (v/v)** acetic acid 3 times.

12.4 Elute peptides off the column with 150uL of **40 % (v/v)** acetonitrile with **0.1 % (v/v)** acetic acid, followed by 150uL of **80 % (v/v)** acetonitrile with **0.1 % (v/v)** acetic acid.

13 Put the sample in the SpeedVac until dry.

#### TMT Labeling

4h

14 The TMT10plex™ Isobaric Label Reagents were prepared as suggested by the manufacture. 4h  
[https://www.thermofisher.com/document-connect/document-connect.html?url=https%3A%2F%2Fassets.thermofisher.com%2FTFS-Assets%2FSLSG%2Fmanuals%2FMAN0016969\\_2162457\\_TMT10plex\\_UG.pdf&title=VXNlciBHdWlkZTogVE1UMTBwbGV4IE1hc3MgVGFnIExhYmVsaW5nIEtpdHMcGmYVW5kIFJlYWdlbnRz](https://www.thermofisher.com/document-connect/document-connect.html?url=https%3A%2F%2Fassets.thermofisher.com%2FTFS-Assets%2FSLSG%2Fmanuals%2FMAN0016969_2162457_TMT10plex_UG.pdf&title=VXNlciBHdWlkZTogVE1UMTBwbGV4IE1hc3MgVGFnIExhYmVsaW5nIEtpdHMcGmYVW5kIFJlYWdlbnRz)

14.1 Resuspend the sample in **100 Milimolar (mM)** TEAB **pH 8.5** and then add assigned TMT tag to the sample.

14.2 Incubate the sample with its associated TMT tag for **01:00:00** at room temperature.

14.3 Quench the reaction with **5 % (v/v)** hydroxylamine.

14.4 Mix each set of 10 samples together in equal amounts.

15 Put the sample in the SpeedVac to dry.

16 Resuspend the sample with **200 Milimolar (mM)** ammonium formate and perform LC-MS analysis.