





Version 2 ▼

Jul 27, 2022

## Demuxlet Cell Preparation Protocol V.2

Woong-Yang Park<sup>1</sup>, Jay Shin<sup>2</sup>, Shyam Prabhakar<sup>3</sup> <sup>1</sup>SMC; <sup>2</sup>RIKEN; <sup>3</sup>GIS



dx.doi.org/10.17504/protocols.io.n2bvjyb55vk5/v2

### Human Cell Atlas Method Development Community

Shvetha Sankaran

**ABSTRACT** 

This protocol details Demuxlet cell preparation.

**ATTACHMENTS** 

Demuxlet\_Cell\_Preparatio n\_Protocol\_21April2022.do

DOI

dx.doi.org/10.17504/protocols.io.n2bvjyb55vk5/v2

PROTOCOL CITATION

Woong-Yang Park, Jay Shin, Shyam Prabhakar 2022. Demuxlet Cell Preparation Protocol. **protocols.io** 

https://dx.doi.org/10.17504/protocols.io.n2bvjyb55vk5/v2

Version created by Nastia Malochka

**KEYWORDS** 

cell preparation, DeMuxlet, human serum

**LICENSE** 

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

May 04, 2022

LAST MODIFIED

Jul 27, 2022



1

#### **OWNERSHIP HISTORY**

May 10, 2022 Shvetha Sankaran

PROTOCOL INTEGER ID

61989

#### **GUIDELINES**

- If cell viability for any sample is <50%, avoid using that sample in the suspension mix.
- If cell viability is between 50-70%, you may try to enrich for viable cells by centrifuging at lower speed (100 g). It may improve the viability to >75%, the cells can then be mixed with the other samples.

MATERIALS TEXT

#### List of reagents and materials

• Fisher Catalog #21870076

**X** Human Serum **Sigma** 

Aldrich Catalog #H4522

▼ Fetal Bovine Serum Sigma

Aldrich Catalog #F2442

⊠L-Glutamine (200 mM) Gibco - Thermo

■ Fischer Catalog #25030081

⊠ Penicillin-Streptomycin (10,000 U/mL) Thermo Fisher

Scientific Catalog #15140122

**⊠** Gibco™ DPBS no calcium no magnesium **Thermo Fisher** 

Scientific Catalog #14190144

BSA (Capricorn Scientific; Cat. No.: BSA-1S)

**⊠** Axygen<sup>™</sup> 1000 μL Universal Pipetter Tips: Wide Bore **Fisher** 

Scientific Catalog #14-222-703

Scientific Catalog #15250061

MACS® SmartStrainers (30 μm) Miltenyi

- Biotec Catalog #130-110-915
- EVE Cell Counting slides (NanoEnTek, Cat. No. 10027-446)
- QIAamp DNA Mini Kit (Qiagen, Cat no. 51306)



2

**SQIAamp DNA Mini Kit** 

(250) Qiagen Catalog #51306

- Inc. Catalog #20030770
- An appropriate single cell reagent kit and instrument for your application, e.g., 10x Genomics
  Chromium Controller

#### SAFETY WARNINGS

Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

1m

### Preparation of Reagents and Media

- 1 Prepare appropriate volume of <u>thawing media</u> (RPMI + 5% HS + 1% Pen/Strep + 1% Glutamine) and keep it at § 4 °C.
- 2 Prepare appropriate volume of <u>wash media</u> (RPMI + 10% FBS + 1% Pen/Strep + 1% Glutamine) and and keep it at **§ 4 °C**.
- 3 Prepare appropriate volume of fresh PBS + 0.04% BSA.

Thawing Frozen PMBCs and Preparing the Suspension Mix

4 /

Warm up thawing media, wash media and PBS + 0.04% BSA in § 37 °C water bath.

5

Transfer **9 mL** of **37 °C** pre-warmed thawing media into each of the 15 mL Falcon tube.

Take the cryovial containing PBMCs out of liquid nitrogen storage, place cryovial on dry ice, and <sup>2m</sup> transfer immediately to the § 37 °C water bath. Thaw for 1 minute- © 00:02:00 until no visible ice crystals remain.

7 / 2m

After thawing for 1 minute-  $\bigcirc$  **00:02:00**, open the cryovial in a biosafety cabinet and add  $\square$  **500**  $\mu$ L -  $\square$  **1** mL of pre-warmed thawing media into the cryovial using the  $\square$  **1** mL wide-bore blue tips.

8

After adding the thawing media, use the **1 mL** wide-bore blue tips to gently transfer the whole suspension from the cryovial into the 15 mL Falcon tube containing **9 mL** of pre-warmed thawing media.

9 💢

Mix the suspension extremely gently by inverting the Falcon tube twice.

10

Centrifuge at **300 x g, 21°C, 00:05:00**.

- 11 Decant the supernatant.
- 12 Leave around **■200 µL** of supernatant behind.
- 13

Using a serological pipette, gently re-suspend the cell pellet (3-5 times) in **5 mL** of pre-warmed wash media.

14

Centrifuge at **300 x g, 21°C, 00:05:00**.

Decant the supernatant. Leave around **■200 µL** of supernatant behind.

16

Using a serological pipette, gently re-suspend the cell pellet (7 times) in  $\blacksquare 3$  mL of pre-warmed PBS + 0.04 % BSA.

Avoid bubbles.

17



5m

Repeat Step 14-16.

Centrifuge at 300 x g, 21°C, 00:05:00 . Decant the supernatant. Leave around  $\blacksquare$ 200  $\mu$ L of supernatant behind. Using a serological pipette, gently re-suspend the cell pellet (7 times) in  $\blacksquare$ 3 mL of pre-warmed PBS + 0.04 % BSA.

Avoid bubbles.

18



5m

After the second wash, centrifuge at  $\$300 \times g$  for \$00:05:00 at  $\$21 \degree C$ . Decant the supernatant until around  $$\square200 \mu L$  of supernatant is left behind.

- 19 Gently re-suspend the cells in  $\blacksquare 800 \ \mu L$  of PBS + 0.04% BSA (pipette ~10 times).
- Filter the cell suspension through the  $\rightarrow$  Macs SmartStrainer to remove clumps or debris. After filtering, keep cells & On ice .



5

## 21 🗸 🎝

Thoroughly mix  $\Box 15 \mu L$  of cell suspension with  $\Box 15 \mu L$  of trypan blue. Load  $\Box 10 \mu L$  of mixture into each of the two chambers of a cell counting slide.

- Let the samples sit in the cell counting slide for © **00:01:00** before performing cell counting on an automated cell counter.
- 23 Make aliquots of  $1.50 \times 10^6$  cells/mL for each sample (  $\blacksquare 100 \, \mu L$  aliquot per sample).
- 24 Keep the remaining cell suspension from individual samples & **On ice** for DNA extraction using a QIAamp DNA Mini Kit (Qiagen, Cat no. 51306) according to the manufacturer's protocol. Perform genotyping using Illumina Global Screening Array-24 v3.0 BeadChip (Cat. No.: 20030770).

# 25

Mix equal volumes (  $\blacksquare 80~\mu L$  ) of cell suspension from each sample (up to 16 samples) to make a pooled suspension with a final concentration of  $1.50 \times 10^6$  cells/mL (total volume: 1280  $\mu L$ ). Keep this pooled suspension § On ice.

# 26

Thoroughly mix  $\Box 15 \mu L$  of the pooled suspension with  $\Box 15 \mu L$  of trypan blue. Load  $\Box 10 \mu L$  of the mixture into each of the two chambers of a cell counting slide.

- Let the samples sit in the cell counting slide for © 00:01:00 before performing cell counting on an automated cell counter.
- Count cells from the pooled suspension using an automated cell counter and verify that the concentration is in the range of  $1.10 \times 10^6$  cells/mL to  $1.50 \times 10^6$  cells/mL.
- Proceed with single cell capturing: Load 40,000 cells per 10x well using an appropriate 10x reagent kit for your application and run the 10x Genomics Chromium Controller according to the manufacturer's protocol.