




Nov 03, 2022

# Immunofluorescent staining of LN

Guido Krähenbühl<sup>1</sup><sup>1</sup>MDC-Berlin1 *Works for me* Share

This protocol is published without a DOI.

 Guido Krähenbühl

## ABSTRACT

Generic Staining Protocol

## PROTOCOL CITATION

Guido Krähenbühl 2022. Immunofluorescent staining of LN. **protocols.io**  
<https://protocols.io/view/immunofluorescent-staining-of-ln-cisquedw>



## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Nov 03, 2022

## LAST MODIFIED

Nov 03, 2022

## PROTOCOL INTEGER ID

72240

1 Thaw Tissue sections at RT

2 wash in PBS for 10 minutes

 00:10:00

10m

3 Circle samples with PAP-Pen

30m

- 4 Block in 10% Serum for 30 minutes  
🕒 00:30:00
- 5 Add primary Antibody for 90 minutes at RT or o/n at 4C
- 6 wash 3x with PBS?
- 7 incubate with secondary AB for 60 minutes
- 8 wash 3x with PBS?
- 9 Add DABCO-Mowiol on samples to increase IF
- 10 Mikroskopie!