

6



Oct 09, 2022

IPTG_induction

Ana Belem García González¹, Georgina Diego¹, Jair Alexis Gardea Sáenz¹, Irán Alessandra Chaparro Rodríguez¹

¹Tecnologico de Monterrey Campus Chihuahua





dx.doi.org/10.17504/protocols.io.4r3l2od1pv1y/v1



ABSTRACT

Protocol to induce protein expression using IPTG.

DOI

dx.doi.org/10.17504/protocols.io.4r3l2od1pv1y/v1

PROTOCOL CITATION

Ana Belem García González, Georgina Diego, Jair Alexis Gardea Sáenz, Irán Alessandra Chaparro Rodríguez 2022. IPTG_induction. **protocols.io** https://dx.doi.org/10.17504/protocols.io.4r3l2od1pv1y/v1

LICENSE

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 08, 2022

LAST MODIFIED

Oct 09, 2022

PROTOCOL INTEGER ID

68367

Protocol 5h

From a relatively fresh plate, grab a colony and grow © Overnight at § 37 °C and \$\to\$225 rpm in \$\boxed{10}\$ mL of LB+CAM.



1

Citation: Ana Belem GarcÃÂa González, Georgina Diego, Jair Alexis Gardea Sáenz, Irán Alessandra Chaparro RodrÃÂguez IPTG_induction https://dx.doi.org/10.17504/protocols.io.4r3l2od1pv1y/v1

- 2 Inoculate □3 mL in □30 mL of LB+CAM and incubate at 8 37 °C and ⊕225 rpm
- 3 When the cells reach to a OD600 between 0.4 and 0.8, the protein expression must be induced (IPTG concentrations can vary from 0.1 to 1M).

5h

4 Incubate a § 37 °C at **225 rpm** for **35:00:00** (Hours can vary).