



VERSION 2

JAN 12, 2023

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.e6nvwkkewvmk/v2

External link:
<https://doi.org/10.7554/eLife.79208>

Collection Citation: Hanqin Li, Oriol Busquets, Steven Poser, Dirk Hockemeyer, Frank Soldner 2023. Nucleofection (Amaxa) and electroporation (Biorad) of hPSCs. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvwkkewvmk/v2> Version created by [Devin E Snyder](#)

MANUSCRIPT CITATION:
 Hanqin Li, Oriol Busquets, Yogendra Verma, Khaja Mohieddin Syed, Nitzan Kutnowski, Gabriella R Pangilinan, Luke A Gilbert, Helen S Bateup, Donald C Rio, Dirk Hockemeyer, Frank Soldner (2022) Highly efficient generation of isogenic pluripotent stem cell models using prime editing eLife 11:e79208

<https://doi.org/10.7554/eLife.79208>

🌐 Nucleofection (Amaxa) and electroporation (Biorad) of hPSCs V.2

Hanqin Li¹, Oriol Busquets², Steven Poser², Dirk Hockemeyer¹, Frank Soldner¹

¹University of California, Berkeley; ²Albert Einstein College of Medicine



Devin E Snyder

ABSTRACT

This collection describes the standard procedure for the delivery of plasmids, mRNA or ribonucleoprotein (RNP) into human pluripotent stem cells (hPSCs) using nucleofection or electroporation.

Collection Overview

Preparing MEF-cultured hPSCs for nucleofection

Preparing Feeder-free hPSCs for nucleofection

Preparing plasmids for nucleofection

Preparing mRNA for nucleofection

In vitro assembling of RNP for nucleofection

Nucleofection of hPSCs

Electroporation of hPSCs

A. Preparation of MEFs-cultured hPSCs for electroporation

B. Preparing plasmids for electroporation

C. Electroporation

FACS-based enrichment of transfected cells

A. Preparation of samples for FACS sorting

B. After FACS

General notes

1. Throughout these protocols, the term hPSC is used to collectively refer to both hiPSCs and hESCs. All described procedures have been tested and work equally well for hiPSCs and hESCs.

MATERIALS

Item	Vendor	Catalog #
DMEM/F12	Thermo Fisher	11320082
DPBS w/o calcium and magnesium	Corning	MT21031CV

License: This is an open access collection distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this collection and it's working

Created: Jan 12, 2023

Last Modified: Jan 12, 2023

COLLECTION integer ID:
75230

Keywords: ASAPCRN

Item	Vendor	Catalog #
Fetal Bovine Serum (FBS)	Corning	35-011-CV
Knockout Serum Replacement	Thermo Fisher	10828-028
L-Glutamine	Sigma	G8540
Penicillin & Streptomycin (100X)	Thermo Fisher	15140163
MEM Non-Essential Amino Acids (100X)	Thermo Fisher	11140050
Heat Stable Recombinant Human FGF2	Thermo Fisher	PHG0360
Collagenase type IV	Thermo Fisher	17104019
2-Mercaptoethanol	Sigma	M3148
mTeSR-plus	STEMCELL Technologies	100-0276
StemFlex	Thermo Fisher	A3349401
Vitronectin (VTN-N) Recombinant Human Protein, Truncated	Thermo Fisher	A14700
Accutase	Thermo Fisher	SCR005
Dispase	STEMCELL Technologies	NC9995391
Y-27632	Chemdea	CD0141
Cas9, purified protein, 40uM	Macrolab, QB3 UC Berkeley	
Synthetic pegRNAs	IDT or Synthego	
Synthetic sgRNAs	Synthego	
P3 primary Cell 4D X kit S	Lonza	V4XP-3032
Countess™ Cell Counting Chamber Slides	Thermo Fisher	C10228
pCMV-PE2	Addgene	132775
4D-Nucleofector TM Core + X Unit	Lonza	AAF-1002B, AAF-1002X
5 ml polystyrene round-bottom tube with cell-strainer cap	Corning	352235
Cell-strainer (70 µm)	Fisher	07201431
Gene Pulser Xcell Eukaryotic System	Bio-Rad	1652661
Gene Pulser Electroporation Cuvettes, 0.4 cm gap	Bio-Rad	1652081

Item	Vendor	Catalog #
Exact N Amp Blood PCR Kit	Sigma	XNAB2-1TK

FILES

Protocol



NAME

Preparing MEF-cultured hPSCs for nucleofection

VERSION 2

CREATED BY



Devin E Snyder

OPEN →

Protocol



NAME

Preparing feeder-free hPSCs for nucleofection

VERSION 1

CREATED BY



Devin E Snyder

OPEN →

Protocol



NAME

Preparing plasmids for nucleofection of hPSCs

VERSION 1

CREATED BY



Devin E Snyder

OPEN →

Protocol



NAME

Preparing mRNA for nucleofection of hPSCs

VERSION 1

CREATED BY



Devin E Snyder

[OPEN](#) →

Protocol



NAME

In vitro assembling of RNP for nucleofection of hPSCs

VERSION 1

CREATED BY



Devin E Snyder

[OPEN](#) →

Protocol



NAME

Nucleofection of hPSCs

VERSION 1

CREATED BY



Devin E Snyder

[OPEN](#) →

Protocol



NAME

Electroporation of hPSCs

VERSION 1

CREATED BY



Devin E Snyder

[OPEN](#) →

Protocol



NAME

FACS-based enrichment of transfected hPSCs

VERSION 1

CREATED BY



Devin E Snyder

[OPEN](#) →