

JUN 16, 2023

Intraperitoneal L-dopa injection in mice

Zachary

natalia.lopezgonzalezdelrey¹, Gaertner¹

¹Northwestern University, Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815



natalia.lopezgonzalezdelrey

ABSTRACT

Intraperitoneal L-dopa injection in mice

OPEN BACCESS

DOI:

dx.doi.org/10.17504/protocol s.io.81wgbxd51lpk/v1

Protocol Citation: natalia.lo pezgonzalezdelrey, Zachary Gaertner 2023. Intraperitoneal L-dopa injection in mice. protocols.io https://dx.doi.org/10.17504/p

rotocols.io.81wgbxd51lpk/v1

License: This is an open access protocol distributed under the terms of the Creative Commons
Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: Jun 12, 2023

Last Modified: Jun 16, 2023

PROTOCOL integer ID: 83209

Keywords: ASAPCRN, Ldopa, Parkinson's disease,

benserazide

Preparation

1h

Mice were acclimated to the room for 01:00:00 1



L-dopa preparation

1h

2 Prepare the L-dopa mix fresh every day. Prepare a master mix of: 12.5mg/kg L-dopa and 12.5mg/kg benserazide dissolved in sterile saline with 0.25% citrate (wt/vol). Calculate grams depending on number of mice to be injected. Using vortex will dissolve the solution easier.

L-dopa injection

5m

3 Mice are injected intraperitoneally with the 12.5mg/kg L-dopa and 12.5mg/kg benserazide solution.

Behaviour

4 One hour post injection, mice were tested using the appropriate behavior protocol.