



APR 14, 2023

Transformation of Arabidopsis Thaliana Protocol

creative-biogene¹

¹Creative Biogene



creative-biogene

ABSTRACT

In this experiment, the flower immersion method was used to transfer the target gene into Arabidopsis thaliana using Agrobacterium mediated.

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.261ge3b7wl47/v1

Protocol Citation: creative-biogene 2023. Transformation of Arabidopsis Thaliana Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261ge3b7wl47/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
 We use this protocol and it's working

Created: Apr 14, 2023

Last Modified: Apr 14, 2023

PROTOCOL integer ID:
 80506

Experiment Summary

- 1 In this experiment, the flower immersion method was used to transfer the target gene into Arabidopsis thaliana using Agrobacterium mediated.

Main Reagents

- 2 YEB liquid medium, LB medium, 0.1 M CaCl₂, 0.05 M MgSO₄, flower immersion buffer (0.5XMS, 5% sucrose, 0.03% Silwet L-77), Rif, Kan.

Main Equipment

- 3 Shaker, centrifuge, culture bowl, greenhouse, trays, plastic film.

Experimental Steps

- 4
1. Preparation and transformation of *Agrobacterium tumefaciens* receptor cells
 - (1) Pick a single colony of GV1301 in 10 ml of YEB or LB liquid medium and incubate with shaking at 28°C until late logarithmic phase.
 - (2) Add 0.5 ml of bacterial solution to 50 ml of fresh YEB or LB liquid medium and incubate at 28°C with shaking until OD₆₀₀=0.5.
 - (3) Transfer to 50 ml centrifuge tube and ice bath for 20 min.
 - (4) Transfer to a 50 ml centrifuge tube and ice bath for 20 min.
 - (5) Centrifuge at 4000 rpm for 10 min at 4°C and collect the organisms.
 - (6) Resuspend the precipitate with 8 ml of pre-cooled 0.1 M CaCl₂ and 0.05 M MgSO₄.
 - (7) Centrifuge at 4000 rpm for 10 min at 4°C and collect the bacteriophage.
 - (8) Take 200 µl of receptor cells, add 0.5 µg of plasmid and mix gently.
 - (9) Freeze in liquid nitrogen for 50 s, place on ice for 40 min, and heat-shock at 42°C for 90 s.
 - (10) Add 800 µl of LB medium and incubate for 1 hour at 28°C.
 - (11) Apply the appropriate amount of cells on the screening medium and blow dry the surface liquid on the ultra-clean table.
 - (12) Incubate at 28°C for two days.
 2. Transformation of *Arabidopsis thaliana* by flower immersion
 - (1) Inoculate *Agrobacterium* clones containing expression vectors in 5 ml YEB medium (containing 100 µg/ml Rif, 100 µg/ml Kan) at 28°C, 200 rpm and incubate with shaking overnight.
 - (2) Transfer to 200 ml YEB medium at a ratio of 1:50 and incubate at 28°C, 200 rpm for 5 hours; centrifuge at 5000 X g for 15 min and collect the bacteriophage; resuspend in floral immersion buffer (0.5XMS, 5% sucrose, 0.03% Silwet L-77) and adjust OD₆₀₀ to 0.8.

(3) Place the Arabidopsis culture bowl upside down on a suitable size container containing flower soaking buffer, soak for 3-5 min, remove the culture bowl upside down in a tray, cover the tray with plastic film, remove the film after 24 hours, and continue the culture in the greenhouse.

Experimental Materials

- 5 Constructed expression vector, *Arabidopsis thaliana*.