

Feb 16, 2021

Immunofluorescence staining of PFA or fresh frozen mouse brain section

Daniel Manrique-Castano¹¹Université Laval

1

Works for me

This protocol is published without a DOI.

Daniel Manrique-Castano
Université Laval

SUBMIT TO PLOS ONE

ABSTRACT

This protocol is suitable for immunofluorescence staining of PFA fixed or Fresh Frozen mouse brain sections.

PROTOCOL CITATION

Daniel Manrique-Castano 2021. Immunofluorescence staining of PFA or fresh frozen mouse brain section.
protocols.io
<https://protocols.io/view/immunofluorescence-staining-of-pfa-or-fresh-frozen-bsf6nbre>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 16, 2021

LAST MODIFIED

Feb 16, 2021

PROTOCOL INTEGER ID

47326

MATERIALS TEXT

Blocking-Permeabilization buffer

 1X PBS (Phosphate-buffered saline) Contributed by users

[M] 0.1 % (v/v)  Triton-X100 Contributed by users


 Tween 20 Bio-rad

[M] 0.05 % (v/v) Laboratories Catalog #170-6606-MSDS

 Glycine Sigma –

[M] 0.3 Molarity (M) Aldrich Catalog #G8898

Primary antibody buffer

 1X PBS (Phosphate-buffered saline) Contributed by users

 Tween 20 Bio-rad

[M] 0.1 % (v/v) Laboratories Catalog #170-6606-MSDS

 Normal Goat Serum Gibco - Thermo

[M] 5 % (v/v) Fischer Catalog #LSPCN5000

or

 Normal Donkey Serum Contributed by

users Catalog #017-000-121

depending on your secondary antibody

species.

Secondary antibody buffer

 1X PBS (Phosphate-buffered saline) Contributed by users

 Normal Goat Serum Gibco - Thermo

[M] 1 % (v/v) Fischer Catalog #LSPCN5000

 Normal Goat Serum Gibco - Thermo

[M] 1 % (v/v) Fischer Catalog #LSPCN5000

or

 Normal Donkey Serum Contributed by

users Catalog #017-000-121

depending on your secondary antibody

species.



SAFETY WARNINGS




Paraformaldehyde (PFA) and DAPI are highly toxic substances. Manipulation should be performed carefully and according to security measurements.

BEFORE STARTING

Please read the whole protocol before starting the procedure.

Rehydration and permeabilization 1h 10m

- 1 If fixed sections are kept in  -80 °C , rehydrate them before starting the staining procedure by incubating them in ^{10m} PBS for  00:10:00 .

For fresh frozen sections kept at  -80 °C , incubate with [M] 4 % volume PFA during  00:15:00 followed by 3 washes in PBS,  00:05:00 each.

- 2 To permeabilize and block for unspecific staining, decant the PBS and incubate the section with **Blocking-Permeabilization buffer** for ⌚ **01:00:00** at 🌡 **Room temperature** in a humid chamber/box to prevent sections from drying out. 1h

If required, encircle brain sections using a hydrophobic pen previous incubation to prevent buffer leakage. Ensure that all sections are well covered with buffer.


Antibody incubation 1h 10m

- 3 When permeabilization is achieved, decant the buffer and incubate the sections in **primary antibody buffer** containing **diluted primary antibodies** ⌚ **Overnight** at 🌡 **4 °C**.

- 4  5m

When incubation time is finished, wash the sections 3 times with **PBS**, ⌚ **00:05:00** each.

- 5 Incubate the cells in **secondary antibody buffer** containing **diluted secondary antibodies** ⌚ **01:00:00** at 1h 🌡 **Room temperature**

- 6  5m

Subsequently, wash the cells 3 times with **PBS**, ⌚ **00:05:00** each.

- 7 To mount the sections, Remove as much of the remaining buffer as possible from the crystals, and with a pipette tip  **Fluoromount-G Electron Microscopy** drop 1-2 small drops of **Sciences Catalog #17984-25**. Subsequently, cover the sections with appropriate glass coverslips.

- 8 Let section air dry in the dark before imaging.