

Version 2

Feb 14, 2021

Recombinant expression and purification of codon-optimized M-MLV and Mashup V.2

Maira Rivera¹, Javiera Reyes¹, Paula Blazquez-Sanchez¹, Cesar A Ramirez-Sarmiento¹¹Institute for Biological and Medical Engineering, Pontificia Universidad Católica de Chile

1 Works for me dx.doi.org/10.17504/protocols.io.bsernbd6

Reclone.org (The Reagent Collaboration Network)

Tech. support email: protocols@reclone.org[Click here to message tech. support](#)

Cesar A Ramirez-Sarmiento

Pontificia Universidad Catolica de Chile

SUBMIT TO PLOS ONE

ABSTRACT

This protocol has been optimized for the recombinant expression of codon-optimized M-MLV and MashUp RT.

The goal of this protocol was to eliminate the use of large volumes for dialysis and potential issues with the protein crashing out of the solution due to the use of concentrators for buffer exchange of this enzyme into storage conditions.

M-MLV RT has been further functionally tested and works as expected. The sequence of the plasmid encoding the codon-optimized M-MLV RT enzyme can be found [here](#).

DOI

dx.doi.org/10.17504/protocols.io.bsernbd6

PROTOCOL CITATION

Maira Rivera, Javiera Reyes, Paula Blazquez-Sanchez, Cesar A Ramirez-Sarmiento 2021. Recombinant expression and purification of codon-optimized M-MLV and Mashup. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bsernbd6>

Version created by Cesar A Ramirez-Sarmiento

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 14, 2021

LAST MODIFIED

Feb 14, 2021

PROTOCOL INTEGER ID

47281

MATERIALS TEXT

MATERIALS

[Sodium phosphate monobasic monohydrate](#) Sigma

Aldrich Catalog #S9638

 **PMSF Sigma**

Aldrich Catalog #P7626

 **Sodium phosphate dibasic Sigma**

Aldrich Catalog #7558-79-4

 **DTT Sigma**

Aldrich Catalog #D0632

 **Imidazole Sigma Catalog #I5513**

 **NaCl Sigma**

Aldrich Catalog #53014

 **HiTrap Heparin HP affinity column Ge Life**

Sciences Catalog #17040701

 **HisTrap FF Crude Column Ge**

Healthcare Catalog #17528601

 **Glycerol Sigma**

Aldrich Catalog #G5516

 **Dextrose Sigma –**

Aldrich Catalog #D9434

 **Tween-20 Sigma**

Aldrich Catalog #P9416

 **Nonidet P40 Substitute Sigma**

Aldrich Catalog # 74385

 **EDTA Sigma**

Aldrich Catalog #ED2SS

Buffer A, pH 8.0

[M]**50 Milimolar (mM)** NaPO₄, pH 8.0

[M]**50 Milimolar (mM)** dextrose

[M]**300 Milimolar (mM)** NaCl

[M]**1 Milimolar (mM)** EDTA

[M]**0.1 % volume** Nonidet P-40

[M]**0.1 % volume** Tween-20

[M]**40 Milimolar (mM)** Imidazole, pH 8.0

Buffer B, pH 8.0

[M]**50 Milimolar (mM)** NaPO₄, pH 8.0

[M]**300 Milimolar (mM)** NaCl

[M]**1 Milimolar (mM)** EDTA

[M]**0.1 % volume** Nonidet P-40

[M]**0.1 % volume** Tween-20

[M]**10 % volume** Glycerol

[M]**150 Milimolar (mM)** Imidazole, pH 8.0

Buffer HA, pH 8.0

[M]50 Milimolar (mM) Tris-HCl, pH 8.0

[M]100 Milimolar (mM) NaCl

[M]0.1 % volume Nonidet P-40

Buffer HB, pH 8.0

[M]50 Milimolar (mM) Tris-HCl, pH 8.0

[M]2000 Milimolar (mM) NaCl

[M]0.1 % volume Nonidet P-40

Storage Conditions

[M]25 Milimolar (mM) Tris-HCl, pH 8.0

[M]150 Milimolar (mM) NaCl

[M]0.1 Milimolar (mM) EDTA



[M]0.1 % volume Nonidet P-40

[M]5 Milimolar (mM) DTT

[M]50 % volume Glycerol



DAY 1 – Plasmid transformation

1d

- 1 Transform  100 ng of plasmid containing codon-optimized M-MLV or Mashup into *E. coli*/BL21(DE3) competent cells using either heat shock or electroporation. 2h
- 2 Spread transformed cells in LB Agar plates supplemented with [M]0.05 mg/ml Kan for Mashup, or [M]0.1 mg/ml Amp for M-MLV. Grow plate overnight at  37 °C. 12h




DAY 2 – Preinoculum

1d

- 3 Select a single colony from the LB agar plate to prepare a preinoculum in  10 mL LB media supplemented with [M]0.05 mg/ml Kan for Mashup or [M]0.1 mg/ml Amp for M-MLV. Grow overnight at  250 rpm, 37°C. 1d




DAY 3 – Protein Overexpression

1d

- 4 Use the full volume of the preinoculum to inoculate  1 L of LB media supplemented with [M]0.05 mg/ml Kan for Mashup or [M]0.1 mg/ml Amp for M-MLV (1% inoculation). Grow at  250 rpm, 37°C until reaching an optical density at 600 nm (OD₆₀₀) = 0.8. 4h
- 5 Upon reaching OD₆₀₀ = 0.8, add [M]0.5 Milimolar (mM) IPTG and incubate overnight at  160 rpm, 18°C. 16h

DAY 4A – Protein Purification by IMAC

3h

- 6 Centrifuge the cell culture  4000 x g, 4°C, 00:20:00. Then, resuspend the cell pellet in  50 mL of Buffer A freshly supplemented with [M]0.5 Milimolar (mM) PMSF and [M]0.2 mg/ml lysozyme. 30m
- 7 Incubate the resuspended cells  80 rpm, Room temperature, 00:30:00. 30m

- 8 Sonicate on ice for **00:04:00** using cycles of **00:00:01** ON and **00:00:06** OFF at 40% amplitude (Qsonica Q125, 125W). 30m
- 9 Centrifuge the unclarified lysate **20000 x g, 4°C, 00:20:00** and collect the supernatant. You might want to collect a small sample for SDS-PAGE afterwards. 30m
- 10 On a **1 mL HisTrap column** preequilibrated with 10 column volumes (c.v.) (here, 10 mL) of **Buffer A**, load the supernatant. Wash with 10-20 c.v. of **Buffer A**. Then, elute with 5 c.v. of **Buffer B**, collecting the eluted fractions every **0.5 mL** in 1.5 ml tubes. 1h
- 11 To quickly pool the fractions containing the protein of interest, prepare a 96-well plate or 1.5 mL tubes with **40 µl** of Bradford reagent and **160 µl** of distilled water. Then, add **10 µl** of each protein fraction and compare against a blank reference sample corresponding to **10 µl** of **Buffer B**. You can determine your protein-containing fractions either by absorbance at 595 nm on a plate reader or visually by comparing the blue coloration of each fraction against the blank reference. Pool your fractions and collect a **10 µl** sample for SDS-PAGE 5m

DAY 4B – Second purification and buffer exchange by Heparin

2h


- 12 Dilute the pooled fractions 3-fold in **50 Milimolar (mM) Tris-HCl pH 8.0**, such that the final concentration of NaCl is 100 mM. 5m
- 13 Load the diluted sample onto a **1 ml HiTrap Heparin column** previously equilibrated with 10 c.v. (here, 10 mL) **Buffer HA**. Then, elute the protein using a 10 c.v. linear gradient against **Buffer HB**, collecting the eluted fractions every **0.5 mL** in 1.5 mL tubes. This linear gradient can be achieved by connecting two containers, one with 5 c.v. **Buffer HA** and the other with 5 c.v. **Buffer HB**, using a syphon or tubing, and withdrawing solution from the **buffer HA** container to the column using a peristaltic pump or by gravity. The protein will elute at high concentrations between 200 and 400 mM NaCl. 1h
- 14 Determine your protein-containing fractions using the Bradford assay. Pool your fractions and determine its protein concentration using the same method. Also, collect a **10 µl** sample for SDS-PAGE. 10m
- 15 For storage, supplement your pooled fraction with **0.2 % volume** Nonidet P-40, **0.2 Milimolar (mM) EDTA** and **10 Milimolar (mM) DTT**. Then, add glycerol up to **100 % volume** to reach **Storage Conditions**. Do consider that a final protein concentration of **0.3 mg/ml** is appropriate for subsequent experiments. 5m
- 16 Generate **200 µl** aliquots of the enzyme and store it at **-20 °C** until required. 30m

IMAC SDS-PAGE Result


10m

17





The unlabelled lane corresponds to the PageRuler protein ladder. His = Pooled fraction from HisTrap IMAC column (20 µl). Numbered lanes correspond to eluted fractions from Heparin column (10 µl). The high molecular weight contaminations in M-MLV are due to disulfide bonds between the free cysteins in the protein.



Effect of reducing agents on the elimination of species with high molecular weight. M-MLV and Mashup were incubated overnight in oxidizing (Ox) or reducing conditions using either BME or DTT as the reducing agent. For the SDS-PAGE, the loading buffer was prepared with (+) or without (-) reducing agents. Condition Ox/- shows some high molecular weight species, which are eliminated upon addition of reducing agents in the protein-containing solution (BME, DTT). For M-MLV, 5 µg were loaded on each lane, whereas 2 µg were loaded for Mashup.