



Mar 03, 2022

# Culturing Primary Cortical Neurons

Haley Geertsma<sup>1</sup><sup>1</sup>University of Ottawa

1

[dx.doi.org/10.17504/protocols.io.b5u6q6ze](https://dx.doi.org/10.17504/protocols.io.b5u6q6ze) Haley Geertsma  
University of Ottawa

This protocol is used to isolate and culture primary cortical neurons from mouse embryos at E14.5-E15.5.

DOI

[dx.doi.org/10.17504/protocols.io.b5u6q6ze](https://dx.doi.org/10.17504/protocols.io.b5u6q6ze)

Haley Geertsma 2022. Culturing Primary Cortical Neurons. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.b5u6q6ze>



protocol ,

Mar 03, 2022

Mar 03, 2022

59006

Day 1 1d 3h 16m

- 1 Coat culture dishes with 0.05mg/mL Poly-D-lysine overnight at 37°C.

1d

Day 2 3h 14m

- 2 Wash culture dishes with sterile water x2 then leave to dry.

1m

- 3 Sedate pregnant mouse with 120mg/kg Euthanyl and confirm sedation prior to proceeding<sup>2m</sup> with the next step.

\*The neurons must be plated within approximately 4 hours of sedation.

- 4 Soak pregnant mouse with 70% ethanol then dissect open the abdominal cavity. 1m
- 5 Isolate uterus with embryos and place in a tube of 1X PBS. 1m
- 6 In a laminar flow hood, dissect the brains of E14.5-15.5 embryos in Hank's Balanced Salt Solution (HBSS) and isolate each brain in 700uL of HBSS. 2h
- 7 Gently pipette up/down x10 to mix cortices then add 20uL trypsin for every embryo and incubate on a rotator at 37°C for 20 minutes. 20m
- 8 Add 300uL Solution A per embryo and carefully pipette up/down x5 then centrifuge at 2500g for 5 minutes at 4°C. 10m  
Solution A: 2.75mL Neurobasal media (unsupplemented) + 150uL trypsin inhibitor + 100uL DNase1
- 9 Remove supernatant and carefully resuspend pellet in 300uL Solution B per embryo by pipetting up/down x10. 5m  
Solution B: 2.55mL Neurobasal media (unsupplemented) + 200uL trypsin inhibitor + 250uL DNase1
- 10 Centrifuge at 2500g for 5 minutes at 4°C then remove supernatant. 5m
- 11 Resuspend pellet in 1mL Neurobasal complete per embryo. 1m
- 11.1 Optional trypan blue cell exclusion assay: Mix 100uL trypan blue + 100uL 1X PBS + 20uL resuspended cells. Count with haemocytometer. 5m
- 12 Plate neurons on pre-coated dishes at desired confluency. 1m

- 13 Change neurobasal medium at 4 days *in vitro*. 5m
- 14 Fix neurons at 7 days *in vitro* by removing media and replacing with 4% paraformaldehyde for 10 minutes at room temperature. 10m
- 15 Wash neurons with 1X PBS x2 to remove traces of paraformaldehyde then store in 1X PBS at 4°C. 5m