

May 29, 2024

Quick protocol for protein extraction from adherent fish-derived fibroblasts

DOI

dx.doi.org/10.17504/protocols.io.q26g71ydkgwz/v1

Joao M Moreno^{1,2}, Vitor C Sousa¹, Romana Santos²

¹cE3c – Centre for Ecology, Evolution and Environmental Changes & CHANGE – Global Change and Sustainability Institute, Faculdade de Ciências da Universidade de Lisboa, Lisboa, Portugal;

²MARE – Centro de Ciências do Mar e do Ambiente (MARE) & ARNET—Aquatic Research Network, Faculdade de Ciências da Universidade de Lisboa, Lisboa, Portugal



Joao M Moreno

cE3c – Centre for Ecology, Evolution and Environmental Chang...

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.q26g71ydkgwz/v1

Protocol Citation: Joao M Moreno, Vitor C Sousa, Romana Santos 2024. Quick protocol for protein extraction from adherent fish-derived fibroblasts. [protocols.io https://dx.doi.org/10.17504/protocols.io.q26g71ydkgwz/v1](https://dx.doi.org/10.17504/protocols.io.q26g71ydkgwz/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 16, 2024

Last Modified: May 29, 2024

Protocol Integer ID: 99955

Keywords: Protein extraction, Protein purification, Cell culture

Funders Acknowledgement:

**Fundação para a Ciência e a
Tecnologia**

Grant ID:

SFRH/BD/143199/2019

**Fundação para a Ciência e a
Tecnologia**

Grant ID: PTDC/BIA-

EVL/4345/2021

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

We present a rapid and efficient protocol for extracting total proteins from adherent fish-derived fibroblasts, specifically optimized for applications in Western blotting and mass spectrometry. The method involves directly dissociating cells from the culture flask into a lysis buffer (e.g. RIPA buffer) followed by centrifugation to collect the total soluble proteins. The use of buffers like RIPA ensures comprehensive solubilization of cellular proteins while preserving their integrity and functionality. This straightforward and reproducible protocol yields high-quality protein extracts suitable for various downstream analytical techniques. Its simplicity and reliability make it an invaluable tool for researchers studying proteomics using fish cell culture.

Materials

Solutions

- 1X Phosphate Buffered Saline (PBS)
- RIPA lysis buffer: 50 mM Tris-HCl (pH 7.4), 150 mM NaCl, 1mM EDTA, 1% NP-40, 1% sodium deoxycholate, 0.1% SDS, sterile-filtered.

Materials

- 2 mL or 1.5 mL microcentrifuge tubes
- Cell scrappers
- Pipettes and pipette tips
- Serological pipettes (variable volume)
- Pasteur pipettes

Equipments

- Refrigerated microcentrifuge
- Flow hood chamber

- 1 Carefully remove all culture media from the flask and add enough ice-cold 1X PBS to wash the cells
- 2 Carefully remove the ice-cold 1x PBS and add ice-cold lysis buffer (RIPA buffer) according to the estimated number of cells:
 1 mL for 10^7 cells (roughly a T75 flask).
- 3 Use a cell scraper to dissociate the cells from the bottom of the flask.
- 4 Resuspend the cells in the lysis buffer and transfer the suspension to a microcentrifuge tube.
- 5 Agitate for 20 minutes at 4°C.
- 6 Centrifuge at 13,000 x g for 20 min at 4°C.
- 7 Carefully transfer the supernatant containing the soluble protein to a new tube and keep on the ice. Discard the pellet.

Note: The protein solution can be kept in the freezer for longer storage periods until further use.