



OCT 16, 2023

Expression and purification for cryo-EM samples of FIP200NTD:ATG13(363-517)-ULK1(MIT) complexes

Minghao Chen¹, Xuefeng Ren¹¹Department of Molecular and Cell Biology, University of California, Berkeley, Berkeley, CA 94720, USA

ASAP Collaborative Research Network



Minghao Chen

OPEN  ACCESS

ABSTRACT

Expression and purification for cryo-EM samples of FIP200NTD:ATG13(363-517)-ULK1(MIT) complexes

DOI:

dx.doi.org/10.17504/protocols.io.e6nvwjxw7lmk/v1

Protocol Citation: Minghao Chen, Xuefeng Ren 2023. Expression and purification for cryo-EM samples of FIP200NTD:ATG13(363-517)-ULK1(MIT) complexes.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.e6nvwjxw7lmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working







Created: Jun 05, 2023


Last Modified: Oct 16, 2023

PROTOCOL integer ID:
82914

Keywords: ASAPCRN


Funders
Acknowledgement:
Aligning Science Across
Parkinson’s (ASAP) initiative
Grant ID: ASAP-000350

Expression		2d 0h 40m
1	Transfect HEK GNTi cells at concentration of 2×10^6 cells/ml	
2	Dilute PEI with Warm Hybridoma-SFM(1X)	
3	In a separate tube, dilute DNA with Hybridoma-SFM(1X)	
4	Add PEI to DNA dilution. Incubate mixture for  00:30:00 at  37 °C	30m
5	Add mixture to cells. Let cells grow for  48:00:00	2d
6	Harvest Cells  500 rpm ,  4 °C ,  00:10:00	10m

7 Wash pellet with cold PBS. Store pellet at  -80 °C until purification.

Purification

2h 50m



8 Cell pellets were lysed at room temperature for 20 min with lysis buffer (25 mM HEPES pH 7.5, 200 mM NaCl, 2 mM MgCl₂, 1 mM TCEP, 10% Glycerol) with 5 mM EDTA, 1% Triton X-100 and protease inhibitor cocktail (Thermo Scientific) for  00:15:00

15m

9 Clarify lysate for  17000 rpm for  00:35:00 at  4 °C .

35m



10 Wash strep-tactin resin (IBA Lifesciences, Germany) into lysis buffer (without Triton). Load clarified lysate onto resin

11 Rock supernatant with equilibrated Strep-Tactin Sepharose resin for  01:00:00 at  4 °C


1h

12 Wash with 5CV lysis buffer (25 mM HEPES pH 7.5, 200 mM NaCl, 2 mM MgCl₂, 1 mM TCEP, 5 mM EDTA, 10% Glycerol)

13 Elute with lysis buffer plus 4 mM desthiobiotin for STREP resin

14 His6-TEV cleavage at  4 °C  Overnight ,

1h

- 15 Concentrate elution and inject onto pre-equilibrated Superose 6 Increase 10/300 GL column (Cytiva)
(25 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl₂, 1 mM TCEP)
- 16 Pool peak fractions, concentrate, snap freeze, and store at  -80 °C