



APR 12, 2024

## 🌐 Preparing fixed cells for immunofluorescence

Harper JW<sup>1,2</sup>, Louis R R Hollingsworth<sup>1</sup>, Chan Lee<sup>1</sup>

<sup>1</sup>Department of Cell Biology, Harvard Medical School Boston, MA 02115, USA;

<sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA

ASAP Collaborative Research Network



Louis R R Hollingsworth  
Harvard Medical School

### ABSTRACT

This protocol provides a general method to prepare fixed cells for immunofluorescence imaging, e.g., by confocal microscopy.

### PROTOCOL REFERENCES

<https://www.protocols.io/view/immunocytochemical-analysis-c7g8zjzw>

<https://www.protocols.io/view/immunofluorescence-of-galectin-3-puncta-after-lyso-x54v9jnmpg3e/v2>

OPEN ACCESS



DOI:

[dx.doi.org/10.17504/protocols.io.kxygyeeol8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygyeeol8j/v1)

**Protocol Citation:** Harper JW, Louis R R Hollingsworth, Chan Lee 2024. Preparing fixed cells for immunofluorescence.

**protocols.io**

<https://dx.doi.org/10.17504/protocols.io.kxygyeeol8j/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** Apr 09, 2024

Last Modified: Apr 12, 2024

PROTOCOL integer ID: 97980

Keywords: ASAPCRN, immunofluorescence, IF

Funders Acknowledgement:  
National Institutes of Health  
Grant ID: R01NS083524

MATERIALS

Cell culture -- sterile materials

- PBS (Thermo J67802.AP)
- 24-well glass bottom plate (Cellvis P24-1.5H-N)
- 0.1% poly-L-lysine solution (Sigma P8920)
- Relevant cell culture media

IF

- PBS (MedChemExpress HY-K1023)
- Triton-X 100 (Sigma X100)
- Tween-20 (Sigma 655205)
- 16% paraformaldehyde (PFA; Electron Microscopy Sciences 15710)
- BSA (GoldBio A-420-10)
- Glycine (Sigma 410225)
- Methanol (Sigma 34860-2L-R)

Prepare-ahead buffers (see protocol for more details)

- PBS (from MedChemExpress HY-K1023, filtered)
- PBS-T (PBS with 0.02% Tween-20)
- PBS-TritonX (PBS with 0.3%-0.4% Triton-X 100)
- PBS-PFA (PBS with 4% PFA; prepared fresh)
- PBS-Block (PBS with 3% BSA and 22.52 mg/mL)
- 100% methanol pre-chilled at -20 degrees celsius

Representative secondary antibodies and conjugates:

A	B	C
Secondary antibodies and conjugates		
Antobody	Vendor, Cat#	Dilution
Anti-mouse 488	Thermo A-11029	1:200
Anti-mouse 568	Thermo A-11004	1:200
Anti-mouse 647	Thermo A-21235	1:200
Anti-rat 488	Thermo A-11006	1:200
Anti-rat 647	Thermo A-21247	1:200
Anti-Rabbit 488	Thermo A-11034	1:200
Anti-Rabbit	Thermo A-11011	1:200

A	B	C
568		
Anti-Rabbit 647	Thermo A-21244	1:200
Phalloidin-405	Thermo A30104	1:200
Phalloidin-488	Thermo A12379	1:200

## Plate coating and cell seeding

- 1 Dilute 0.1% poly-L-lysine solution (Sigma P8920) in PBS (0.01%). Make enough solution for 0.5 mL per well, plus some excess.
- 2 Add 0.5 mL of the diluted poly-L-lysine solution (0.01% in PBS) to each well of a 24-well glass bottom plate (Cellvis P24-1.5H-N). Place the plate back in the cell culture incubator (37 °C, 5% CO<sub>2</sub>) for 1 – 24 h to adhere the solution to the plate.
- 3 Wash the wells 3x with sterile PBS or water.
- 4 Add 1 mL of cell culture medium to each well, then plate cells at an appropriate density such that they are not more than 70% confluent on the day of the experiment.

### Note

Use DMEM for iBMDMs and HEK293(T)s. For these cell types, seed the cells at approximately ~15% confluence for iBMDM or ~30% for HEK293Ts to have an ideal cell density for imaging the next day.

## Fixation and permeabilization

- 5 Prepare fresh solutions for IF including: PBS, PBST (PBS with 0.02% Tween-20), PBS-Triton (PBS with 0.3%-0.4% Triton-X 100) or 100% methanol pre-chilled at -20 degrees celsius, PBST-Block (see step 5.1 below), and PBS-PFA (see step 5.2 below).

#### Note

Different antibodies work with different fixation and permeabilization protocols. Additionally, certain organellar populations are better preserved with alternative fixation/permeabilization protocols. For example, we've found that lysosomes and other vesicles are better resolved with methanol permeabilization instead of PBS-TritonX (see note at step 8 for that particular alternative). Because of the difference in permeabilization temperatures, we recommend using different plates if you intend to conduct experiments in parallel that require different permeabilization protocols.

**5.1** To prepare PBST-Block, dissolve 3% BSA and 22.52 mg/mL glycine in PBST (for example, add 1.5 g BSA and 1.126 g glycine to 50 mL PBST). Glycine binds free aldehyde groups (from PFA) to prevent antibody crosslinking and thus high background, whereas BSA reduces nonspecific antibody binding.

**5.2** To prepare PBS-PFA, freshly dilute an ampule of 16% paraformaldehyde (PFA; Electron Microscopy Sciences 15710) in PBS to 4% (add the 10 mL ampule to 30 mL of PBS) in a falcon tube. Warm the solution to 37 degrees celsius prior to fixation.



#### Note

PFA is toxic, so do this step carefully in a chemical hood. Additionally, be careful when breaking the glass ampule.

**6** Conduct relevant cellular treatment(s). Stagger any treatments so that all of the wells of the plate are ready to harvest simultaneously.

**7** Aspirate the cell culture medium and add 1 mL of pre-warmed PBS-PFA to each well. Incubate for 10 mins **10m** with gentle rocking at room temperature.

#### Note

Certain cell types can be extremely labile, and so use caution when aspirating and adding fixation medium. Alternatively, you can add 8% PBS-PFA directly to the cell culture medium (1:1 volumetric ratio for a final concentration of 4%), which can prevent tricky cells from lifting off of the plate.

**8** For Triton-X 100 permeabilization, slowly add 1 mL of PBS-Triton directly to each well then immediately aspirate the solution (this can help prevent the formation of aggregates). Replace with 1 mL PBS-Triton and

incubate at RT for 10 min.

PFA-containing waste should be collected and disposed of properly.

#### Note

For methanol permeabilization, instead first wash the plate 3x with PBS, prior to adding 1 mL cold methanol (-20 degrees) to each well and placing the plate at -20 degrees celsius for 10 min.

- 9 Wash cells 3x with 1 mL PBST. For the first wash, aspirate immediately and replace with PBST, and then conduct the subsequent washes over 5-10 mins. 10m

## Blocking and antibodies

19h

- 10 Aspirate PBST, and add 0.5 mL PBST-Block to each well. Incubate at RT for 1 h. 1h

- 11 Dilute primary antibodies in pre-chilled PBST-Block on ice. Aspirate the wells and add 200 µL of the antibody solution to each well. Rock gently at 4 degrees celsius overnight. 18h

#### Note

Note that some organelles, like mitochondria, have significant autofluorescence in the 488 channel. Additionally, with multiple primary/secondary combinations, bleed through can occur. Control wells lacking particular primary antibodies (minus-one controls) should always be included at least once to assess bleed-through and autofluorescence with a given set of cells, antibodies, and laser settings.

A useful resource for checking dyes and fluorophores is FPbase: <https://www.fpbases.org/>

- 12 The next day, wash wells 3x with 1 mL PBST. For the first wash, aspirate immediately and replace with PBST, and then conduct the subsequent washes over 5-10 mins.

- 13** Dilute fluorescent secondaries (see materials) 1:200 in PBST-Block. Phalloidin can also be added at this step. Aspirate the wells and add 200  $\mu$ L of secondary antibody solution to each well. Rock gently at RT for 60-90 mins. 1h 30m

**Note**

Wrap the plate in aluminum foil to minimize fluorophore quenching over subsequent incubation and wash steps.

If possible, spread the secondary antibodies out as far as possible to avoid bleed-through (e.g., 488 and 647).

- 14** Wash cells 3x with 1 mL PBST. For the first wash, aspirate immediately and replace with PBST, and then conduct the subsequent washes over 5-10 mins. 10m

- 15** Aspirate the PBST and add Hoechst 33342 (0.5 mL, 5  $\mu$ M in PBS) to each well, if desired. Hoechst labels nuclei in the ~405 channel. Incubate at RT for 10 min. 10m

- 16** Wash cells 3x with 1 mL PBS. For the first wash, aspirate immediately and replace with PBST, and then conduct the subsequent washes over 5-10 mins. 10m

- 17** Keep cells at RT if imaging with the same day (ideal), or transfer the plate to 4 degrees celsius for imaging at a later time. Image the plate within ~3 - 5 d with an appropriate confocal microscope.

**Note**

Imaging settings highly depend on the cells and antibodies used in the experiment, in addition to the cellular phenomenon. We thus recommend consulting with an imaging core during imaging optimization.