



Sep 22, 2021

Preparation of Enhanced Orthogonal Aminoacyl-tRNA-Synthetase

In 1 collection

Anne Zemella¹, Theresa Richter^{1,2}, Lena Thoring¹, Stefan Kubick¹¹Cell-free and Cell-based Bioproduction, Fraunhofer Institute for Cell Therapy and Immunology (IZI), Branch Bioanalytics and Bioprocesses (IZI-BB), Potsdam, Germany;²University of Potsdam, Potsdam, Germany

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bqnumvew

Springer Nature Books

Satyavati Kharde
Springer Nature

ABSTRACT

This is part 3.1 of the "A Combined Cell-Free Protein Synthesis and Fluorescence-Based Approach to Investigate GPCR Binding Properties" collection of protocols:

<https://www.protocols.io/view/a-combined-cell-free-protein-synthesis-and-fluores-bqntmven>

Collection Abstract: Fluorescent labeling of de novo synthesized proteins is in particular a valuable tool for functional and structural studies of membrane proteins. In this context, we present two methods for the site-specific fluorescent labeling of difficult-to-express membrane proteins in combination with cell-free protein synthesis. The cell-free protein synthesis system is based on Chinese Hamster Ovary Cells (CHO) since this system contains endogenous membrane structures derived from the endoplasmic reticulum. These so-called microsomes enable a direct integration of membrane proteins into a biological membrane. In this protocol the first part describes the fluorescent labeling by using a precharged tRNA, loaded with a fluorescent amino acid. The second part describes the preparation of a modified aminoacyl-tRNA-synthetase and a suppressor tRNA that are applied to the CHO cell-free system to enable the incorporation of a non-canonical amino acid. The reactive group of the non-canonical amino acid is further coupled to a fluorescent dye. Both methods utilize the amber stop codon suppression technology. The successful fluorescent labeling of the model G protein-coupled receptor adenosine A2A (Adora2a) is analyzed by in-gel-fluorescence, a reporter protein assay, and confocal laser scanning microscopy (CLSM). Moreover, a ligand-dependent conformational change of the fluorescently labeled Adora2a was analyzed by bioluminescence resonance energy transfer (BRET).

For **Introduction and Notes**, please see: <https://www.protocols.io/view/a-combined-cell-free-protein-synthesis-and-fluores-bqntmven/guidelines>

DOI

dx.doi.org/10.17504/protocols.io.bqnumvew

EXTERNAL LINK

https://link.springer.com/protocol/10.1007/978-1-4939-9121-1_4

PROTOCOL CITATION

Anne Zemella, Theresa Richter, Lena Thoring, Stefan Kubick 2021. Preparation of Enhanced Orthogonal Aminoacyl-tRNA-Synthetase. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bqnumvew>

Zemella A., Richter T., Thoring L., Kubick S. (2019) A Combined Cell-Free Protein Synthesis and Fluorescence-Based Approach to Investigate GPCR Binding Properties. In: Tiberi M. (eds) G Protein-Coupled Receptor Signaling. Methods in Molecular Biology, vol 1947. Humana Press, New York, NY.
https://doi.org/10.1007/978-1-4939-9121-1_4

COLLECTIONS ⓘ



A Combined Cell-Free Protein Synthesis and Fluorescence-Based Approach to Investigate GPCR Binding Properties

KEYWORDS

Cell-free protein synthesis, G protein-coupled receptor, Protein modification, Non-canonical amino acids, Amber suppression, Confocal laser scanning microscopy

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 12, 2020

LAST MODIFIED

Sep 22, 2021

OWNERSHIP HISTORY

Dec 12, 2020	Lenny Teytelman	protocols.io
Jul 05, 2021	Emma Ganley	protocols.io
Jul 06, 2021	Emma Ganley	protocols.io
Aug 24, 2021	Satya K	
Aug 26, 2021	Satyavati Kharde	Springer Nature
Aug 26, 2021	Satyavati Kharde	Springer Nature
Sep 22, 2021	Satya K	
Sep 22, 2021	Satyavati Kharde	Springer Nature

PROTOCOL INTEGER ID

45492

PARENT PROTOCOLS

Part of collection

[A Combined Cell-Free Protein Synthesis and Fluorescence-Based Approach to Investigate GPCR Binding Properties](#)

2.1 Materials for Preparation of Enhanced Orthogonal Aminoacyl-tRNA-Synthetase

1. Coding sequence for the modified tyrosyl-tRNA-synthetase (eAzFRS, including the mutations Thr37, Ser182, Ala183, and Arg265 [11, 12] and a C-terminal Strep-Tag) from *E.coli*.
2. *E.coli* expression system (RTS 500 *E.coli*/HY Kit, biotechrabbit).
3. 100 mM Isopropyl β -D-1-thiogalactopyranoside (IPTG)
4. Gravity flow Strep-Tactin® superflow mini-column (0.2 ml).
5. Strep-Tactin® Purification Buffer Set: 10× Washing Buffer (1 M Tris-Cl, pH 8.0, 1.5 M NaCl, 10 mM EDTA), 10× Elution Buffer (1 M Tris-Cl, pH 8.0, 1.5 M NaCl, 10 mM EDTA, 25 mM Desthiobiotin) and 10× Regeneration Buffer (1 M Tris-Cl, 1.5 M NaCl, 10 mM EDTA, 10 mM HABA (hydroxyl-azophenyl-benzoic acid)).
6. Zeba™ Spin Desalting Columns (7 K MWCO, 0.5 ml).
7. Amicon® Ultra Centrifugal Filters (10 K device, 0.5 ml).
8. Synthetase storage buffer: 50 mM HEPES pH 7.6, 10 mM KOAc, 1 mM MgCl₂, 4 mM DTT.
9. Thermomixer with a microtiter plate adapter and a RTS 500 adapter.

SAFETY WARNINGS

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

3.1 Preparation of Enhanced Orthogonal Aminoacyl-tRNA-Synthetase

20h

- 1 For prokaryotic cell-free synthesis, the eAzFRS gene should be cloned into a vector containing a T7 promotor, ribosomal binding site, and T7 terminator such as pIX3.0, pIVEX2.3d, and pIVEX2.4d vectors or alternatively containing a T5 promotor such as pQE2 vectors as used in this protocol. eAzFRS is synthesized in a cell-free system using an *E. coli* lysate in a dialysis mode.

A typical 1.1 ml reaction is composed of **0.525 mL E.coli lysate** , **0.225 mL reaction mix** ,

0.27 mL amino acids without methionine , **30 μ L methionine** , **11 μ L IPTG** for the induction of the protein expression pQE2 vector, **39 μ L template** containing **110 μ g plasmid DNA** .

- 2 The surrounding feeding mixture contains **7990 μ L feeding mix** , **110 μ L IPTG** , **2650 μ L amino acids without methionine** and **300 μ L methionine** (see Note 2).

- 3 Fill the reaction solution into the reaction compartment (marked through the red lid).

- 4 Fill the feeding mix into the feeding chamber (marked through the colorless lid).

- 5 Insert the prepared chamber into the RTS 500 adapter in a thermomixer. The reaction time is **20:00:00** at **30 °C** and a shaking speed of **900 rpm** .

20h



- 6 For the separation of aggregated proteins from soluble eAzFRS a centrifugation step at **16000 x g, 4°C, 00:10:00** is recommended.



Equilibrate two Strep-Tactin columns with **400 µl 10× washing buffer** and add **500 µl supernatant of the cell-free reaction** to each column.

8 

After the supernatant has completely entered the column, wash each column 5× with **200 µl washing buffer** (see **Note 3**).

9 Elute the protein 6× with **100 µl elution buffer** and collect the fractions.

10 Elution fractions containing the target protein are pooled.

11 Regenerate the column with 3× **1 mL 1× regeneration buffer** and remove the regeneration buffer 2× with **800 µl 1× washing buffer**. Store the column in **2 mL washing buffer** at **4 °C**.

12  

The combined elution fractions are applied to Zeba™ Spin Desalting Columns to exchange the elution buffer of the strep-tag purification to a synthetase storage buffer. Therefore, remove the storage solution of the Zeba™ Spin Desalting Column by centrifugation at **1500 x g, 00:01:00**. Add **300 µl synthetase storage buffer** to the resin bed and centrifuge at **1500 x g, 00:01:00**. Repeat this **step 2×**.

Place the column in a new collection tube and apply **100 µl pooled synthetase solution** to each column.

Centrifuge at **2000 x g, 00:02:00** and collect the synthetase.

13  

The concentration of the synthetase can be performed with Amicon® Ultra Centrifugal Filters. Add up to **500 µl synthetase solution** to the concentrator and centrifuge at **14000 x g, 4°C, 00:10:00**. Collect the concentrated sample and determine the concentration by NanoDrop measurement using the molecular mass (48.6 kDa) and the extinction coefficient (54.3) (see **Note 4**).

14 

The synthetase can be stored at **-80 °C** after shock freezing in liquid nitrogen.