

# Rat Brain Tissue RNA Extraction/cDNA Synthesis for qPCR

Kokila Shankar<sup>1</sup>, Olivier George<sup>1</sup>

<sup>1</sup>University of California, San Diego

1 Works for me This protocol is published without a DOI.

George Lab  
Tech. support email: [olgeorge@ucsd.edu](mailto:olgeorge@ucsd.edu)

Kokila Shankar

## ABSTRACT

This protocol is to obtain cDNA from rat brain tissue micropunches for downstream qPCR. Briefly, RNA will be extracted from brain tissue and will be converted to cDNA through reverse transcription reaction. Protocols have been adapted from Zymo ([https://files.zymoresearch.com/protocols/\\_r2060\\_r2061\\_r2062\\_r2063\\_direct-zol\\_rna\\_microprep.pdf](https://files.zymoresearch.com/protocols/_r2060_r2061_r2062_r2063_direct-zol_rna_microprep.pdf)) and New England Biolabs (<https://www.neb.com/protocols/2013/01/23/first-strand-cdna-synthesis-protocols-e6560>).

## PROTOCOL CITATION

Kokila Shankar, Olivier George 2020. Rat Brain Tissue RNA Extraction/cDNA Synthesis for qPCR.  
**protocols.io**  
<https://protocols.io/view/rat-brain-tissue-rna-extraction-cdna-synthesis-for-7wdhpa6>

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Oct 02, 2019

## LAST MODIFIED

Jul 23, 2020

## PROTOCOL INTEGER ID

28325

## MATERIALS

NAME	CATALOG #	VENDOR
<a href="#">ProtoScript II First Strand cDNA Synthesis Kit - 150 rxns</a>	E6560L	New England Biolabs
<a href="#">Molecular Biology Grade Water</a>	10154604	Fisher Scientific
<a href="#">Ethanol (95 - 100%), molecular grade</a>		
<a href="#">TRIZOL reagent</a>	15596-026	Invitrogen - Thermo Fisher
<a href="#">Microcentrifuge tubes (1.5 ml)</a>	C-2170	Denville Scientific Inc.
<a href="#">Direct-zol RNA Microprep Kit</a>	R2062	Zymo Research

## STEPS MATERIALS

NAME	CATALOG #	VENDOR
<a href="#">TRIZOL reagent</a>	15596-026	Invitrogen - Thermo Fisher
<a href="#">Ethanol (95 - 100%), molecular grade</a>		
<a href="#">Molecular Biology Grade Water</a>	10154604	Fisher Scientific

## MATERIALS TEXT

Ensure all reagents are nuclease-free (not just sterile)

#### EQUIPMENT

NAME	CATALOG #	VENDOR
C1000 Touch Thermal Cycler	1851148	<a href="#">BioRad Sciences</a>
Nanodrop 2000C	TSC-ND2000C	<a href="#">Thermo Scientific</a>

#### SAFETY WARNINGS

Wear lab coat, gloves, and goggles as the protocol involves use of Trizol. Change gloves frequently to avoid contamination.

#### BEFORE STARTING

Make sure all surfaces/tools/equipment are sterile and/or have been cleaned with RNase decontaminant spray.  
Make sure centrifuge has been pre-cooled to 4°C. Thaw tissue samples on ice right before starting protocol.



1 Remove tissue samples from  **-80 °C** and thaw on ice

2 If samples are stored in RNALater ICE, remove from tube and add  **300 µl**



#### TRIZOL reagent

by Invitrogen - Thermo Fisher  
Catalog #: 15596-026

to each sample. Pipet homogenize until tissue is fully dissolved and leave at RT while working on other samples. Spin samples at 13000 RPM at  **4 °C** for  **00:05:00** and transfer supernatant to new tube.








If samples were stored directly in Trizol, ensure tissue is fully dissolved in Trizol, and spin down as directed.

3 Add  **300 µl**



Ethanol (95 - 100%), molecular grade

to each sample, mix well (invert tube, do not pipet), and transfer mixture to spin columns (with collection tube underneath). Spin samples at 13000 RPM at room temperature (  **21 °C** ) for  **00:00:30** and transfer column to new collection tube.


4 Add  **400 µl** Direct-zol RNA Prewash buffer to each column. Spin samples at 13000 RPM at room temp (  **21 °C** ) for  **00:00:30** and discard flowthrough. Repeat this step (2 total washes).

5 

Add **700 µl** RNA Wash buffer to each column. Spin at 13000 RPM at room temp ( **21 °C** ) for **00:02:00** .  
Spin one more time for **00:01:00** to ensure all liquid is removed from column.


6 

Add **8 µl**




**Molecular Biology Grade Water**  
by Fisher Scientific  
Catalog #: 10154604

directly to column matrix (be careful not to pierce column) and spin at 13000 RPM at room temp ( **21 °C** ) for **00:00:30** . IMMEDIATELY TRANSFER COLUMNS TO ICE. Quantify RNA concentration and purity (260/280) using




**Nanodrop 2000C**  
Thermo Scientific TSC-ND2000C




Make sure all RNA samples stay on ice and as sterile as possible to avoid RNase contamination. Move immediately to next step to limit RNA degradation.

7 On ice, in separate 200 uL PCR strip tubes, mix up to **1 µg** RNA (up to **6 µl** ) with **2 µl** d(T)23VN and



**Molecular Biology Grade Water**  
by Fisher Scientific  
Catalog #: 10154604

(if needed) to make an **8 µl** total reaction. Mix samples well and make sure nothing is sticking to the sides of the tube. Make 1 additional tube to serve as negative control.







Master mix of d(T)23VN and/or H2O should be made to avoid pipetting in and out of stock tube

8 Denature RNA at **70 °C** for **00:05:00** and cool down at **4 °C** for **00:05:00** . Use preset "cDNAincubate" function on




C1000 Touch Thermal Cycler  
PCR machine  
BioRad 1851148

. Immediately put on ice.

- 9 To each sample tube, add  **10 µl** M-MuLV Reaction Mix and  **2 µl** M-MuLV Enzyme Mix. To the negative control tube, add  **10 µl** M-MuLV Reaction Mix and  **2 µl**







Molecular Biology Grade Water  
by Fisher Scientific  
Catalog #: 10154604

. Total reaction volume for all samples should be  **20 µl** .




Ensure this step is done on ice

- 10 Incubate samples at  **42 °C** for  **01:00:00** and then inactivate enzyme at  **80 °C** for  **00:05:00** . Use preset "NEB cDNA" program on



C1000 Touch Thermal Cycler  
PCR machine  
BioRad 1851148



At this point, samples can either be used directly in qPCR assay or stored for later use at  **-20 °C**