



Nov 18, 2020

# Dip-C (Part 1: Chromosome Conformation Capture, for Fixed Nuclei)

Longzhi Tan<sup>1</sup><sup>1</sup>Stanford University

1 Works for me [dx.doi.org/10.17504/protocols.io.bpt7mnrn](https://dx.doi.org/10.17504/protocols.io.bpt7mnrn)



Longzhi Tan  
Stanford University

DOI

[dx.doi.org/10.17504/protocols.io.bpt7mnrn](https://dx.doi.org/10.17504/protocols.io.bpt7mnrn)

## PROTOCOL CITATION

Longzhi Tan 2020. Dip-C (Part 1: Chromosome Conformation Capture, for Fixed Nuclei). **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bpt7mnrn>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Nov 18, 2020

## LAST MODIFIED

Nov 18, 2020

## PROTOCOL INTEGER ID

44639

Digestion 40m

1 Thaw a tube of fixed nuclei **On ice**.

2 Prepare 0.5% SDS ( **50 µl** per sample; recipe below for **100 µl** ):

- **5 µl**

[Sodium dodecyl sulfate solution BioUltra for molecular biology 10% in H2O Sigma](#)

**Aldrich Catalog #71736**

(final: **0.5 Mass Percent** )


- **95 µl** water





- Vortex to mix


3 Resuspend nuclei in **50 µl** 0.5% SDS.


4 Incubate at **62 °C** for **00:10:00**.





10m

- 5 Add:
  -  **145 µl** water
 

 [Triton X-100, 10% solution](#)
  -  **25 µl** [Sigma Catalog #93443](#)
- 6 Rotate at  **37 °C** for  **00:15:00** . 15m
- 7 Add restriction enzyme and its buffer:
 

 [NEBuffer 2 \(10X\) New England](#)

  -  **25 µl** [Biolabs Catalog #B7002S](#)

 [Mbol \(25,000 units/ml\) - 2,500 units New England](#)
  -  **20 µl** [Biolabs Catalog #R0147M](#)
- 8 Rotate at  **37 °C**  **Overnight** . 15m

Ligation 40m

- 9 Centrifuge at  **1000 x g, 4°C, 00:05:00** .
- 10 Make Ligation Buffer (  **2 mL** per sample; recipe below for  **1 mL** ):
 

 [T4 DNA Ligase Reaction Buffer - 6.0 ml New England](#)



  -  **100 µl** [Biolabs Catalog #B0202S](#)

 [BSA, molecular biology grade, 20 mg/ml New England](#)
  -  **5 µl** [Biolabs Catalog # B9000S](#)
  -  **865 µl** water
  - Vortex to mix.
- 11 Remove supernatant leaving ~  **50 µl** .
- 12 Resuspend in  **1 mL** Ligation Buffer.
- 13 Centrifuge at  **1000 x g, 4°C, 00:05:00** .
- 14 Remove supernatant leaving ~  **50 µl** .



15 Resuspend in  **1 mL** [Ligation Buffer](#).

16  [T4 DNA Ligase \(1 U/μL\)](#) **Thermo**  
Add  **10 μl** [Fisher Catalog #15224025](#)

17 Pipette to mix.

18 Incubate at  **16 °C** for  **04:00:00**, occasionally inverting the tube.

4h

19  [Falcon 40 μm Cell](#)  
Optionally filter with [Strainer Corning Catalog #352340](#)  
 [Corning™ Falcon™ Test Tube with 35μm Cell Strainer Snap](#)  
[Cap Corning Catalog #352235](#)  
the flow cytometer.

or

to avoid clogging


20 Aliquot if needed.

21 Centrifuge at  **1000 x g, 4 °C, 00:05:00**.

22 Remove supernatant.



23 Store at  **-80 °C**.

Flow Sorting 40m

24 On the day of flow sorting, thaw a tube of ligated nuclei  **On ice**.

25  [PBS, pH 7.4](#) **Thermo**  
Resuspend in  **1 mL** [Fisher Catalog #10010023](#)

26 Make 300 μM DAPI:

-  [PBS, pH 7.4](#) **Thermo**  
 **100 μl** [Fisher Catalog #10010023](#)

- **2.1 µl** 14.3 mM (5 mg/mL) DAPI (stock made by dissolving **DAPI (4',6-Diamidino-2-Phenylindole, Dihydrochloride) Thermo Fisher Catalog #D1306**

in **2 mL** water

and stored at **4 °C** )

- Vortex to mix.

27 Add **1 µl** 300 uM DAPI (final: **0.5 Nanomolar (nM)** ).

28 Pipette to mix.

29 **96 well LoBind PCR plates Semi-**

Flow sort single nuclei into **skirted Eppendorf Catalog #0030129504**

either

dry or containing lysis buffer (which requires lysis by incubation before storing; see Part 2 for details).

30 Proceed directly to Part 2, or store at **-80 °C** .