


FEB 10, 2023

🌐 Preparation and transformation of chemically hyper-competent Escherichia coli

Andreas Sagen¹

¹University of Oslo



Andreas Sagen

University of Oslo, The National Institute of Occupational H...

ABSTRACT

Based on a protocol produced by Liu, et al.

Liu, Chang, W., Pan, L., Liu, X., Su, L., Zhang, W., Li, Q., & Zheng, Y. (2018). An Improved Method of Preparing High Efficiency Transformation Escherichia coli with Both Plasmids and Larger DNA Fragments. *Indian Journal of Microbiology*, 58(4), 448–456. <https://doi.org/10.1007/s12088-018-0743-z>

MATERIALS

LAF
Scale
Incubator
Heat cycler / Heat-bath
Ice

OPEN ACCESS

Protocol Citation: Andreas Sagen 2023. Preparation and transformation of chemically hyper-competent Escherichia coli. **protocols.io** <https://protocols.io/view/preparation-and-transformation-of-chemically-hyper-competent>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol

Created: Feb 03, 2023








Last Modified: Feb 10, 2023

PROTOCOL integer ID:
76363

Keywords: transformation, e. coli, Modified Inoue

Preparation of CRM transformation buffer

1 In a sterile flask, add  200 mL distilled water

- 2 Measure  2.5 mL ( 1 Molarity (M)) PIPES,  2.177 g Manganese(II) chloride tetrahydrate,  7.5 mL ( 1 Molarity (M)) Calcium chloride,  4.660 g Potassium chloride and  110 μ L (1 μ g/mL) Lactoferrin B.

Materials:


 PIPES **Sigma-aldrich Catalog #P1851**

 Manganese(II) chloride tetrahydrate **Sigma-aldrich Catalog #M3634** ^c


 Calcium chloride **Sigma-aldrich Catalog #C3881**

 Potassium chloride **Sigma-aldrich Catalog #P4504**

 Lactoferrin B **Sigma-aldrich Catalog #L9507**

- 3 Add measured reagents and mix for  00:05:00


5m

- 4 Adjust pH to  6.7 with Potassium hydroxide solution



Materials:

 Potassium hydroxide solution **Supelco Catalog #P4494**


- 5 Fill flask with distilled water to  250 mL









- 6 Filter sterilize solution with a filter (0.2 μ m) and store refrigerated ( 4 °C)

Preparation of chemically hyper-competent cells

- 7 Prepare a culture of *E. coli* on an LB agar plate. Pick a single colony, and inoculate in  500 mL S. O. C. broth in a  1000 mL flask.

Incubate at  18 °C with shaking at 100 rpm overnight, until OD₆₀₀ reaches 0.6


- 8 Aliquot entire culture volume into  50 mL canonical tubes

- 9 Place tubes on ice for  00:10:00 10m
- 10 Centrifuge tubes with  5000 rcf, 4°C, 00:10:00 10m
- 11 Discard supernatant and resuspend with  16 mL CRM transformation buffer
- 12 Incubate cells on ice for  00:10:00 10m
- 13 Centrifuge tubes with  5000 rcf, 4°C, 00:10:00 10m
- 14 Discard supernatant and resuspend with  8 mL CRM transformation buffer. Pool into two tubes
- 15 Centrifuge tubes with  5000 rcf, 4°C, 00:10:00 . Meanwhile, prepare  50 mL DMSO-CRM transformation buffer by diluting  3.5 mL DMSO in  46.5 mL CRM transformation buffer 10m



Materials:

 DMSO MP Biomedicals Catalog #196055

- 16 Discard supernatant and resuspend with  10 mL DMSO-CRM per tube

17 Incubate cells on ice for  00:30:00

30m

18 Aliquot  100 μL cell suspension into sterile  500 μL screw cap reaction tubes (Sarstedt #72.704.200).

Note

While creating aliquots, keep original tubes, and aliquots on ice, until snap-freeze take place

19 Snap freeze tubes in liquid nitrogen using a floating foam tube rack (Southern labware #HS2166)

20 Transfer aliquots storage box, and place in an ultra-low temperature freezer or vapor-phase nitrogen tank

Note

Store tubes in 50 mL canonical tubes, or similar containers

Transformation

1h 32m 30s

21 Quickly thaw a single reaction tube with  100 μL hyper-competent cells in CRM-DMSO









22 Mix 1-5 μL plasmid (ligation product)

Note

Do not exceed 5% of the volume competent cells

Note

Use up-to  25 ng per  50 μ L of competent cells

- 23** Incubate cells on ice for  00:30:00 30m
- 24** Heat-shock cells at  42 °C for  00:00:30, followed by  00:02:00 at  4 °C 2m 30s
immediatly
- 25** Add  500 μ L prewarmed S. O. C. medium and incubate for  37 °C at 200 rpm for  01:00:00 1h
- 26** Add desired amount of suspension on LB plates with ampicillin (100 μ g/mL) and incubate overnight