

Oct 11, 2022

# Tuning the expression levels of native genes

In 1 collection

Carolyn N Bayer<sup>1</sup>, Ana Gabriela Veiga Sepulchro<sup>1</sup>, Maja Rennig<sup>1</sup>, Morten Norholm<sup>1</sup>

<sup>1</sup>Technical University of Denmark

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.kqdg3px1ql25/v1](https://dx.doi.org/10.17504/protocols.io.kqdg3px1ql25/v1)

Carolyn N Bayer  
Technical University of Denmark

## ABSTRACT

This protocol collection describes how to use our optimised *tetA*<sup>OPT</sup> dual selection marker in *E. coli* K12 and Nissle. This dual selection marker can be used for positive selection based on tetracycline resistance and counterselection based on NiCl<sub>2</sub> sensitivity. *tetA* can be used to engineer all stages of the central dogma of molecular biology. On the DNA-level *tetA*<sup>OPT</sup> can be used to create scarless knockouts across the *E. coli* genome with an efficiency above 90%, whereas recombinant gene integrations can be achieved with approximately 50% efficiency. On the expression level, *tetA*<sup>OPT</sup> enables advanced genome engineering of both gene translation and transcription.

DOI

[dx.doi.org/10.17504/protocols.io.kqdg3px1ql25/v1](https://dx.doi.org/10.17504/protocols.io.kqdg3px1ql25/v1)

## PROTOCOL CITATION

Carolyn N Bayer, Ana Gabriela Veiga Sepulchro, Maja Rennig, Morten Norholm 2022. Tuning the expression levels of native genes. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.kqdg3px1ql25/v1>

## COLLECTIONS ⓘ

**tetA dual selection protocols**

## KEYWORDS

Recombineering, counterselection, selection markers, synthetic biology, genome engineering, tetA

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Mar 07, 2022

## LAST MODIFIED

Oct 11, 2022

## PROTOCOL INTEGER ID

59144

## PARENT PROTOCOLS

Part of collection

[tetA dual selection protocols](#)

## Ordering of oligonucleotides

### 1 order the following oligonucleotides:

**1.1** 2 primers annealing in the *tetA* cassette. These two primers have to each include 50 bp overhangs. The primer annealing upstream of *tetA* should contain homology to the 50 first basepairs of the gene of interest. The primer annealing downstream of *tetA* should contain homology to the promoter region of the gene of interest.

### 1.2 For introducing a TIR library while keeping the native promoter:

1 degenerated oligonucleotide. The primer needs to contain the randomization of the 6 nucleotides upstream of the start codon and also change the two codons downstream of the start to all synonymous codons. This primer needs to harbor 50 bp homology to the promoter region and 50 bp downstream of the randomized region.

#### For introducing a new promoter:

For short promoters, the promoter sequence can be included in the oligonucleotide that contains 50 bp homology to the regions up- and

downstream of *tetA*.

Promoter constructs over 100 bp need to be amplified via PCR. Create 2 primers, that bind in the promoter constructs and each contain 50 bp homology to the regions up- and downstream of *tetA*. Make sure to purify the PCR product over an agarose gel.

If an oligonucleotide is used for the removal of *tetA*, the primer needs to anneal to the lagging strand. Use [modest.biosustain.dtu.dk](https://modest.biosustain.dtu.dk) to create a sample "MAGE" oligonucleotide. Select your locus of interest in the dropdown menu "gene" and delete "A" in position 1. This will generate a MAGE oligonucleotide that will delete the A of the start codon in your gene that you want to delete. Align this oligonucleotide in your sequence software to see which strand it aligns to. This strand represents the lagging strand. Now create your custom oligonucleotide that aligns to the same strand as the MAGE sample oligonucleotide.

#### preculture and PCR - day 1

1d

- 2 Setup a preculture of the strain with pSIM19 (recombineering plasmid) in LB medium supplemented with Spectinomycin **10.05 mg/mL** and incubate at **250 rpm, 30°C** overnight. From now on the strain has to be kept at **30 °C** to maintain pSIM19 inside the cells.
- 3 Prepare a PCR product of the *tetA* cassette using a proof-reading polymerase and purify it.

#### Recombineering: tetA integration - day 2-4














3d










- 4 Prepare:

3d

4.1 Cold sterile water



4.2 Cold Glycerol **15 % volume**

- 4.3 Pre-chilled centrifuge and tabletop centrifuge at  **4 °C**
- 4.4 LB agar supplemented with  **0.05 mg/mL** tetracycline
- 4.5 M9 agar supplemented with  **50 micromolar (μM)** NiCl<sub>2</sub>
- 5 Inoculate  **50 mL** LB-Medium supplemented with Spectinomycin (  **0.05 mg/mL** ) with  **500 μL** of the preculture from step
  - 5.1 Incubate at  **250 rpm, 30°C** until cultures reached an OD600 of 0.5
- 6 Induce expression by transferring the culture to a shaking water bath at  **150 rpm, 42°C, 00:20:00**
- 7 Transfer culture to prechilled  **50 mL** falcon tubes and put on ice for  **00:15:00** 15m
- 8 Spin the culture down at  **4000 x g, 4°C, 00:15:00** and discard the supernatant 15m
- 9 Add  **1 mL** of ice cold water, resuspend and transfer to a  **1.5 mL** tube
- 10 Spin at  **11000 x g, 4°C, 00:00:30** in a tabletop centrifuge 30s

- 11 Wash pellet twice with  **1 mL** ice cold water 
  - 12 Resuspend the pellet in  **600 µL** cold glycerol ( **15 % volume** )
  - 13 Unused cells can be stored at  **-80 °C**
- This is not possible for *E. coli* Nissle
- 14 Electroporate  **50 µL** of cells with  **200 ng** of purified PCR product from step 3
  - 15 Recover cells  **800 rpm, 30°C, 01:00:00** in a tabletop shaker using SOC medium.
  - 16 plate cells on LB agar supplemented with  **0.05 mg/mL** tetracycline. Cell might need up to 2 days to grow.





Recombineering: *tetA* removal - day 4



1d

- 17 Select a colony from the LB tetracycline plate and start a preculture in LB medium supplemented with Spectinomycin  **0.05 mg/mL** . Incubate at  **250 rpm, 30°C** overnight.

By integrating *tetA* the expression of the gene of interest was interrupted. Possible changes in the phenotype can be screened for in this step.

- 18 prepare cells following steps 5-13

19 Electroporate  **50  $\mu$ L** of the prepared cells with either  **2  $\mu$ L** of a  **100 micromolar ( $\mu$ M)** oligonucleotide or  **200 ng** of a gel-purified PCR products.



20 Recover cells for  **800 rpm, 30°C, 01:00:00 in a tabletop shaker** . Afterwards transfer the cells into  **5 mL** LB medium supplemented with Spectinomycin

21 Incubate at  **250 rpm, 30°C** overnight



Important! Cells need to lose tetA transporter in the membrane to get resistant to  $\text{NiCl}_2$

Plating - day 5-7

2d

22 Wash  **1 mL** of the recovered cells twice with sterile water. Centrifuge at  **11000 rpm, 20°C, 00:00:30**


30s

23 Make a dilution series and plate  **100  $\mu$ L** of the 1:10 - 1:1000 dilution on M9 agar supplemented with  **50 micromolar ( $\mu$ M)**  $\text{NiCl}_2$

If the promoter was changed, the inducer can be added to the M9 plates


24 incubate the plates at  **30 °C** for  **48:00:00** to  **72:00:00**

5d

Cells can be incubated at  **37 °C** if it is the last step of recombination and pSIM19 is no longer needed in the cells

Screening - day 8

25 Screen for positive colonies by colony PCR to identify the correct recombinants. Restreak



correct colony on LB agar. In case of the TIR library, select colonies and assay the effect on the gene of interest.