



Feb 05, 2021

Library construction of metabarcoding at DNBSEQ-G400 with MGIEasy universal DNA library prep Kit

In 2 collections

Xiaohuan Sun¹, Yuehua Hu², Zewei Song¹¹BGI-Shenzhen; ²CAS Key Laboratory of Tropical Forest Ecology

1 Works for me dx.doi.org/10.17504/protocols.io.bn95mh86

Xiaohuan Sun

ABSTRACT

Library construction steps following MGIEasy universal DNA library prep set user manual started from end repair and A-tailing to single strand circularization in chapter 3. 500 ng of gDNA is used in this Library Construction Protocol. The steps were described below.

<https://en.mgitech.cn/Uploads/Temp/picture/20191023/5db00efef2bd7.pdf>

DOI

dx.doi.org/10.17504/protocols.io.bn95mh86

PROTOCOL CITATION

Xiaohuan Sun, Yuehua Hu, Zewei Song 2021. Library construction of metabarcoding at DNBSEQ-G400 with MGIEasy universal DNA library prep Kit. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bn95mh86>

COLLECTIONS ⓘ



Protocols for Efficient and stable metabarcoding sequencing data using DNBSEQ-G400 sequencer validated by comprehensive community analyses



Protocols for Efficient and stable metabarcoding sequencing data using DNBSEQ-G400 sequencer validated by comprehensive community analyses

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 31, 2020

LAST MODIFIED

Feb 05, 2021

OWNERSHIP HISTORY

Oct 31, 2020 Hongling Zhou

Nov 02, 2020 Xiaohuan Sun

PROTOCOL INTEGER ID

44061

Protocols for Efficient and stable metabarcoding sequencing data using DNBSEQ-G400 sequencer validated by comprehensive community analyses

Protocols for Efficient and stable metabarcoding sequencing data using DNBSEQ-G400 sequencer validated by comprehensive community analyses

- 1 End-repair and A-tailing steps are performed with ERAT Enzyme Mix in MGIEasy Universal DNA Library Prep Set (Cat. No: 1000006985).
- 2 Adaptor ligation step following the instruction below by using 5 mL MGIEasy DNA Adaptors.
- 3 Cleanup of Adapter-ligated DNA by using DNA clean beads and freshly prepared 80% ethanol, final DNA was elute by 40 mL of TE buffer.
- 4 19 µL of purified Adapter-ligated DNA sample with 25 mL PCR enzyme mix and 6 mL PCR primer mix were used for PCR amplification. We amplified samples using the following cycling conditions: 98 °C for 00:03:00^{9m 5s}; 7 cycles of 98 °C for 00:00:20, 60 °C for 00:00:15, and 72 °C for 00:00:30; and then a final extension at 72 °C for 00:05:00.
- 5 PCR products are cleaned by DNA clean beads. 32 mL of TE buffer is used to elute the DNA.
- 6 Quantify the purified PCR products with dsDNA Fluorescence Assay Kits such as Qubit® dsDNA HS Assay Kit or Quant-iT™ PicoGreen® dsDNA Assay Kit. The required yield for PCR products is ≥ 1 pmol. For pooled sequencing, please follow instructions provided by MGIEasy DNA Adaptors User Manual. The total yield after pooling should be 1 pmol, with a total volume ≤ 48 µL.
- 7 PCR products are stored at 95 °C for 00:03:00 for denaturation. After the reaction is complete, immediately^{5m} place the PCR tube on ice for 00:02:00.
- 8 Single strand circularization is prepared with single strand circularization reaction mixture consisting of 11.6 µL^{30m} splint buffer and 0.5 µL DNA rapid ligase. The mixture is incubated at 37 °C for 00:30:00. The mixture should be immediately placed on ice after the completeness of the reaction.
- 9 The circularized products are therefore used for sequencing at DNBSEQ-G400.