

MAR 27, 2024

🌐 Basic pipet technique training

Ying-Yu Hu¹

¹Dalhousie University

Marine Microbial Macroecology Lab
Tech. support email: ruby.hu@dal.ca



Ying-Yu Hu
Dalhousie University

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.36wgqjo4yvk5/v1

Protocol Citation: Ying-Yu Hu 2024. Basic pipet technique training. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.36wgqjo4yvk5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jan 08, 2023

ABSTRACT

Here we adopt the molybdate colorimetric reaction from our Total particulate phosphorus (TPP) measurement as a training material for our new lab members to learn basic pipetting technique.

CITATION

YY Hu, AJ Irwin, ZV Finkel (2022). Improving quantification of particulate phosphorus. *Limnology and Oceanography: Methods*.

LINK

[10.1002/lom3.10517](https://doi.org/10.1002/lom3.10517)

CITATIONS

YY Hu, AJ Irwin, ZV Finkel. Improving quantification of particulate phosphorus [10.1002/lom3.10517](https://doi.org/10.1002/lom3.10517)

Last Modified: Mar 27, 2024


PROTOCOL integer ID: 74966

Keywords: orthophosphate

Funders Acknowledgement:

Simons Foundation
Grant ID: 549937
Simons Foundation
Grant ID: 723789

PROTOCOL MATERIALS


 Ammonium molybdate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #09878-100G**

Step 6

 Ascorbic acid **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A5960-100G**

Step 7

SAFETY WARNINGS

 Acidic reagents need to be neutralized before dumped into the sink.

Preparing standard working solutions

2h

1 Date:


Trainee:

2 Need the following materials:

Supply	Check
Microtube rack	
17 X 2 mL microtubes	
MilliQ water container	
10 uL pipet	
100 uL pipet	
1000 uL pipet	
10 uL pipet tip	
100 uL pipet tip	
1000 uL pipet tip	
Vortex mixer	
Sharpies	
Primary stock solution of KH2PO4	
MilliQ	

3 How to aspire and dispense with pipet

- (a) When aspirating solution, ensure the pipette to be held vertically.
- (b) When dispensing, ensure you hold the pipette at an angle (10-45°). Working to these angles ensures the desired liquid amount is drawn into the tip properly and that all of the liquid is fully dispensed without leaving any residue in the tip.



<https://www.americanlaboratory.com/914-Application-Notes/240482-Ten-Tips-for-Proper-Pipetting>

Note

Always pay attention to the pipet tip, check if there is air bubble in the tip or slightly leaking on the tip.

4 Standard working solutions

Concentration of KH2PO4 primary stock	Unit

KH2PO4	Primary (uL)	Check	MilliQ (uL)	Check
S1	0		500	
S2	2		498	
S3	5		495	
S4	10		490	
S5	25		475	
S6	50		450	
S7	75		425	
S8	100		400	

Prepare working reagents

2h

5 Need the following materials:

Supply	Check
50 mL Falcon tube (3)	
Transfer pipet (1)	
Scoopula (1)	
Anti-static weighting plate (2)	
MilliQ	
Kimwipe	
Chamber balance	
Graduated cylinder (1X25 mL)	
Goggles	

6 [M] 2.5 % ammonium molybdate reagent

⊗ Ammonium molybdate **Sigma Aldrich Catalog #09878-100G**

6.1 In an anti-static weight plate, weigh ~0.25 g ammonium molybdate, transfer into the labelled falcon tube.

Actual weight (g)

6.2 Measure 10 mL MilliQ (use transfer pipet for the final topping)

6.3 Rinse the residue on the weighing plate into the falcon tube with the MilliQ.

6.4 Transfer all remaining MilliQ into the falcon tube

7 [M] 10 % ascorbic acid reagent (avoid light exposure)

⊗ Ascorbic acid **Sigma Aldrich Catalog #A5960-100G**

7.1 In an anti-static weight plate, weigh ~1 g ascorbic acid, transfer into the labelled falcon tube.

Actual weight (g)

7.2 Measure 9 mL MilliQ (use transfer pipet for the final topping)

7.3 Rinse the residue on the weighing plate into the falcon tube with the MilliQ.


7.4 Transfer all remaining MilliQ into the falcon tube


8 Sulfuric acid reagent:

8.1 Measure 17 mL MilliQ (use transfer pipet for the final topping) and transfer into a Falcon tube

8.2 Carefully add  1 mL  18 M concentrated sulfuric acid into the water

9 Vortex ammonium molybdate and ascorbic acid solutions.


10 Use graduated cylinder, measure and transfer 6 mL  2.5 % ammonium molybdate reagent into the sulfuric acid tube.

11 Use graduated cylinder, measure and transfer 6 mL  10 % ascorbic acid reagent into the sulfuric acid tube.

12 Vortex the reagent.

Colorimetric measurement

2h

13 Preheat incubator/shaker to  37 °C

Equipment	
SHAKING INCUBATOR	NAME
71L	TYPE
Corning® LSE™	BRAND
6753	SKU

- 14Add

🧪 500 µL

 reagent to each standard working solution.
- 15Vortex each tube.
- 16Incubate at

🌡️ 37 °C

 for

🕒 01:00:00

 while shaking at 200 rpm

1h
- 17Load microplate with 250 uL reactant from each tube, duplicate.

18 Layout in the microplate:

	1	2	3	4
A	S1	S1	S1	S1
B	S2	S2	S2	S2
C	S3	S3	S3	S3
D	S4	S4	S4	S4
E	S5	S5	S5	S5
F	S6	S6	S6	S6
G	S7	S7	S7	S7
H	S8	S8	S8	S8

19 Turn on microplate reader

Equipment

Varioskan LUX Multimode Microplate Reader

NAME

Thermo Fisher

BRAND

VL0L00D0

SKU

20 Setup layout and method

A	B
Shake duration	00:00:05
Shaking type	Continuous
Shaking force	High
Shaking speed [rpm]	600
Wavelength [nm]	820
Measurement Time [ms]	100

21 Place microplate onto the plate holder

22 Read the plate

23 Export data to Excel sheet

Calculation

24 In the Excel sheet, subtract the average absorbance at 820 nm of the blank standard replicates from the absorbance at 820 nm of all other standard working solutions.

This gives the corrected absorbance from each well.

25 Calculate the concentration of each standard working solution

Unit for each content:

Primary concentration g/L

Primary volume uL

MilliQ volume uL

Standard working solution g/L

Phosphate in standard working solution uM

26 Y-axis: Corrected absorbance

X-axis: Phosphate (uM)

Calculate slope (S), intercept (b), coefficient of determination (R^2), standard deviation (S_y) of the absorbance from the four replicates of the blank, and $LOD = 3.3 \times S_y / S$

Parameter	Value
S	
b	
R2	
Sy	

Parameter	Value
LOD	

Review

- 27 What information does LOD provide?
- 28 Is it necessary to be aware of the order when mixing water with concentrated sulfuric acid?

Clean-up

- 29 Neutralize acidic chemicals with baking soda in a plastic beaker
- 30 Rinse microtubes and microplate with tap water
- 31 Dispose microtubes and microplate into recyclable garbage bin
- 32 Dispose pipet tips into sharp garbage container
- 33 Rinse graduated cylinder, Falcon tubes and MilliQ container with DI water, place them into a basket to air-dry

