

MAR 01, 2024

🌐 DNA Extraction with ZymoBIOMICS MiniPrep Kit

Carlos Goller¹

¹North Carolina State University

BIT Metagenomics

MICROnet



Carlos Goller

ABSTRACT

DNA extraction of bacterial DNA Using ZymoBIOMICSTM DNA Miniprep Kit (D4300) using lysis tubes with DNA/RNA Shield.

IMAGE ATTRIBUTION

Bacteria from Dmitry Mirolyubov from The Noun Project

MATERIALS

ZymoBIOMICSTM DNA Miniprep Kit

PROTOCOL MATERIALS



DNA/RNA Shield **Zymo Research Catalog #R1100-50**

Step 1

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.81wgbxjn1lpk/v1

Protocol Citation: Carlos Goller 2024. DNA Extraction with ZymoBIOMICS MiniPrep Kit. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.81wgbxjn1lpk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Feb 27, 2024

Last Modified: Mar 01, 2024




PROTOCOL integer ID: 95833

Keywords: DNA, bacteria, isolate, microbe

Funders Acknowledgement:
Biotechnology Program
Grant ID: NC State University

D


14m











- 1
- Add  Sample to a ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm). Add  750 µL ZymoBIOMICSTM 10m
Lysis Solution to the tube and cap tightly. Note: For samples stored and lysed in  DNA/RNA Shield **Zymo Research Catalog #R1100-50** Lysis Tubes, do not add ZymoBIOMICSTM Lysis Solution and proceed to Step 2.








A	B
Sample Type	Maximum Input
Feces	200 mg
Soil	250 mg
Liquid Samples and Swab Collections	250 ul
Cells (isotonic buffer, el.g., PBS)	50-100 mg (wet weight) (109 bacterial and 108 yeast cells)
Samples in DNA/RNA Shield	<1 ml

Table from ZymoBIOMICS Protocol D4300

- 2
- Obtain DNA/RNA Shield Lysis Tubes with samples/swab heads from instructors. Secure in a bead beater 40m
fitted with a 2 ml tube holder assembly and process using optimized beat-beating conditions (speed and

time) for different devices. We will use a Vortex Genie with 2ml Bashing Bead tubes for  00:40:00 .

- 3 Centrifuge the ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm) in a microcentrifuge at ≥  10000 rpm, 20°C, 00:01:00 . 1m
- 4 Transfer up to  400 µL supernatant to the Zymo-Spin™ III-F Filter in a Collection Tube and centrifuge  8.000 x g, 00:01:00 . Discard the Zymo-Spin™ III-F Filter. 1m
- 5 Add  1200 µL of ZymoBIOMICS™ DNA Binding Buffer to the filtrate in the Collection Tube from Step 4. Mix well.
- 6 Transfer  800 µL of the mixture from Step 5 to a Zymo-Spin™ IICR Column in a Collection Tube and centrifuge at  10.000 x g, 00:01:00 . 1m
- 7 Discard the flow through from the Collection Tube and repeat Step 6.
- 8 Add  400 µL ZymoBIOMICS™ DNA Wash Buffer 1 to the ZymoSpin™ IICR Column in a **new** Collection Tube and centrifuge at  10.000 x g, 00:01:00 . Discard the flow-through. 1m
- 9 Add  700 µL ZymoBIOMICS™ DNA Wash Buffer 2 to the ZymoSpin™ IICR Column in a Collection Tube and centrifuge at  10.000 x g, 00:01:00 . Discard the flow-through. 1m

- 10** Add  200 µL ZymoBIOMICSTM DNA Wash Buffer 2 to the ZymoSpinTM IICR Column in a Collection Tube and centrifuge at  10.000 x g, 00:01:00 . 1m
- 11** Transfer the Zymo-SpinTM IICR Column to a clean 1.5 ml microcentrifuge tube and add 100 µl (50 µl minimum) ZymoBIOMICSTM DNase/RNase Free Water directly to the column matrix and incubate for  00:01:00 . Centrifuge at  10.000 x g, 00:01:00 to elute the DNA5, 6. 2m
- 12** Place a Zymo-SpinTM III-HRC Filter in a new Collection Tube and add  600 µL ZymoBIOMICSTM HRC Prep Solution. Centrifuge at  8.000 x g, 00:03:00 . 3m
- 13** Transfer the eluted DNA (Step 11) to a prepared Zymo-SpinTM III-HRC Filter in a clean 1.5 ml microcentrifuge tube and centrifuge at exactly  16.000 x g, 00:03:00 . The filtered DNA is now suitable for PCR and other downstream applications. 3m