



APR 17, 2023

🌐 Lakes ABPS Protocol - Optimized protocol for the extraction of fish DNA from freshwater sediments (Thomson-Laing et al., 2022)

Georgia Thomson-Laing¹

¹Cawthron Institute, Nelson, New Zealand



Grayson Huston

ABSTRACT

DNA was extracted from lake sediment samples by an alkaline lysis method with ethanol precipitation adapted from method described by Kuwae et al. (2020); Sakata et al. (2020b); Sakata et al. (2020a).

In a comparison of multiple sedDNA extraction methods, the ABPS (Alkaline buffer - power soil) protocol yielded the highest concentrations of target genets across a range of lake sediments. This protocol was further optimized (65C incubation temperature, pooling of multiple PowerSoil extractions) to overcome technical challenges related to co-precipitation of organic content in lake-surface sediments.

The optimized ABPS protocol is called the "Lakes ABPS protocol"

This protocol has proven to be successful at detecting fish sedDNA from surface sediments in multiple systems for multiple species.

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.yxmvm24d6g3p/v1

Protocol Citation: Georgia Thomson-Laing 2023. Lakes ABPS Protocol - Optimized protocol for the extraction of fish DNA from freshwater sediments (Thomson-Laing et al., 2022). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvm24d6g3p/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working




Created: Apr 17, 2023


Last Modified: Apr 17, 2023

PROTOCOL integer ID:
80636



Alkaline extraction


1h 51m

- 1 **INTO** a 50 mL tube add:
-  10 g of sediment sample
 -  6 mL sodium hydroxide (0.33M)
 -  3 mL Tris-EDTA (pH 8)

2 **VORTEX** for  00:01:00

51m

INCUBATE at  65 °C for  00:50:00


3 **ALLOW** samples to cool to  Room temperature


1h

CENTRIFUGE at  15000 x g for  01:00:00

Ethanol precipitation

2h

4 **TRANSFER**  7.5 mL of supernatant to a new 50 mL tube


ADD  7.5 mL of Tris HCl (1M, pH 6.7) to neutralize

5 **ADD**  1.5 mL sodium acetate (3M, pH 5.2)

1h

ADD  30 mL of molecular grade 100% ethanol

INCUBATE samples at  -20 °C for  01:00:00

6 **CENTRIFUGE** samples at  10000 x g for  01:00:00

1h


DISCARD supernatant

RETAIN precipitated pellet

DNeasy PowerSoil extraction

7 **EXTRACT** the total pellet using multiple DNeasy PowerSoil DNA Isolation Kit extractions following the manufacturer's instructions

Note

 0.25-0.5 g of pellet per extraction

POOL resultant DNA elutes

DNA is now ready for downstream applications