

Jul 31, 2020

DNeasy Protocol

 [Nature Communications](#)

DOI

dx.doi.org/10.17504/protocols.io.bi8xkhxn

John E Gorzynski¹

¹Stanford University



John E Gorzynski

Stanford University

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bi8xkhxn>

External link: <https://doi.org/10.1038/s41467-022-32397-8>

Document Citation: John E Gorzynski 2020. DNeasy Protocol. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bi8xkhxn>

Manuscript citation:

Parikh, V.N., Ioannidis, A.G., Jimenez-Morales, D. *et al.* Deconvoluting complex correlates of COVID-19 severity with a multi-omic pandemic tracking strategy. *Nat Commun* **13**, 5107 (2022). <https://doi.org/10.1038/s41467-022-32397-8>

License: This is an open access document distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: July 31, 2020

Last Modified: July 31, 2020

Document Integer ID: 39927

Keywords: protocol

Troubleshooting



DNeasy protocol-- nasal swabs in VTM and DNA RNA Shield

Date _____

Preheat Nuclease Free water to 95C for elution

4X sets of labeled tubes (2X from the 2ml tubes provided in the kit and an additional 2X 1.5ml eppendorf tubes)

1. To **200ul** of sample add **200ul**ATL Buffer and Vortex
2. Incubate 56C **10** min
3. Add **200ul**100% EtOH
4. Pipette mixture into a DNeasy spin column and place in a 2ml collection tube
5. Centrifuge **>6,000** rcf (g) for **1** min at RT
6. Discard collection tube and replace with a fresh collection tube
7. Add **500ul**AW1 into spin column
8. Centrifuge **>60,000** rcf (g) for **1** min at RT
9. Discard collection tube and replace with a fresh collection tube
10. Add **500ul**AW2 into spin column
11. Centrifuge **20,000** rcf for **1** min
12. Discard collection tube and replace with a clean 1.5ml low-bind eppendorf tube
13. Add **30ul**Nuclease Free Water into spin column
14. Incubate at RT for 1 min
15. Centrifuge **>6,000** rcf for **1** min
16. Quantify using Qubit and Nanodrop
17. Store at -20C