

Jul 11, 2018

Testing antibacterial activity of bacteriophages on hard surfaces

DOI

<https://dx.doi.org/10.17504/protocols.io.nwmdfc6>

Elisabetta Caselli¹

¹University of Ferrara



Elisabetta Caselli

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.nwmdfc6>

Protocol Citation: Elisabetta Caselli 2018. Testing antibacterial activity of bacteriophages on hard surfaces. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.nwmdfc6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 19, 2018

Last Modified: May 18, 2026

Protocol Integer ID: 10925

Keywords: bacteriophages, bacterial contamination, antibacterial activity of bacteriophage, testing antibacterial activity, bacteriophage, hard surface

Abstract

Testing antibacterial activity of bacteriophages on hard surfaces

Attachments



Phages on surfaces.d...

14KB

Guidelines

To test the killing/decontaminating activity of bacteriophages against specific bacterial targets.

Materials

MATERIALS

 Luria-Bertani (LB) broth, makes 1L **Amresco Catalog #K488**

Safety warnings

-  Perform experiments in sterile conditions under a biohazard flow chamber.
Biosafety level=2

Before start

The following materials are needed:

- individual pure bacterial cultures
- bacteriophage preparation (either commercial or home-prepared)
- bacterial growth media (liquid and agarized)
- Rodac contact agar plates

