

Sep 19, 2019

Unsucessful atemps on transfromation of the Ichthyosporean *Sphaeroforma artica*

DOI

<https://dx.doi.org/10.17504/protocols.io.z6ef9be>

Elena Casacuberta¹, Cristina Aresté², Omayra Dudin¹

¹Institute of Evolutionary Biology; ²IIBB

Protist Research to Opti...

Multicellgenomelab

1 more workspace



Elena Casacuberta

Institute of Evolutionary Biology

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.z6ef9be>

Protocol Citation: Elena Casacuberta, Cristina Aresté, Omayra Dudin 2019. Unsucessful atemps on transfromation of the Ichthyosporean *Sphaeroforma artica*. protocols.io <https://dx.doi.org/10.17504/protocols.io.z6ef9be>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 18, 2019

Last Modified: May 18, 2026

Protocol Integer ID: 22438

Keywords: electroporation on spaheroforma artica, testing electroporation, ichthyosporean sphaeroforma artica, different electroporators neon, spaheroforma artica, transfromation, unsucessful atemps on transfromation



Abstract

Testing electroporation on *Sphaeroforma artica* with two different electroporators

Neon (invitrogen)

Lonza (nucleofector)

Materials

MATERIALS

⊗ 100ml Sorbitol solution [1M] **G-Biosciences Catalog #R038**

⊗ Glycerol **Merck Millipore (EMD Millipore) Catalog #4095**

⊗ Marine Broth 2216 **Becton Dickinson (BD) Catalog #279110**

⊗ Neon™ Transfection System 10 µL Kit **Thermo Fisher Scientific Catalog #MPK1096**

⊗ Lipofectamine 2000 **Life Technologies**

Sphaeroforma artica cells



- 1 Cells were seeded 24 hours prior transfection
- 2 200.000 cells per transfection condition were tested
- 3 Tested techniques and settings are as follows (with 1 to 5 micrograms of DNA endogenous promoter (actin) plus fluorescent reporter (venus or mCherry))

	<i>S.arctica</i>	Neon	1000-2500V, 10-40 ms, 1-3 pulses	n/a (not successful)
		Lipofectamina		n/a (not successful)
		LONZA®	16 preset codes P3/P4/P5 buffer	n/a (not successful)

- 4 cells were recovered immediately in MB medium and inspected for positive transfection by fluorescence after 24, 48 and up to 5 days.