

Nov 06, 2017

Version 2

DNA extraction of Nephila clavipes using the Wizard Genomic DNA Purification kit (Promega) V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.knxcvfn>



Luiz Filipe Bartoleti¹

¹University of Campinas



Luiz Filipe Bartoleti

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.knxcvfn>

Protocol Citation: Luiz Filipe Bartoleti 2017. DNA extraction of Nephila clavipes using the Wizard Genomic DNA Purification kit (Promega) . protocols.io <https://dx.doi.org/10.17504/protocols.io.knxcvfn>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: November 06, 2017

Last Modified: May 18, 2026

Protocol Integer ID: 8631

Keywords: dna extraction, spider, kit, araneae, Nephila clavipes, dna extraction of nephila clavipe, whole genome dna from nephila clavipe, wizard genomic dna purification kit, nephila clavipe, araneidae, whole genome dna, dna, extraction, purification, genome

Abstract

Protocol used to extract whole genome DNA from Nephila clavipes (Araneae: Araneidae) specimens.

Attachments



Protocol for DNA ext...

13KB

- 1 Macerate 1-8 legs (around 50 mg of tissue) using liquid nitrogen.
- 2 Add 600 μ L of Nuclei Lysis solution and homogenize.
 600 μ L Nuclei Lysis solution
- 3 Take the microtubes to the dry bath at 65°C for one hour.
 01:00:00 Dry bath
 65 °C
- 4 Add 200 μ L of Protein Precipitation solution and homogenize.
 200 μ L Protein Precipitation solution
- 5 Centrifuge the microtubes at 13000 rpm for 10 minutes.
 00:10:00 Centrifugation
- 6 Transfer the supernatant to a new microtube and add 600 μ L of Isopropyl alcohol.
 600 μ L Isopropyl alcohol
- 7 Keep the tubes at -20°C for at least one hour.
 -20 °C
 01:00:00 Keeping tubes at -20°C
- 8 Centrifuge the microtubes at 13000 rpm for 6 minutes at 4°C.
 4 °C
 00:06:00 Centrifugation
- 9 Discard the supernatant and add 600 μ L of absolut etanol.
 600 μ L Etanol
- 10 Homogenize gently for one minute.
 00:01:00 Homogenizing
- 11 Centrifuge the microtubes at 13000 rpm for 3 minutes at 4°C.
 4 °C
 00:03:00 Centrifugation
- 12 Discard the supernatant.



13 Let the samples dry.

14 Rehydrate with 50 μ L of Tris-EDTA buffer.

 50 μ L Tris-EDTA buffer