



Jun 24, 2022

Determination of MIC according to the EUCAST E.DEF 11.0 method

Khalid El Moussaoui¹¹Université de Liège[Other](#)[Share](#)dx.doi.org/10.17504/protocols.io.kxygxz83kv8j/v1Khalid El Moussaoui
Université de Liège

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

Determination of MIC according to the EUCAST E.DEF 11.0 method.

DOI

dx.doi.org/10.17504/protocols.io.kxygxz83kv8j/v1

PROTOCOL CITATION

Khalid El Moussaoui 2022. Determination of MIC according to the EUCAST E.DEF 11.0 method. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kxygxz83kv8j/v1>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 24, 2022

LAST MODIFIED

Jun 24, 2022

PROTOCOL INTEGER ID

65232

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

- 1 With the disinfected tweezers, remove the filter (attention: do not take the filters of protections having a slightly bluer tint) and place them in the support provided for this purpose. Close the holder.
- 2 Pour the contents of TUBE A (sterile water + Tween) onto the Petri dish containing the starting sample and wipe with a swab until the suspension is relatively cloudy.
- 3 Aspirate the suspension with the syringe.
- 4 Filter this suspension by injecting it into the filter (0,45 µm) holder and collect the filtered suspension in suspension in TUBE A.
- 5 Measure the density of TUBE A and adjust it to 0.5 McFarland. Then take 500 µL of this suspension and place it in TUBE B (to be done 2x since we are working in duplicates). TUBE B contains sterile demineralized water (Thermo Fisher T3339).

- 6 Place an inoculation head on the tube and inoculate with the Sensititre Aim.
- 7 Incubate at room temperature for 7 days and read the result at 492 nm.