


JUN 22, 2023

TIMSTOF Compass 4.0 Methods - LCMS and MALDI - OrsburnLab 2023

Benjamin Orsburn¹

¹The Johns Hopkins University School of Medicine

LCMSMethods.org Version 2



lcmsmethods

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.x54v9p8zzg3e/v1

Protocol Citation: Benjamin Orsburn 2023. TIMSTOF Compass 4.0 Methods - LCMS and MALDI - OrsburnLab 2023. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.x54v9p8zzg3e/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jun 22, 2023

Last Modified: Jun 22, 2023

PROTOCOL integer ID:
83870

DISCLAIMER

No guarantees are implicit or implied in the use of these methods. For research use only.

ABSTRACT

The methods included here are in support of research submissions from the www.OrsburnLab.org group in 2023 from a TIMSTOF Flex system with EasyNLC1200, EvoSep One and MALDI (1) front end sources.

Methods include:

pasefRIQ used for SCoPE-MS (DefaultTIMSTOF40_TMT_highsensitivity_pasefSCOPE diaPASEF methods optimized for the EvoSep One - diaPASEF-shortgradient-1800V the same - optimized for EasyNLC 1200 using a 10um emitter - TIMSTOF_40_shortgradient_DIA ddaPASEF method for normal high input (100ng - 500ng proteomics loads) - DefaultDDA_TIMST40_5_6_2023 ddaPASEF methods optimized for the EvoSep One Whisper for ultralow input immunoenriched samples - DefaultDDA_TIMS40_high_sensitivity

A MALDI-MS method for TIMSTOF Flex to test the effectiveness of direct transfer of MALDI instrument methods between different TIMSTOF imaging mass spectrometers:
20220925_SCZ_GelMethod_AngelLab_transfer.m

ATTACHMENTS

[TIMSTOF40_methods.zip](#)

1