



MAR 19, 2024

Neuronal co-culture

Nisha Mohd Rafiq¹, Pietro De Camilli^{2,3,4,5}

¹University of Tuebingen;

²Department of Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA.;

³Department of Cell biology, Yale University School of Medicine, New Haven, Connecticut 06510, USA.;

⁴Program in Cellular Neuroscience, Neurodegeneration and Repair.;

⁵Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA.

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.dm6gpze38lzp/v1

Protocol Citation: Nisha Mohd Rafiq, Pietro De Camilli 2024. Neuronal co-culture. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gpze38lzp/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 12, 2024

ASAP Collaborative Research Network



Berrak Ugur
Yale University

ABSTRACT

This protocol describes the co-culturing of iPSC-derived dopaminergic (DA) neurons and iPSC-derived medium spiny neurons (MSNs) in a microfluidic compartmentalization device.

ATTACHMENTS

[Rafiq_co-culture.docx](#)

Last Modified: Mar 19, 2024

PROTOCOL integer ID: 96913

Keywords: ASAPCRN

Funders Acknowledgement:
ASAP
Grant ID: ASAP-000580

MATERIALS

Preparation of device for seeding cells:

Prepare enough amount of NB/B27 medium.

1. For 500 mL

- NB/B27 medium (see Method section of paper), add:

A	B
Neurobasal medium	484 mL
B27 supplement without vitamin A	10 mL
GlutaMAX	5 mL
Penicillin-Streptomycin	1 mL

Storage: NB/B27 medium can be stored for 5 days at 4 °C or for up to one month at -20 °C .

- Warm NB/B27 medium at 37 °C .
- Make NB/B27 complete medium by adding:



A	B
BDNF	20 ng/ml
Ascorbic acid	0.2 mM
GDNF	20 ng/ml
db-cAMP	0.5 mM
TGFβ3	1 ng/ml
DAPT	10 uM
Y-27632	10 uM

Neuronal co-culture device set-up


1

Note

The OMEGA4 device has 2 pairs of interconnected chambers, where each pair of chambers is joined via a series of microfluidic channels.

Coat chambers with  200 µL per well with  0.1 undetermined Poly-L-Ornithine (PLO) in PBS.

2

Incubate plates overnight at  37 °C .






3


Wash the chambers thrice with PBS.



4

Coat chambers with  200 µL per well  10 undetermined Laminin plus  2 undetermined fibronectin, both diluted in PBS.

5

Incubate plates overnight at  37 °C . Do not store coated plates. Proceed with preparation of plates for seeding cells.



Preparation of device for seeding cells

15m

6

Prepare enough amount of NB/B27 medium.

1. For  500 mL

- NB/B27 medium (see Method section of paper), add:

A	B
Neurobasal medium	484 mL
B27 supplement without vitamin A	10 mL
GlutaMAX	5 mL
Penicillin-Streptomycin	1 mL

Storage: NB/B27 medium can be stored for 5 days at 4 °C or for up to one month at -20 °C .

- Warm NB/B27 medium at 37 °C .
- Make NB/B27 complete medium by adding:

A	B
BDNF	20 ng/ml
Ascorbic acid	0.2 mM
GDNF	20 ng/ml
db-cAMP	0.5 mM
TGFβ3	1 ng/ml
DAPT	10 uM
Y-27632	10 uM

7 Discard coating reagents and add 200 µL per well of NB/B27 complete medium. 15m

Keep the plate at 37 °C for 00:15:00 before seeding cells.

8 Replate cultured iPSC-derived dopaminergic neurons (day 30, see Method section of paper) on one side of the two-chamber microfluidic compartmentalization device (OMEGA4, eNuvio) at a cell concentration of 3×10^5 . Only the axons of DA neurons can migrate through the microfluidic channels connected to the adjacent chamber.

9 Feed neurons with fresh NB/B27 media every 3 days. Add 10 undetermined Laminin to NB/B27 media every 10 days before feeding the neurons.



- 10 After an additional 25 days in the co-culture device, thaw frozen iPSC-derived medium spiny neurons (MSN) from BrainXell and plate on the other half of the device (where only the axons of DA neurons are present) at a cell concentration of 3×10^5 cells.
- 11 Fix the DA-MSN co-cultures till 7-10 days later for immunofluorescence (see Method section of paper).