



Version 4

Nov 15, 2020

# NEBNext Ultra End Prep Mixture E7442 V.4

Isabel Gautreau<sup>1</sup><sup>1</sup>New England Biolabs

1

Works for me

[dx.doi.org/10.17504/protocols.io.bpiumkew](https://dx.doi.org/10.17504/protocols.io.bpiumkew)**New England Biolabs (NEB)**Tech. support phone: +1(800)632-7799 email: [info@neb.com](mailto:info@neb.com)Isabel Gautreau  
New England Biolabs

## EXTERNAL LINK

<https://www.neb.com/protocols/2012/11/15/nebnext-end-prep-e7370>

## DOI

[dx.doi.org/10.17504/protocols.io.bpiumkew](https://dx.doi.org/10.17504/protocols.io.bpiumkew)

## EXTERNAL LINK

<https://www.neb.com/protocols/2012/11/15/nebnext-end-prep-e7370>

## PROTOCOL CITATION

Isabel Gautreau 2020. NEBNext Ultra End Prep Mixture E7442. **protocols.io**<https://dx.doi.org/10.17504/protocols.io.bpiumkew>

Version created by Isabel Gautreau

## EXTERNAL LINK

<https://www.neb.com/protocols/2012/11/15/nebnext-end-prep-e7370>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Nov 09, 2020

## LAST MODIFIED

Nov 15, 2020

## OWNERSHIP HISTORY

Nov 09, 2020



Lenny Teytelman

protocols.io

Nov 09, 2020



Isabel Gautreau

New England Biolabs

## PROTOCOL INTEGER ID

44340

## MATERIALS TEXT

## MATERIALS

[NEBNext End Prep Enzyme Mix](#) **New England****Biolabs Catalog #E7646**

Step 1

[NEBNext End Repair Reaction Buffer New England](#)

Biolabs Catalog #E7372

Step 1

STEP MATERIALS

[NEBNext End Prep Enzyme Mix New England](#)

Biolabs Catalog #E7646

Step 1

[NEBNext End Repair Reaction Buffer New England](#)

Biolabs Catalog #E7372

Step 1

[NEBNext End Prep Enzyme Mix New England](#)

Biolabs Catalog #E7646

Step 1

[NEBNext End Repair Reaction Buffer New England](#)

Biolabs Catalog #E7372

Step 1

## NEBNext End Prep

1

**Starting Material:** 5 ng-1 µg fragmented DNA

Mix the following components in a sterile, nuclease-free tube:

Component	Volume
End Prep Enzyme Mix	3 µl
End Repair Reaction Buffer (10X)	6.5 µl
Fragmented DNA	55.5 µl
Total Volume	65 µl

[NEBNext End Prep Enzyme Mix New England](#)

Biolabs Catalog #E7646

[NEBNext End Repair Reaction Buffer New England](#)

Biolabs Catalog #E7372

2 Mix by pipetting, followed by a quick spin to collect all liquid from the sides of the tube.

3 Place in a thermocycler, with the heated lid on, and run the following program:

**30 minutes @ 20°C**

**30 minutes @65°C**

**Hold at 4°C**

4 Proceed directly to NEBNext Ultra Ligation Module (NEB #E7445).