



Sep 16, 2020

Protocols for Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons

Yingchao Xue^{1,2}, Xiping Zhan³, Shisheng Sun⁴, Senthilkumar S. Karuppagounder^{5,6,7}, Shuli Xia^{2,5}, Valina L Dawson^{5,6,7,8,9}, Ted M Dawson^{5,6,7,8,10}, John Laterra^{2,5,8,11}, Jianmin Zhang¹, Mingyao Ying^{2,5}

¹Department of Immunology, Research Center on Pediatric Development and Diseases, Institute of Basic Medical Sciences, Chinese Academy of Medical Sciences and School of Basic Medicine, Peking Union Medical College, State Key Laboratory of Medical Molecular Biology;

²Hugo W. Moser Research Institute at Kennedy Krieger; ³Department of Physiology and Biophysics, Howard University;

⁴College of Life Sciences, Northwest University; ⁵Department of Neurology, Johns Hopkins University School of Medicine;

⁶Neuroregeneration and Stem Cell Programs, Institute for Cell Engineering, Johns Hopkins University School of Medicine;

⁷Adrienne Helis Malvin Medical Research Foundation; ⁸Department of Neuroscience, Johns Hopkins University School of Medicine;

⁹Department of Physiology, Johns Hopkins University School of Medicine;

¹⁰Department of Pharmacology and Molecular Sciences, Johns Hopkins University School of Medicine;

¹¹Department of Oncology, Johns Hopkins University School of Medicine

1 Works for me This collection is published without a DOI.

Neurodegeneration Method Development Community
Tech. support email: ndcn-help@chanzuckerberg.com



Anita Broellochs
protocols.io

ABSTRACT

Proneural transcription factors (TFs) drive highly efficient differentiation of pluripotent stem cells to lineage-specific neurons. However, current strategies mainly rely on genome-integrating viruses. Here, we used synthetic mRNAs coding two proneural TFs (Atoh1 and Ngn2) to differentiate induced pluripotent stem cells (iPSCs) into midbrain dopaminergic (mDA) neurons. mRNAs coding Atoh1 and Ngn2 with defined phosphosite modifications led to higher and more stable protein expression, and induced more efficient neuron conversion, as compared to mRNAs coding wild-type proteins. Using these two modified mRNAs with morphogens, we established a 5-day protocol that can rapidly generate mDA neurons with >90% purity from normal and Parkinson's disease iPSCs. After in vitro maturation, these mRNA-induced mDA (miDA) neurons recapitulate key biochemical and electrophysiological features of primary mDA neurons and can provide high-content neuron cultures for drug discovery. Proteomic analysis of Atoh1-binding proteins identified the nonmuscle myosin II (NM-II) complex as a new binding partner of nuclear Atoh1. The NM-II complex, commonly known as an ATP-dependent molecular motor, binds more strongly to phosphosite-modified Atoh1 than the wild type. Blebbistatin, an NM-II complex antagonist, and bradykinin, an NM-II complex agonist, inhibited and promoted, respectively, the transcriptional activity of Atoh1 and the efficiency of miDA neuron generation. These findings established the first mRNA-driven strategy for efficient iPSC differentiation to mDA neurons. We further identified the NM-II complex as a positive modulator of Atoh1-driven neuron differentiation. The methodology described here will facilitate the development of mRNA-driven differentiation strategies for generating iPSC-derived progenies widely applicable to disease modeling and cell replacement therapy.

EXTERNAL LINK

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6344911/>

THIS COLLECTION ACCOMPANIES THE FOLLOWING PUBLICATION

Xue Y, Zhan X, Sun S, et al. Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons. *Stem Cells Transl Med.* 2019;8(2):112–123. doi:10.1002/sctm.18-0036

ATTACHMENTS

[Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons.pdf](#)

COLLECTION CITATION

Yingchao Xue, Xiping Zhan, Shisheng Sun, Senthilkumar S. Karuppagounder, Shuli Xia, Valina L Dawson, Ted M Dawson, John Latterra, Jianmin Zhang, Mingyao Ying 2020. Protocols for Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons. **protocols.io**
<https://protocols.io/view/protocols-for-synthetic-mrnas-drive-highly-efficient-9e5h3g6>

MANUSCRIPT CITATION please remember to cite the following publication along with this collection

Xue Y, Zhan X, Sun S, et al. Synthetic mRNAs Drive Highly Efficient iPS Cell Differentiation to Dopaminergic Neurons. *Stem Cells Transl Med.* 2019;8(2):112–123. doi:10.1002/sctm.18-0036

EXTERNAL LINK

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6344911/>

KEYWORDS

ipsc, Atoh1, Ngn2, phosphosite modification, midbrain dopaminergic, differentiation, nonmuscle myosin II, NM-II

LICENSE

— This is an open access collection distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 17, 2019

LAST MODIFIED

Sep 16, 2020

OWNERSHIP HISTORY

Nov 17, 2019  Liz Brydon **Protocols.io**

Sep 16, 2020  Anita Broellochs **protocols.io**

COLLECTION INTEGER ID

29885

GUIDELINES

All results reported should represent at least three independent replications. Perform statistical analysis using Prism software. Post hoc tests should include the Student's *t* test and the Tukey multiple comparison tests as appropriate. For neurophysiological recordings, the recorded data should first be visualized with Clampfit 9.2 and exported to MATLAB for further analysis and plotting. Visualize recording traces with Igor Pro 6.0. All data to be represented as mean \pm SEM.



BEFORE START

Obtain approval to work with human stem cells from an appropriate Institutional Review Board.



SAFETY WARNINGS

Please refer to the Safety Data Sheets (SDS) for safety and environmental hazards.



FILES







Cell culture
Version 1
by Anita Broellochs, protocols.io







mRNA Synthesis and Transfection
Version 1
by Anita Broellochs, protocols.io



- 



iPSC Differentiation
Version 1
by Anita Broellochs, protocols.io
- 



Quantitative Real-Time Polymerase Chain Reaction (qRT-PCR)
Version 1
by Anita Broellochs, protocols.io
- 

Western Blot
Version 1
by Anita Broellochs, protocols.io
- 

Immunofluorescence and Cell Counting
Version 1
by Anita Broellochs, protocols.io
- 

Electrophysiological Recordings
Version 1
by Anita Broellochs, protocols.io
- 

High-Performance Liquid Chromatography (HPLC)
Version 1
by Anita Broellochs, protocols.io
- 

6-Hydroxydopamine (6-OHDA) Treatment and Neurite Tracing in miDA Neurons
Version 1
by Anita Broellochs, protocols.io
- 

Immunoprecipitation (IP) and Mass Spectrometry
Version 1
by Anita Broellochs, protocols.io