

DEC 30, 2022

OPEN ACCESS

dx.doi.org/10.17504/protocol s.io.dm6gpjdy1gzp/v1

Protocol Citation: Kabir Umar. Idris Abdullahi Nasir. Abdurrahman E Ahmad, Abdullahi Hassan Kawo, Abduqadir Magaji Magashi, Abubakar Umar Anka 2022. MecA detection Protocol. protocols.io https://dx.doi.org/10.17504/p

rotocols.io.dm6gpjdy1gzp/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: Dec 29, 2022

Last Modified: Dec 30, 2022

PROTOCOL integer ID:

74561

MecA detection Protocol

Kabir Umar¹, Idris Abdullahi Nasir¹, <u>Abdurrahman E Ahmad</u>², Abdullahi Hassan Kawo³, Abdugadir Magaji Magashi³, Abubakar Umar Anka²

¹Department of Medical Laboratory Science Ahmadu Bello University Zaria

²Department of Medical Laboratory science Ahmadu Bello University;

³Department of Microbiology Bayero Univeristy Kano



Kabir Umar

ABSTRACT

Polymerase chain reaction (PCR) was conducted for the detection of MecA gene from the Staphylococcus aureus Isolated from both healthcare workers and burns and wound patients using the primers listed in table 3.1. A total PCR reaction volume of 25µl containing primers 2 µl, Pre-mix 15 µl, sample 2 µl, and water 6 µl was used, and a cycling condition, with Initial denaturation at 94 °C for 5 min, and 30 cycles of denaturation at 94 °C for 30 sec, annealing at 55 °C for 30 sec and elongation at 72 °C for 1 min and final extension 72 °C 10 min was adopted. PCR amplification was conducted in a thermal cycler (MyCycler Thermal Cycler, Bio-Rad, Portugal) after which agarose gel electrophoresis was run to detect the presence (Band size) or absence of MecA gene as described by Poulsen et al., 2003

MecA gene primers