



Sep 06, 2022

## Primary neuronal cultures

Miguel Da Silva Padilha<sup>1</sup>, Irina Dudanova<sup>1</sup>, F. Ulrich Hartl<sup>1</sup>, Itika Saha<sup>2</sup>, Mark S. Hipp<sup>3,4</sup>

<sup>1</sup>Molecular Neurodegeneration Group, Max Planck Institute of Neurobiology, 82152 Martinsried, Ger many;

<sup>2</sup>Department of Cellular Biochemistry, Max Planck Institute of Biochemistry, Am Klopferspitz 18, 821 52 Martinsried, Germany;

<sup>3</sup>Department of Biomedical Sciences of Cells and Systems, University Medical Center Groningen, University of Groningen, Antonius Deusinglaan, 1, 9713 AV Groningen, The Netherlands;

<sup>4</sup>School of Medicine and Health Sciences, Carl von Ossietzky University Oldenburg, 26129 Oldenburg, Germany



dx.doi.org/10.17504/protocols.io.ewov1ojwklr2/v1

## Felix Kraus

## **ABSTRACT**

This protocol describes the preparation of primary neuronal cultures from E15.5 CD-1 wild type mouse embryos. Experiments involving animal models must be performed in accordance with relevant institutional guidelines and regulations.

DOI

dx.doi.org/10.17504/protocols.io.ewov1ojwklr2/v1

EXTERNAL LINK

https://www.biorxiv.org/content/10.1101/2022.02.18.481043v1.full

PROTOCOL CITATION

Miguel Da Silva Padilha, Irina Dudanova, F. Ulrich Hartl, Itika Saha, Mark S. Hipp 2022. Primary neuronal cultures. **protocols.io** 

https://protocols.io/view/primary-neuronal-cultures-cf7ztrp6

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Itika Saha, Patricia Yuste-Checa, Miguel Da Silva Padilha, Qiang Guo, Roman Körner, Hauke Holthusen, Victoria A. Trinkaus, Irina Dudanova, Rubén Fernández-Busnadiego, Wolfgang Baumeister, David W. Sanders, Saurabh Gautam, Marc I. Diamond, F. Ulrich Hartl, Mark S. Hipp bioRxiv 2022.02.18.481043; doi: https://doi.org/10.1101/2022.02.18.481043



LICENSE
This is an open access protocol distributed under the terms of
the Creative Commons Attribution License, which permits unrestricted use,
distribution, and reproduction in any medium, provided the original author and sourc
are credited
CREATED
Sep 05, 2022

Sep 06, 2022

LAST MODIFIED

PROTOCOL INTEGER ID 69593

- 1 Sacrifice pregnant female mice by cervical dislocation.
- 2 Remove the uterus from the abdominal cavity and place into a 10 cm sterile Petri dish on ice containing dissection medium, consisting of Hanks' balanced salt solution (HBSS) supplemented with 0.01 M HEPES, 0.01 M MgSO4 and 1% penicillin/streptomycin.
- 3 Isolate each embryo, decapitate the heads, remove the brains from the skull and immerse in ice-cold dissection medium.
- 4 Dissect cortical hemispheres, and remove meninges under a dissection microscope.
- 5 Collect the cortices in a 15 mL sterile tube and digest with 0.25% trypsin containing 1 mM ethylenediaminetetraacetic acid (EDTA) and 15  $\mu$ L 0.1% DNAse I for 20 min at 37 °C.
- 6 Stop digestion by removing the supernatant and washing the tissue twice with Neurobasal medium (Invitrogen) containing 5% Fetal Bovine Serum.
- 7 Resuspend the tissue in 2 mL Neurobasal medium and triturate to achieve a single cell suspension.



- 8 Spin cells at 130 x g, remove the supernatant, and resuspend the cell pellet in Neurobasal medium with 2% B-27 supplement (Invitrogen), 1% L-glutamine (Invitrogen) and 1% penicillin/streptomycin (Invitrogen).
- 9 Plate cells at desired density in dishes or coverslips coated with 1 mg/mL poly-D-lysine (Sigma) and 1  $\mu$ g/mL laminin (Thermo Fisher Scientific).