

Upload image


Jul 24, 2020

SOLUTION- 03 - Dextran solution 5%

Marco Cosentino¹, Elisa Storelli¹, Alessandra Luini¹, Massimiliano LM Legnaro¹, Emanuela Rasini¹, Marco Ferrari¹, Franca Marino¹

¹Center for Research in Medical Pharmacology, University of Insubria (Varese, Italy)

1 Works for me dx.doi.org/10.17504/protocols.io.biyukfww



DOI

dx.doi.org/10.17504/protocols.io.biyukfww

DOCUMENT CITATION

Marco Cosentino, Elisa Storelli, Alessandra Luini, Massimiliano LM Legnaro, Emanuela Rasini, Marco Ferrari, Franca Marino 2020. SOLUTION- 03 - Dextran solution 5%. **protocols.io**
dx.doi.org/10.17504/protocols.io.biyukfww

LICENSE

— This is an open access document distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 24, 2020

LAST MODIFIED

Jul 27, 2020

DOCUMENT INTEGER ID

39668

Dextran solution 5% composition g/ml:

Dextran 5 g

Add 100 ml of NaCl 0.15 M.

Storage: 4 °C Fridge 1 - (Room TS08).

Dextran code: 31392, Sigma