

Aug 09, 2022

SAMOSA-Tag - combining in-situ adenine methyltransferase footprinting with transposase-mediated single-molecule sequencing

Scott Nanda¹, Ke Wu¹, Siva Kasinathan², vijay.ramani¹

¹Gladstone Institute for Data Science & Biotechnology;

²Division of Allergy, Immunology, and Rheumatology, Department of Pediatrics, Stanford University

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.3byl4bob2vo5/v1


Ke Wu

ABSTRACT

Here we describe a protocol for SAMOSA-Tag: **T**agmentation-assisted **S**ingle-molecule **A**denine **M**ethylated **O**ligonucleosome **S**equencing **A**ssay - a transposase-mediated strategy that combines our method SMRT-Tag (<https://www.protocols.io/view/FB06B8000CA311ED84710A58A9FEAC02>) for producing low-input PacBio sequencing libraries with *in situ* footprinting to enable multimodal profiling of single chromatin fibers. SAMOSA-Tag facilitates the joint analysis of genetic and epigenetic variations (i.e. somatic variation, nucleosome repeat length, CTCF occupancy, CpG methylation).

In brief, nuclei are methylated *in situ* using the non-specific EcoGII m6dAase, tagged using hairpin-loaded Tn5-adaptors, gap-repaired following DNA purification, and then sequenced on the PacBio Sequel II platform. All steps from live cells to sequencing-ready libraries can be performed in one to two days. We have demonstrated the utility of SAMOSA-Tag by profiling chromatin accessibility, CpG methylation, and CTCF occupancy on low cell numbers in mammalian cell lines.

1 Transposome design & assembly



2 Tagmentation of high-molecular weight gDNA



3 Optimized gap repair & exonuclease clean-up



4 PacBio circular consensus sequencing



An overview of SMRT-Tag schematic design.

DOI

[dx.doi.org/10.17504/protocols.io.3byl4bob2vo5/v1](https://doi.org/10.17504/protocols.io.3byl4bob2vo5/v1)

PROTOCOL CITATION

Scott Nanda, Ke Wu, Siva Kasinathan, vijay.ramani 2022. SAMOSA-Tag - combining in-situ adenine methyltransferase footprinting with transposase-mediated single-molecule sequencing. **protocols.io**
<https://protocols.io/view/samosa-tag-combining-in-situ-adenine-methyltransferase-cd59s896>



KEYWORDS

Long-read sequencing, low-input genomics, chromatin, epigenetics

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 25, 2022

LAST MODIFIED

Aug 09, 2022

PROTOCOL INTEGER ID

67489

GUIDELINES

Sensitive multimodal profiling of native DNA by transposase-mediated single-molecule sequencing

Arjun S Nanda^{1,2*}, Ke Wu^{1*}, Sivakanthan Kasinathan^{3,4*#}, MeganSOstrowski¹, Andrew S Clugston^{5,6}, Ansuman T Satpathy^{4,7,8}, E Alejandro Sweet-Cordero^{5,6}, Hani Goodarzi^{2,6,8,9}, Vijay Ramani^{1,2,6,9#}

¹Gladstone Institute for Data Science & Biotechnology, Gladstone Institutes, San Francisco, CA 94158

²Department of Biochemistry & Biophysics, UCSF, San Francisco, CA 94143

³Division of Allergy, Immunology, and Rheumatology, Department of Pediatrics, Stanford University, Stanford, CA 94305

⁴Gladstone-UCSFInstituteForGenomic Immunology, Gladstone Institutes, San Francisco, CA 94143

⁵Department of Pediatrics, UCSF, San Francisco, CA 94143

⁶Helen-Diller Cancer Center, San Francisco, CA 94158

⁷Department of Pathology, Stanford University, Stanford, CA 94305

⁸Parker Institute for Cancer Immunotherapy, San Francisco, CA 94129

⁹Bakar Computational Health Sciences Institute, San Francisco, CA 94143

*These authors contributed equally to this work.

#correspondence to skas@stanford.edu or vijay.ramani@gladstone.ucsf.edu

MATERIALS TEXT

- LoBind microfuge tubes (1.5 mL)
- PCR tubes (0.2 mL)
- Tn5^{R27S,E54K,L372P} purified enzyme (QB3 Core)

- Barcoded SMRT-Tag Adaptor sequence (IDT, see supp. table)
- Cell suspension. We have tried human K562 cells, murine embryonic stem cells (mESCs), osteosarcoma cells (OS152)
- Distilled, deionized or RNase-free H₂O (dH₂O; e.g. Invitrogen, cat. no. 10977023)
- 1 M Hydroxyethyl piperazineethanesulfonic acid (HEPES; Gibco, cat. no. 15630080)
- 2 M Potassium Chloride (KCl; Invitrogen, cat. no. AM9640G)
- 1 M Magnesium Chloride (MgCl₂; G-Bioscience, cat. no. R004)
- Triton X-100 (Sigma-Aldrich, cat. no. X100)
- Glycerol (Invitrogen, cat. no. 15514011)
- cOmplete EDTA-free Protease Inhibitor Cocktail (Roche, cat. no. 4693132001)
- 1 M Tris-HCl pH 7.5 (Invitrogen, cat. no. 15567027)
- 1 M Tris-HCl pH 8.0 (Invitrogen, cat. no. 15568025)
- Elution buffer (1x EB; Qiagen, cat. no. 19086)
- 5 M Sodium chloride (NaCl; Invitrogen, cat. no. AM9760G)
- 1 M Spermidine (Sigma-Aldrich, cat. no. S0266)
- Dimethylformamide (DMF; Sigma-Aldrich, cat. no. 319937)
- Phosphate-buffered saline pH 7.4 (1x PBS; Gibco, cat. no. 10010023)
- Dimethyl sulfoxide (DMSO; VWR, cat. no. 97063-136)
- EcoGII Methyltransferase (25,000U/mL New England Biolabs, cat. no. M0603B-HC1)
- S-Adenosyl methionine (SAM; 32mM New England Biolabs, cat. no. B9003S)
- 5% Digitonin (Thermo Scientific, cat. no. BN2006)
- RNase A (10mg/mL Fisher Scientific, cat. no. EN0531)
- 0.5 M Ethylenediaminetetraacetic acid (EDTA; Invitrogen, cat. no. 15575020)
- 10% Sodium dodecyl sulfate (SDS; Invitrogen, cat. no. 15553027)
- Proteinase K (20 mg/ml Ambion, cat. no. AM2548)
- dNTP mix (100mM Thermo Fisher Scientific, cat. no. R0181)
- 10x Taq DNA Ligase Reaction Buffer (New England Biolabs, cat. no. B0208S)
- Phusion High-Fidelity DNA Polymerase (2,000U/mL New England Biolabs, cat. no. M0530L)
- Taq DNA Ligase (40,000U/mL New England Biolabs, cat. no. M0208S)
- 10x NEBuffer 2 (New England Biolabs, cat. no. B7002S)
- Exonuclease III (100,000U/mL New England Biolabs, cat. no. M0206L)
- Ethanol (Decon Labs, cat. no. V1016)
- SPRI beads (e.g. Agencourt AMPure XP, Beckman Coulter, cat. no. A63880)
- AMPure PB beads (Pacific Biosciences, cat. no. 100-265-900)
- FlowMi cell strainer (Sigma-Aldrich, cat. no. BAH136800040)
- Qubit 1X dsDNA HS Assay Kit (Invitrogen, cat. no. Q33231)
- Bioanalyzer High Sensitivity DNA Reagent Kit (Agilent Technology, cat. no. 5067-4627)
- Bioanalyzer DNA 12000 Reagent Kit (Agilent Technology, cat. no. 5067-1508)

Supplementary Table - SMRT-Tag Adaptor sequence:

A	B
Barcode Name	Sequence
SMRT-A_bc-none	/5Phos/CTG TCT CTT ATA CAC ATC TAT CTC TCT CTT TTC CTC CTC CTC CGT TGT TGT TGT TGA GAG AGA TAG ATG TGT ATA AGA GAC AG

SMRT-A_bc001	/5Phos/CTG TCT CTT ATA CAC ATC TTT CTT CCG ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT CGG AAG AAA GAT GTG TAT AAG AGA CAG
SMRT-A_bc003	/5Phos/CTG TCT CTT ATA CAC ATC TTT CCA CAC ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT GTG TGG AAA GAT GTG TAT AAG AGA CAG
SMRT-A_bc006	/5Phos/CTG TCT CTT ATA CAC ATC TTT GTC GCA ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT TGC GAC AAA GAT GTG TAT AAG AGA CAG
SMRT-A_bc010	/5Phos/CTG TCT CTT ATA CAC ATC TTT AGC TGC ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT GCA GCT AAA GAT GTG TAT AAG AGA CAG
SMRT-A_bc011	/5Phos/CTG TCT CTT ATA CAC ATC TTC CTA AGG ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT CCT TAG GAA GAT GTG TAT AAG AGA CAG
SMRT-A_bc012	/5Phos/CTG TCT CTT ATA CAC ATC TTC CGT TGT ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT ACA ACG GAA GAT GTG TAT AAG AGA CAG
SMRT-A_bc013	/5Phos/CTG TCT CTT ATA CAC ATC TTC GAA TCG ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT CGA TTC GAA GAT GTG TAT AAG AGA CAG
SMRT-A_bc014	/5Phos/CTG TCT CTT ATA CAC ATC TTC ACT GTG ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT CAC AGT GAA GAT GTG TAT AAG AGA CAG
SMRT-A_bc015	/5Phos/CTG TCT CTT ATA CAC ATC TTG CAG GAT ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT ATC CTG CAA GAT GTG TAT AAG AGA CAG
SMRT-A_bc016	/5Phos/CTG TCT CTT ATA CAC ATC TTA TGG CGT ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT ACG CCA TAA GAT GTG TAT AAG AGA CAG
SMRT-A_bc017	/5Phos/CTG TCT CTT ATA CAC ATC TTA CCG ACT ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT AGT CGG TAA GAT GTG TAT AAG AGA CAG
SMRT-A_bc018	/5Phos/CTG TCT CTT ATA CAC ATC TTA CAA GCC ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT GGC TTG TAA GAT GTG TAT AAG AGA CAG
SMRT-A_bc019	/5Phos/CTG TCT CTT ATA CAC ATC TCT GAC CAA ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT TTG GTC AGA GAT GTG TAT AAG AGA CAG
SMRT-A_bc020	/5Phos/CTG TCT CTT ATA CAC ATC TCC TCT CTA ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT TAG AGA GGA GAT GTG TAT AAG AGA CAG

SMRT-A_bc021	/5Phos/CTG TCT CTT ATA CAC ATC TCC TGT AAC ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT GTT ACA GGA GAT GTG TAT AAG AGA CAG
SMRT-A_bc022	/5Phos/CTG TCT CTT ATA CAC ATC TCC GCA TAA ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT TTA TGC GGA GAT GTG TAT AAG AGA CAG
SMRT-A_bc023	/5Phos/CTG TCT CTT ATA CAC ATC TCA AGT GGA ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT TCC ACT TGA GAT GTG TAT AAG AGA CAG
SMRT-A_bc024	/5Phos/CTG TCT CTT ATA CAC ATC TGT GCA TTC ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT GAA TGC ACA GAT GTG TAT AAG AGA CAG
SMRT-A_bc025	/5Phos/CTG TCT CTT ATA CAC ATC TGG CTT CAT ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT ATG AAG CCA GAT GTG TAT AAG AGA CAG
SMRT-A_bc026	/5Phos/CTG TCT CTT ATA CAC ATC TGG AAC TAC ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT GTA GTT CCA GAT GTG TAT AAG AGA CAG
SMRT-A_bc027	/5Phos/CTG TCT CTT ATA CAC ATC TGA CGT TAG ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT CTA ACG TCA GAT GTG TAT AAG AGA CAG
SMRT-A_bc028	/5Phos/CTG TCT CTT ATA CAC ATC TGA GTG TCT ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT AGA CAC TCA GAT GTG TAT AAG AGA CAG
SMRT-A_bc029	/5Phos/CTG TCT CTT ATA CAC ATC TGA AGA AGG ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT CCT TCT TCA GAT GTG TAT AAG AGA CAG
SMRT-A_bc030	/5Phos/CTG TCT CTT ATA CAC ATC TAA CAC CTC ATC TCT CTC TTT TCC TCC TCC TCC GTT GTT GTT GTT GAG AGA GAT GAG GTG TTA GAT GTG TAT AAG AGA CAG

SAFETY WARNINGS

Digitonin and DMF are both toxic and extra care should be taken. Use full PPE including a mask, lab coat, and gloves while handling any amount of either chemical.

Reagent Setup

1 Nuclear Lysis Buffer (prepared on ice):

A	B	C	D
Reagent (final conc.)	Stock conc.	Volume	Note
20 mM HEPES	1M, pH 7.5	200µL	
10 mM KCl	2M	50µL	
1 mM MgCl ₂	1M	10µL	
0.1% Triton X-100	10%	100µL	
20% glycerol	100%	2mL	
1x Protease Inhibitor	25x	400µL	Add 25x Roche protease inhibitor immediately before use
dH ₂ O	-	up to 10mL	
Total Volume		10mL	

Buffer M (prepared on ice):

A	B	C
Reagent (final conc.)	Stock conc.	Volume
15mM Tris-HCl, pH 8.0	1M	150µL
15mM NaCl	5M	30µL
60mM KCl	2M	300µL
0.5mM spermidine	1M	5µL
dH ₂ O	-	up to 10mL
Total Volume		10mL

Nuclei Storage Buffer (prepared on ice):

A	B	C
Reagent (final conc.)	Stock conc.	Volume
20mM HEPES pH 7.5	1M	200µL
150mM NaCl	5M	300µL
0.25mM spermidine	1M	2.5µL
0.1X Protease Inhibitor	25X	40µL
dH ₂ O	-	up to 10mL
Total Volume		10mL

5x TD Premix (prepared at RT):

A	B	C	D
Reagent (final conc.)	Stock conc.	Volume	Note
10mM Tris-HCl, pH 7.5	1M	10µL	
5mM MgCl ₂	1M	5µL	
10% DMF	100%	-	Add directly into Omni-ATAC Reagent Mix at step 20
dH ₂ O	-	185µL	
Total Volume		200µL	

SMRT-Tn5 adaptor assembly:

- we recommend following the "Annealing SMRT-Tag adaptors" and "Assembling SMRT-Tn5 transposomes (Tn5 loaded with SMRT-Tag adaptors)" procedures in SMRT-Tag protocol:
<https://www.protocols.io/view/FB06B8000CA311ED84710A58A9FEAC02>

Nuclei isolation from fresh cells 45m

- 2 Harvest fresh cell culture(s) in a conical centrifuge tube at room temperature and count cells. This protocol can accommodate up to a total of 2 million mammalian cells (e.g. human K562, murine embryonic stem cells) that can be divided into 10-30k nuclei per sample for up to 16 samples to be sequenced.
- 3 Centrifuge a 1 to 2 million cell suspension at 300xg at 4°C for 5 min, and aspirate the supernatant.

[NOTE] For robust cell lines including cancer cell lines, 10 minutes of centrifugation can increase recovery.
- 4 Resuspend the cell pellet in 1mL ice-cold Nuclear Lysis Buffer. Pipet up and down gently 10x using wide-bore tips.
(Use approximately 1mL Nuclear Lysis Buffer for up to 10 million cells.)
- 5 Incubate on ice for exactly 5 min.

- 6 Centrifuge at 600xg at 4°C for 5 min, and aspirate the supernatant.

[NOTE] For robust cell lines including cancer cell lines, 10 minutes of centrifugation can increase recovery.

- 7 Wash nuclei with 1mL Buffer M.

[cell count 0] Determine # of nuclei from a 10 µL aliquot (1:1 dilution) via hemocytometer / cell counter (10 µL nuclei suspension + 10 µL trypan blue)

- 8 Centrifuge again at 600xg at 4°C for 5 min, and aspirate the supernatant.

[NOTE] For robust cell lines including cancer cell lines, 10 minutes of centrifugation can increase recovery.

- 9 Proceed immediately to nuclei freezing or EcoGII methylation.

(Optional) Nuclei freezing 10m

- 10 Isolated nuclei (from step 8) may be slowly frozen by resuspending in nuclei freezing media containing 9 volumes of Nuclei Storage Buffer and 1 volume of DMSO, and further divided into 1 mL per cryogenic vial (~1-2 million nuclei), then stored at -80°C. We have found a cooling rate close to -1°C/min is sufficient for preserving nuclei integrity.

(Optional) Thawing previously frozen nuclei 20m

- 11 If frozen nuclei are used as input, thaw ~1-2 million nuclei in a water bath at 37°C for 1-2 min until most of the solution is a slurry. Expect 30-50% recovery after each wash step. This protocol is designed for up to 1-2 million frozen nuclei that can be further divided into 10-30k nuclei per sample for up to 16 samples to be sequenced but can be adapted for lower numbers of nuclei per replicate, or fewer replicates, as needed.

[CRITICAL] Nuclei are thawed in 10% DMSO and need to be washed rapidly in order to prevent lysis.

- 12 Use a P1000 with wide-bore pipette tips to gently transfer nuclei solution to a 1.5mL LoBind tube.

- 13 Centrifuge at 600xg at 4°C for 5 min, and aspirate the supernatant. Leave a tiny volume around the pellet is fine.

[NOTE] For robust cell lines including cancer cell lines, 10 minutes of centrifugation can increase recovery.

- 14 Resuspend pelleted nuclei in 1mL Buffer M by gentle pipetting using a P1000 with wide-bore tips.

[cell count 1] Determine # of nuclei from a 10 μ L aliquot (1:1 dilution) via hemocytometer / cell counter (10 μ L nuclei suspension + 10 μ L trypan blue)

- 15 Centrifuge again at 600xg at 4°C for 5 min, and aspirate the supernatant.

[NOTE] For robust cell lines including cancer cell lines, 10 minutes of centrifugation can increase recovery.

- 16 Proceed immediately to EcoGII methylation.

in-situ SAMOSA (or EcoGII Methylation)

40m

- 17 Resuspend the nuclei (either freshly isolated nuclei from step 8 or previously-frozen nuclei from step 15) in 400 μ L buffer containing Buffer M + 1mM SAM (e.g. 12.5 μ L of 32mM SAM + 387.5 μ L Buffer M).

[cell count 2] Determine # of nuclei from a 10 μ L aliquot (1:1 dilution) via hemocytometer / cell counter (10 μ L nuclei suspension + 10 μ L trypan blue)

- 18 Split solution into 2 x 200 μ L reactions, one for methylated (+M) and the other for unmethylated (-M) nuclei.

[NOTE] If processing a sample where unmethylated data has been previously obtained, or samples are precious, the unmethylated condition can be skipped.

[PAUSE] Determine if you have enough nuclei for downstream processing, assuming 50% loss after EcoGII methylation. (i.e 8 samples x 2 methylation states x 30k per rxn = ~ 500K nuclei required)

- 19 Add 10 μ L of EcoGII (25U/ μ L) to +M rxn. For -M rxn, do not add EcoGII. Mix well by slowly pipetting with wide-bore tips or gently flicking to mix.

- 20 Incubate both reactions at 37°C for 30 min in a thermomixer with 300rpm shaking every 2 min. At the 15 min time point, supplement with 1 μ L of 32mM SAM stock in both reactions.

- 21 Terminate the reactions by placing them back on ice.

Tagmentation of SAMOSA-treated nuclei 1h 30m

- 22 Prepare Omni-ATAC Reagent Mix on ice. Prepare a little extra (5 samples worth) for sample dilution.

A	B	C
Reagent (final conc.)	Stock conc.	Volume (per sample)
1x TD premix	5x	5µL
0.33x PBS	1x	16.5µL
10% DMF	100%	5µL
0.01% digitonin	1%	0.5µL
0.1% Tween-20	10%	0.5µL
dH2O	-	up to 50µL
Total Volume		50µL

Each tagmented sample of 10-30k nuclei requires a uniquely barcoded SMRT-Tag adaptor. This protocol is designed for up to 16 samples given at least a total of 500k nuclei recovered after EcoGII methylation.

- 23 Centrifuge 2 x 200 µL reactions after EcoGII methylation at 600xg at RT for 10 min, and carefully pipette to remove the supernatant. Leave a tiny volume around the pellet is fine.

[CRITICAL] Nuclei may be smeared out across the LoBind tube wall instead of nicely pelleted at the bottom.

- 24 Resuspend nuclei gently in 250 µL Omni-ATAC Reagent Mix per reaction.

- 25 Using a P1000 with a wide-bore tip, set the pipette to 300 µL and aspirate up the entire volume per reaction. Insert the tip into a FlowMi cell strainer, and then slowly push the reaction through the strainer into a new LoBind tube. There will be bubbles ejecting upward, so make sure filtered tips are used.

[cell count 3] Determine # of nuclei post filtration from a 10 µL aliquot (1:1 dilution) via hemocytometer / cell counter (10 µL nuclei suspension + 10 µL trypan blue) for both +M and -M reactions

[CRITICAL] QC check - if nuclei still look clumped together, make a note - your count is likely NOT accurate, and may require adjusting the Tn5 amount or nuclei input per sample. Nuclei

clumping tends to underestimate the actual count number.

- 26 Split each of the +M and -M reactions into multiple tubes as needed such that ~ 10K - 30K nuclei are in each tube.

[CRITICAL] The majority of samples to be sequenced should be +M which are used to profile chromatin accessibility.

- 27 To generate an average of 3kb libraries, add 2 μL of loaded SMRT-Tn5 adaptor stock (~ 9.4 μM monomer) per sample and fill the total volume up to 50 μL with Omni-ATAC reagent mix.

[OPTIONAL] QC check - Prepare 1 untagged negative control for each of +M and -M reactions to compare against tagged samples. A total of two negative controls: +M -tag and -M -tag. Replace SMRT-Tn5 adaptor volume with Omni-ATAC reagent mix.

- 28 Incubate all tagmentation reactions at 55°C for 1 hr and hold at 4°C.

Tagmentation termination 1h 30m

- 29 Treat the tagmented nuclei with 10 μL of RNase A (10 $\mu\text{g}/\mu\text{L}$) and incubate at 37°C for 15 min on a thermomixer with 300rpm mixing every 2 min.

[CRITICAL] RNase pre-treatment is critical to ensure that only DNA is extracted from the reaction.

- 30 Prepare the Termination Lysis Reagent Mix about 10 min before the tagmentation is complete.

A	B	C
Reagent (final conc.)	Stock conc.	Volume (per sample)
Proteinase K (6.66mg/mL)	20mg/mL	2.5 μL
3.33% SDS	10%	2.5 μL
167mM EDTA	0.5M	2.5 μL
Total Volume		7.5 μL

[CRITICAL] Equilibrate Proteinase K at RT before adding to the mixture. Do NOT keep it on ice, or the SDS mixture will precipitate with EDTA. Vortex EDTA stock solution prior to use. If the mixture appears cloudy, warm it up at RT until it appears homogenous again.

- 31 After RNase treatment, add 7.5 μL of termination lysis reagent mix to each sample and mix well.

- 32 Incubate at 60°C on a thermomixer with 1000rpm continuous shaking for at least 1 hr, up to 2 hr for improved lysis.

gDNA purification with SPRI cleanup

1h


- 33 Remove samples from thermomixer and incubate on ice for 1 min to cool down samples to RT.
- 34 Add 2x volume of resuspended, room-temperature SPRI beads. (i.e. For a 50 µL tagmentation reaction, sample volume should be 50 µL tagmentation reaction + 10 µL RNase A + 7.5 µL termination lysis MM ~ 67.5 µL, so 2X ~ 135 µL). Mix well using a P200 with wide-bore tips.
- 35 Incubate the bead-mixed samples in a thermomixer at 23°C for 30 min with interval mixing @350 rpm (1 min on, 3 min off) to keep the beads resuspended.
- 36 Spin down quickly. Place on magnet and allow to clear before carefully withdrawing supernatant.
- 37 Wash beads twice with 80% ethanol for 30 sec. Add ethanol to the wall opposite of beads. Withdraw supernatant and after a quick spin, remove the remaining liquid with a P20 pipette. Beads take ~ 30 s to 1 min to dry. Do not over-dry the beads - dried beads appear to be fragmented with cracks in light brown.
- 38 Remove from the magnet stand, and gently resuspend the bead pellet in 20 µL 1x EB (or 10 mM Tris-HCl, pH 8.5) to elute the extracted DNA. Once beads are resuspended, mix well.
- 39 Incubate the samples in a Thermomixer at 37°C for 15 min with interval mixing @350 rpm (1 min on, 3 min off) to promote increased DNA elution.
- 40 Spin samples down quickly. Place on magnet and allow samples to clear before carefully transferring the supernatant to a new LoBind tube.

[Pause] Purified gDNA samples can be stored at 4°C for up to two weeks.

- 41 Determine concentrations of purified DNA using Qubit 1x High Sensitivity DNA Assay. Expect at least 160ng of DNA recovered from 10-30k nuclei.

- 42 **[OPTIONAL]** Determine the fragment size distribution using Agilent 2100 Bioanalyzer High Sensitivity DNA Assay.

Reference fig. 1: Extracted gDNA from tagmented nuclei



The Bioanalyzer trace shows major peaks around 8-10kb along with a trailing tail <1kb.

[OPTIONAL] If BioAnalyzer trace shows a trailing tail <1kb (e.g. reference fig. 1), we recommend doing an additional 0.6x Ampure PB cleanup to enrich for fragments > 500bp.

- 43 **[Optional]** QC check - Visualize the SMRT-Tn5 transposome efficiency by running 2 untagged controls against tagged samples on a 0.5% agarose gel coupled with the NEB 1kb-Extend ladder. Expect to see the majority of bands clustered towards high molecular weight >48.5kb for untagged samples, whereas a long smear should be observed for tagged samples. Untagged samples will not be sequenced as they do not contain SMRT-Tag adaptors.

Gap Repair

2h

- 44 Prepare the Gap-Repair Reagent Mix.

A	B	C
Reagent (final conc.)	Stock conc.	Volume (per sample)
tagmented sample (up to 160ng)	-	12 μ L *
dNTP mix (0.8mM)	8mM	2 μ L
1x Taq DNA Ligase Rxn Buffer	10x	2 μ L
NEB Phusion (0.1U)	2U/ μ L	1 μ L
NEB Taq DNA Ligase (4U)	40U/ μ L	2 μ L
dH2O	-	1 μ L
Total Volume		20 μ L

*Normalize each tagmented sample to 160ng in 12 μ L with 1x EB.

- 45 Incubate the gap-repair reactions at 37°C for 1hr and hold at 4°C.
- 46 Add 2x volume of resuspended, room-temperature SPRI beads to clean up the reactions. Mix well using a P200 with wide-bore tips.
- 47 Incubate the bead-mixed samples in a thermomixer at 23°C for 15 min with interval mixing @350 rpm (1 min on, 3 min off) to keep the beads resuspended.
- 48 Spin down quickly. Place on magnet and allow to clear before carefully withdrawing supernatant.
- 49 Wash beads twice with 80% ethanol for 30 sec. Add ethanol to the wall opposite of beads. Withdraw supernatant and after a quick spin, remove the remaining liquid with a P20 pipette. Beads take ~ 30 s to 1 min to dry. Do not over-dry the beads - dried beads appear to be fragmented with cracks in light brown.

- 50 Remove from the magnet stand, and gently resuspend the bead pellet in 12 µL 1x EB (or 10 mM Tris-HCl, pH 8.5) to elute the extracted DNA. Once beads are resuspended, mix well.
- 51 Incubate the samples in a Thermomixer at 37°C for 10 min with interval mixing @350 rpm (1 min on, 3 min off) to promote increased DNA elution.
- 52 Spin samples down quickly. Place on magnet and allow samples to clear before carefully transferring the supernatant to a new LoBind tube.
- 53 **[Optional]** QC check - take 1 µL aliquot and determine concentrations of gap-repaired samples using Qubit 1x High Sensitivity DNA Assay. Expect a 70% or above recovery from tagged sample input in step 44.

Exonuclease Digestion

2h

- 54 Prepare the Exo-Digest Reagent Mix.

A	B	C
Reagent (final conc.)	Stock conc.	Volume (per sample)
gap-repaired sample	-	12µL*
NEBuffer 2 (1x)	10x	1.5µL
Exonuclease III (100U)	100U/µL	1µL
dH2O	-	0.5µL
Total Volume		15µL

*Diluted up to 12µL with 1x EB for volume <12µL.

- 55 Incubate the exo-digested reactions at 37°C for 1hr and hold at 4°C.
- 56 Add 2x volume of resuspended, room-temperature SPRI beads to clean up the reactions. Mix well using a P200 with wide-bore tips.

Incubate the bead-mixed samples in a thermomixer at 23°C for 15 min with interval mixing

- 57 @350 rpm (1 min on, 3 min off) to keep the beads resuspended.
- 58 Spin down quickly. Place on magnet and allow to clear before carefully withdrawing supernatant.
- 59 Wash beads twice with 80% ethanol for 30 sec. Add ethanol to the wall opposite of beads. Withdraw supernatant and after a quick spin, remove the remaining liquid with a P20 pipette. Beads take ~ 30 s to 1 min to dry. Do not over-dry the beads - dried beads appear to be fragmented with cracks in light brown.
- 60 Remove from the magnet stand, and gently resuspend the bead pellet in 12 µL 1x EB (or 10 mM Tris-HCl, pH 8.5) to elute the extracted DNA. Once beads are resuspended, mix well.
- 61 Incubate the samples in a Thermomixer at 37°C for 10 min with interval mixing @350 rpm (1 min on, 3 min off) to promote increased DNA elution.
- 62 Spin samples down quickly. Place on magnet and allow samples to clear before carefully transferring the supernatant to a new LoBind tube.

Library QC and Sequencing


- 63 Determine library concentration using Qubit 1x High Sensitivity DNA Assay.

Determine the size distribution using Agilent 2100 Bioanalyzer High Sensitivity DNA Assay.

[CRITICAL] Convert the BioA trace from its default unit (fluorescent unit) to "molarity per length" using the R package "bioanalyzeR". This step is critical for the evaluation of accurate library sizing required for the optimal loading on PacBio Sequel II.

See more instructions on package installation and usage of "bioanalyzeR":
<https://stanford.edu/~jwfoley/bioanalyzeR.html>

Reference fig. 2: SAMOSA-Tag library



The BioAnalyzer trace shows an average size of ~6kb. After unit conversions to "molarity per length" using "bioanalyzerR", the average library size shifts down to ~3kb which is used as the final library size for Sequel II loaded with Sequel II Binding Kit 2.1.

- 64 Pool barcoded libraries equimolarly and sequence on PacBio Sequel II 8M SMRTcells. We recommend using Sequel II Binding Kit 2.1 for the best loading result, as well as a 30h movie time, 2h pre-extension, and 4h immobilization time.