

JUN 19, 2023

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.kxygx3mkog8j/v1

Protocol Citation: Nathan Heath Patterson, ellie.l.pingry, Katerina V Djambazova, Melissa Farrow, Raf Van De Plas, Jeff Spraggins 2023. Kidney Functional Tissue Unit (FTU) Segmentation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygx3mkog8j/v1>

MANUSCRIPT CITATION: Patterson, N. H.; Neumann, E. K.; Sharman, K.; Allen, J.; Harris, R.; Fogo, A. B.; De Caestecker, M.; Caprioli, R. M.; Van De Plas, R.; Spraggins, J. M. Autofluorescence Microscopy as a Label-Free Tool for Renal Histology and Glomerular Segmentation. <https://doi.org/10.1101/2021.07.16.452703>.

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
 We use this protocol and it's working

Kidney Functional Tissue Unit (FTU) Segmentation

Katerina V

Nathan Heath Patterson¹, ellie.l.pingry¹, Djambazova¹,


Raf Van De

Melissa Farrow¹, Plas², Jeff Spraggins¹

¹Vanderbilt University; ²Delft University of Technology

VU Biomolecular Multimodal Imaging Center / Spraggins Research Group

KPM



Katerina V Djambazova

Vanderbilt University

ABSTRACT

This protocol explains how to apply a segmentation model to autofluorescence microscopy images to find kidney functional tissue units (FTUs). Currently the model allows for segmentation of glomeruli, proximal tubules, thick ascending limb, distal tubules, and collecting ducts.

MATERIALS

For performance, a workstation with at least 32 GB of memory is advised and a GPU with 12 GB memory is required for any GPU-based prediction.

Created: Jun 16, 2023

Last Modified: Jun 19, 2023

PROTOCOL integer ID:
83532

Apply kidney FTU segmentation to autofluorescence images

- 1 Configure the segmentation prediction by creating the dataset [.yaml file](#). Add the path to the saved model and images in the configuration file.
- 2 Create the python environment needed to run wsimap. See instructions at the [wsimap GitHub](#).
- 3 Run the model using the command line from the wsimap directory like so:

```
python scripts/config/instance-predict-from-config-newmodels-multi.py  
"/path/to/configuration-file-from-step2.yaml"
```
- 4 Collect .geojson segmentations in output folder specified in .yaml file from Step 2. These can be visualized in QuPath.