



APR 01, 2024

# 🌐 CODA (part 1): setting up environment and preparing sample dataset | HuBMAP | JHU-TMC

Kyu Sang Han<sup>1</sup>, Pei-Hsun Wu<sup>1</sup>, Joel Sunshine<sup>2</sup>, Sashank Reddy<sup>2</sup>, Ashley Kiemen<sup>2</sup>, Denis Wirtz<sup>1,2</sup>

<sup>1</sup>Johns Hopkins University; <sup>2</sup>Johns Hopkins Medicine

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

TMC - Johns Hopkins University

OPEN  ACCESS



Kyu Sang Han  
Johns Hopkins University

## ABSTRACT

CODA workflow part 1. setting up environment and preparing dataset

## DOI:

[dx.doi.org/10.17504/protocols.io.q26g71rpkgwz/v1](https://dx.doi.org/10.17504/protocols.io.q26g71rpkgwz/v1)

**Protocol Citation:** Kyu Sang Han, Pei-Hsun Wu, Joel Sunshine, Sashank Reddy, Ashley Kiemen, Denis Wirtz 2024. CODA (part 1): setting up environment and preparing sample dataset | HuBMAP | JHU-TMC. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g71rpkgwz/v1>

## MANUSCRIPT CITATION:

Kiemen, A.L., Braxton, A.M., Grahn, M.P. *et al.* CODA: quantitative 3D reconstruction of large tissues at cellular resolution. *Nat Methods* **19**, 1490–1499 (2022). <https://doi.org/10.1038/s41592-022-01650-9>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Mar 28, 2024

**Last Modified:** Apr 01, 2024

**PROTOCOL integer ID:** 97507

**Keywords:** CODA, deeplearning, tissuelabeling, segmentation

**Funders Acknowledgement:**

Institute of Arthritis and  
Musculoskeletal and Skin  
Diseases  
Grant ID: U54AR081774  
National Cancer Institute  
Grant ID: U54CA143868

## Software requirements

- 1** MATLAB [MATLAB \(mathworks.com\)](#)  
Image processing toolbox [Image Processing Toolbox - MATLAB \(mathworks.com\)](#)  
Deep learning toolbox [Deep Learning Toolbox - MATLAB \(mathworks.com\)](#)  
MATLAB - Resnet50 model [MATLAB resnet50 \(mathworks.com\)](#)
- 2** Aperio ImageScope [Aperio ImageScope \(leicabiosystems.com\)](#)
- 3** FIJI ImageJ [Fiji Downloads \(imagej.net\)](#)

## Download Source code

- 4 Codes are available at the following GitHub: [CODA Github Repository](#)

## Download Sample dataset

- 5 Here, we discuss application to a sample dataset “lungs”, containing 150 serial histological images. Download the sample dataset (serial images and sample annotations) here: [Lung Sample Dataset on Google Drive](#)
- 6 Images are .ndpi format and were scanned at 20x magnification (approximately 0.5 micron / pixel resolution), spaced 10 micron apart. Save the images in a local drive folder (e.g. `\\Users\\Ashley\\Documents\\lungs`).
- 7 Filenames for each image should be created such that tissue sections are read consecutively by Matlab. Therefore, include zero-padding in numerical indices.

CORRECT FILENAMES: `lungs_001.ndpi, lungs_002.ndpi, ..., lungs_011.ndpi`

INCORRECT FILENAMES (no zero padding): `lungs_1.ndpi, lungs_2.ndpi, ..., lungs_11.ndpi`

## Create downsampled copies of high-resolution images

- 8 The function **create\_downsampled\_tif\_images** will create downsampled copies of the .ndpi files by directly loading each high-resolution images in tiles and down sampling it to the desired pixel resolutions.
- 9 First, decide the resolution of the images you want to create. Here, we create images of 1 micron / pixel, 2 microns / pixel, and 10 micron / pixel resolution:

```
ds=[1 2 10];
```

**10** Next, decide on the name of the output folders for each of the downsampled images you create. Here, we will save the images downsampled to 1 micron / pixel in a folder named “10x”, the images downsampled to 2 micron / pixel in a folder named “5x,” and the images downsampled to 10 micron / pixel in a folder named “1x.”

```
subfolders=["10x" "5x" "1x"];
```

**11** Finally, call the function: **create\_downsampled\_tif\_images(pth,ds,subfolders);**

**12** Using this function, you will make two subfolders within the original folder containing the .ndpi images. One subfolder named “10x” containing the 20x images downsampled by a factor of 2. The other subfolder named “1x” containing the 20x images downsampled by a factor of 20. Most calculations will be performed on these tif images. Note: here we use 10x and 1x for example, but other resolutions could be created as desired.

```
pth10x=[pth,'10x'];
```

```
pth1x=[pth,'1x'];
```

**\*\*Note: If this code fails due to memory constraints on your computer, try python Openslide.**