


MAY 16, 2023

## Using Amira to generate cell body masks in fluorescent light sheet microscopy images of EASI-FISH samples

Grace Park<sup>1</sup>, Aubrey Weigel<sup>1</sup>, Alyson Petruncio<sup>1</sup>,  
Claire Managan<sup>1</sup>, Mark Eddison<sup>1</sup>, Gudrun Ihrke<sup>1</sup>, Gudrun Ihrke<sup>1</sup>

<sup>1</sup>HHMI

Janelia



Grace Park

OPEN ACCESS

**Protocol Citation:** Grace Park, Aubrey Weigel, Alyson Petruncio, Claire Managan, Mark Eddison, Gudrun Ihrke, Gudrun Ihrke 2023. Using Amira to generate cell body masks in fluorescent light sheet microscopy images of EASI-FISH samples.

protocols.io

<https://protocols.io/view/using-amira-to-generate-cell-body-masks-in-fluorescent-light-sheet-microscopy-images-of-easi-fish-samples>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Feb 08, 2023

**Last Modified:** May 16, 2023

**PROTOCOL integer ID:**  
76674

### ABSTRACT

This protocol describes how [CellMap Project Team](#) created the masks for cell bodies used in EASI-FISH (Expansion-Assisted Iterative Fluorescence In Situ Hybridization) spot-counting analysis. The [Radial Symmetry-FISH \(RS-FISH\)](#) software enables single molecule spot detection in two and three dimensions with high precision. RS-FISH can be installed and run as a Fiji plugin. This protocol features Amira 3D 2021.1 classic segmentation workroom and annotation techniques in the protocol: [Using Amira to manually segment organelles in vEM for machine learning](#).

## Loading the Data

- 1 This protocol describes how the masks for cell bodies were generated to assist the RS-FISH

spot-counting analysis. Please reference the [Using Amira to manually segment organelles in vEM for machine learning](#) protocol to find details of the annotation techniques and Amira software.

The raw images are in .dzi format which can be directly opened in Amira by dragging and dropping the file from the folder to the Amira software (Fig 1.1).




Fig 1.1

Since the .dzi file is large in size, Amira will display a pop-up window with loading policy options to choose from (Fig 1.2). Select the third option, "Read complete volume into memory." and hit "ok" to proceed. If the raw image is split into two datasets (i.e., left and right), load both .dzi files onto Amira. The data may take a couple of minutes to load.




Fig 1.2

The "Multi-Channel-Field" consists of four different channels (Fig 1.3). The mask annotation will use channel 0 with the cell bodies as the image source. Other channels may be deleted (Fig 1.4).




Fig 1.3




Fig 1.4




Fig 1.5

Once the appropriate .czi files and channels are successfully loaded onto the "Project" workroom, the image can be viewed from the "Segmentation" workroom (Fig 1.5). Adjust the

toolbar from the "Display Control" to visualize the cell bodies. The optimal display setting may vary between each cell body or dataset.

## Creating Mask

2




Fig 2.1

The critical classifiers when identifying cell bodies are the dark, round nucleus in the center and the "tail" that proceeds to form the axon (Fig 2.1).

Instead of outlining the cell bodies, the mask will remain inside the plasma membrane and follow the cell body from the first frame of the image to the last, including all abnormal projections.

Some cell bodies may be very dim. If the morphology of the cell body is unclear, conservatively follow the cell boundaries.




Fig 2.2




Fig 2.3

The "Paint" tool is used to annotate the cell bodies (Fig 2.2). Depending on the image's contrast, the "Masking" feature can be turned on (Fig 2.3). This will only select the voxels within the set threshold. The "halo" around the cell body is not part of the membrane and should not be annotated. To find more details on how to use the "Paint" tool, refer to the [Amira Annotation Protocol](#), Section 3.




Fig 2.4

Follow the cell body all the way, and cut off when the diameter reaches a consistent size. The uniform, elongated structure attached to the cell body is the axon which is not part of the mask.

Use the up and down arrow keys to navigate through the image frames. Make an annotation on one frame, then skip five slices to make another annotation. "Interpolation" automatically interpolates all frames between each selection. Review the result after selection, and make any needed changes before adding the selection to the material.

Change your orientation frequently to capture the three-dimensional cell body in all planes. Since the annotation is done on one plane at a time, when the orientation is changed, the annotation may appear jagged and irregular on other planes.

Use the "Smoothing" function to even out the membrane edges. The kernel size determines the intensity of smoothing. The bigger the size is, the more intense the smoothing effect becomes. Start with a smaller kernel size, 2, then examine the result in all three planes before increasing the kernel size. In Cellmap, the typical smoothing kernel size ranges from 2-5.




Fig 2.5

Once the cell body is annotated in all frames and smoothed, the annotation will appear as a round sphere with a peak on one end in the 3D viewer (Fig 2.5).

Then with the "Pick & Move" tool, select any region inside the cell body that is enclosed by the plasma membrane. Add the selection to the material to complete the mask. If there is a break in the annotation, the selection will overfill and select every unassigned voxel. Scroll through the frames to locate the problematic frame with the break in the membrane and enclose it.

Before completing the mask, there are a few common mistakes to double-check for. The halo around the cell body should not be included in the annotation (Fig 2.6).






Fig 2.6

There should be no holes in the mask (Fig 2.7)



Fig 2.7

The elongated section of the cell body that proceeds to form the axon (Fig 2.8) should not be part of the mask.




Fig 2.8

Each cell should have its own material. So, for example, if there are 11 cell bodies in the dataset, there should be 11 materials (Fig 2.9)




Fig 2.9

## Exporting Mask

- 3 The completed mask annotation needs to be exported as a 3D tiff file. In the "Project" workroom, right-click on the label source with the mask annotation. Then from the pop-up window, hit "Export Data As" then choose 3D tiff from the "Save as type" dropdown menu.