



Oct 07, 2020

Data Processing of Nano-DESI IMS data

yang1832¹, Julia Laskin¹¹Purdue University

1

Works for me

dx.doi.org/10.17504/protocols.io.bm2sk8ee

NanoDESI MSI Julia Laskin group Purdue



yang1832

ABSTRACT

Scope:

Process a mass spectrometry imaging data into a final imzML file.

DOI

dx.doi.org/10.17504/protocols.io.bm2sk8ee

PROTOCOL CITATION

yang1832, Julia Laskin 2020. Data Processing of Nano-DESI IMS data. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bm2sk8ee>



LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 05, 2020

LAST MODIFIED

Oct 07, 2020

PROTOCOL INTEGER ID

42802

ABSTRACT

Scope:

Process a mass spectrometry imaging data into a final imzML file.

- 1 Import the raw data folder into PeakByPeak Imaging software under "Process Raw folder"
- 2 Set "Number of Processes" to 10.
- 3 Choose "File Format" as imzML.
- 4 Run.

