



APR 06, 2023

OPEN  ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.5qpvor6xxv4o/v1

Protocol Citation: Nathan Heath Patterson, Jamie Allen, Angela Kruse, allison.b.esselman, ellie.l.pingry, Danielle Gutierrez, Jeff Spraggins 2023. Image Registration of MALDI IMS to Microscopy. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvor6xxv4o/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
 We use this protocol and it's working

Created: Apr 05, 2023

Last Modified: Apr 06, 2023

PROTOCOL integer ID:
 80064

Keywords: HuBMAP, BIOMIC, Vanderbilt, Multimodal Imaging, IMS, Image Registration

Image Registration of MALDI IMS to Microscopy

Forked from [Image Registration of MALDI IMS to Microscopy](#)

Nathan Heath Patterson¹, Jamie Allen¹, Angela Kruse¹, allison.b.esselman¹, ellie.l.pingry¹, Danielle Gutierrez¹, Jeff Spraggins¹

¹Vanderbilt University

VU Biomolecular Multimodal Imaging Center

Human BioMolecular Atlas Program (HuBMAP) Method Development Core



Angela Kruse
 Vanderbilt University

ABSTRACT

Scope:

To process to register MALDI IMS images to different types of microscopy images using IMS MicroLink and WSIREG.

- 1 Install Napari with IMS MicroLink plugin using <https://github.com/NHPatterson/napari-imslink>
- 2 Use MicroLink to align laser ablation marks from the post autofluorescence image to the IMS pixel map.
- 3 Install the WSIREG plugin with Napari <https://github.com/NHPatterson/wsireg>

- 4 Use WsiReg to align other microscopy modalities. These modalities are aligned to the post autofluorescence image.
- 5 When complete, all images will be sampled in the same coordinates as the IMS pixel map.