



4 ▼

Feb 14, 2022

Dephosphorylation of 5'-ends of DNA using rSAP in Restriction Enzyme Reaction (M0371) V.4

New England Biolabs¹¹New England Biolabs

1

dx.doi.org/10.17504/protocols.io.be59jg96**New England Biolabs (NEB)**Tech. support phone: +1(800)632-7799 email: info@neb.com**New England Biolabs**
New England Biolabs

This protocol explains methods for dephosphorylation of 5'-ends of DNA using rSAP in restriction enzyme reaction (M0371).

DOI

dx.doi.org/10.17504/protocols.io.be59jg96<https://www.neb.com/protocols/2013/06/10/protocol-for-dephosphorylation-of-5-ends-of-dna-in-restriction-enzyme-reaction-m0371>

New England Biolabs 2022. Dephosphorylation of 5'-ends of DNA using rSAP in Restriction Enzyme Reaction (M0371). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.be59jg96>
New England Biolabs



rSAP, Dephosphorylation

protocol ,

Apr 15, 2020

Feb 14, 2022

35745

MATERIALS

[Shrimp Alkaline Phosphatase \(rSAP\) - 500 units](#) **New England**

Biolabs Catalog #M0371S

Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

1 

Digest **1 µg – 5 µg plasmid DNA** in a 20 µl reaction as follows:

A	B
DNA	≥ 1 µl
Restriction Enzyme Buffer (10X)	2 µl
Restriction Endonuclease	1 µl
H ₂ O, purified	to 20 µl**

Scale larger reaction volumes proportionally.

2 

Incubate at **37 °C** for **01:00:00** or follow manufacturer's recommendations.

3 

Add **1 unit** of rSAP for every **1 pmol** of DNA ends (about **1 µg** of a 3 kb plasmid).

4 

Incubate at **37 °C** for 30-60 minutes.

5 Stop reaction by heat-inactivation of rSAP and restriction enzyme (follow manufacturer's recommendations).

If restriction enzyme cannot be heat-inactivated, DNA purification is required before ligation.

