

Jun 18, 2024

Transgenic mouse colony maintenance



In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.yxmvmeyx5g3p/v1

Ariadna Laguna¹, Iria Carballo-Carbajal¹, Miquel Vila²

¹vall d'hebron research institute; ²VHIR-CIBERNED-ASAP

ASAP Collaborative Rese...

Vilalab Public

1 more workspace



Miquel Vila

VHIR-CIBERNED-ASAP

OPEN ACCESS



DOI: dx.doi.org/10.17504/protocols.io.yxmvmeyx5g3p/v1

Protocol Citation: Ariadna Laguna, Iria Carballo-Carbajal, Miquel Vila 2024. Transgenic mouse colony maintenance. **protocols.io** https://dx.doi.org/10.17504/protocols.io.yxmvmeyx5g3p/v1

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's

working

Created: June 18, 2024

Last Modified: June 18, 2024

Protocol Integer ID: 102012

Funders Acknowledgement: Ministry of Science and

Innovation

Grant ID: PID2020-116339RB-

100



Abstract

Methods for the generation and maintenance of transgenic tyrosinase mice



- 1 Pronuclear microinjection of the human tyrosinase complementary DNA fused to the rat tyrosine hydroxylase promoter into C57BI6-SJL (RRID:MGI:5656718) mouse zygotes
- Selection of positive lines and backcrossing for 8-10 generations using 2 C57BL/6J mice (Charles River; RRID:MGI:3028467)
- 3 Maintainance in heterozygosis by crossing heterozygous males with C57BL76J females
- 4 Housed two to five per cage with ad libitum access to food and water during a 12 hours (h) light/dark cycle (light-on 8 a.m.)