

JAN 09, 2023

OPEN BACCESS

DOI:

<u>dx.doi.org/10.17504/protocol</u> <u>s.io.x54v9dr84g3e/v1</u>

Protocol Citation: Kaia Kukk 2023. Deglycosylation of N-glycosylated proteins using PNGase F . **protocols.io** https://dx.doi.org/10.17504/protocols.io.x54v9dr84g3e/v1

License: This is an open access protocol distributed under the terms of the <u>Creative Commons</u>
Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: Dec 09, 2022

Last Modified: Jan 09, 2023

PROTOCOL integer ID:

73770

1 12,5 µl of yeast lysate containing sufficient amount of target protein for detecting by Western analysis was pipetted into a tube. Two samples were prepared in parallel, one for negative control without *N*-glycosidase and the other with PNGase F.

Obeglycosylation of N-glycosylated proteins using PNGase F

Kaia Kukk¹

¹Latvian Biomedical Research and Study Centre



Kaia Kukk

Latvian Biomedical Research and Study Centre

ABSTRACT

The described protocol was used to confirm that NADPH-cytochrome P450 reductase from *Helianthus annuus* was *N*-glycosylated when recombinantly expressed in *Pichia pastoris* (*Komagataella phaffii*).

MATERIALS

PNGase F from Elizabethkingia meningoseptica **Sigma Aldrich Catalog**#P7367-50UN

2 0,5 µl of 10% sodium dodecyl sulfate (SDS) and 1 µl of 1M dithiothreitol (DTT) were added. 3 The mixtures were incubated at 95 °C for 5 min. 4 The mixtures were cooled to room temperature after which 2 µl of 0,5 M Tris-HCl, pH 8 and 2 µl of 10% Triton X-100 were added. 5 2 μl of PNGase F (500 U/ml) was added to one of the samples. Instead of N-glycosidase, water was added to the negative control sample. 6 The samples were incubated overnight at 37 °C. 7 SDS-PAGE sample buffer was added and the samples were heated for 5 min at 99 °C. SDS-PAGE and Western analysis were carried out to see whether the molecular weight of the target protein had decreased.