



Apr 20, 2022

Counting kidneys

Gabriel J Barrero¹, Abraham Palmer¹, Oksana Polesskaya¹

1	п	п	\cap	C	г
			◡	J	L



dx.doi.org/10.17504/protocols.io.6qpvr6k4pvmk/v1

CGORD



Dissecting and counting kidneys in rats.

DOI

dx.doi.org/10.17504/protocols.io.6qpvr6k4pvmk/v1

Gabriel J Barrero, Abraham Palmer, Oksana Polesskaya 2022. Counting kidneys. **protocols.io**

https://dx.doi.org/10.17504/protocols.io.6qpvr6k4pvmk/v1

____ protocol,


Apr 20, 2022

Apr 11, 2022


60610

- 1 Rats are euthanized as per IACUC guidelines Spray fur with ethanol. Pick up skin above the sternum, and make a cut through skin and muscle layers.
 - Don't cut the chest cavity (don't go above the diaphragm).
- 2 Cut skin and muscles from the sternum down to the lower belly. Mostly you will see liver and colon.






3 Lift liver lobes and push aside fat tissue and intestines. You will see dark-colored organs – kidney and spleen.



Gross dissection



Representative renal anomalies observed in rats. A. Grossly normal renal development. B. Unilateral renal agenesis (URA) PloS one. 10. e0118147. 10.1371/journal.pone.0118147

4 Count the number of observed kidneys and make a record in the file.