




FEB 13, 2023

🌐 Insert Alignment ELYRA7

Rafael Camacho¹

¹Centre for Cellular Imaging, Core Facilities, Gothenburg University

CCI_Gothenburg



Rafael Camacho

Centre for Cellular Imaging, Core Facilities, Gothenburg Uni...

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

Alignment of the stage insert to ensure that the sample is perpendicular to the optical axis. This protocol is based on material provided by ZEISS' application specialists.

What is needed before to start:

- This protocol was designed for the ELYRA7
- You will need 2 microscopy slides, one of standard size and one of about 2x2cm

OPEN ACCESS

Protocol Citation: Rafael Camacho 2023. Insert Alignment ELYRA7.

protocols.io

<https://protocols.io/view/insert-alignment-elyra7-cmu4u6yw>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jan 17, 2023

Last Modified: Feb 13, 2023

PROTOCOL integer ID:
75388

Turn on the microscope

1 Turn on the ELYRA 7

1.1 Turn on **Main Switch**



1.2 Turn on Components



1.3 Turn on the Microscope Computer

Note

In the ELYRA 7 if the microscope computer is not on, then the microscope's touch screen will not light up.

1.4 Open **ZEN black**, and choose *start system*


Find sample plane via the eye pieces

2

Note


This protocol assumes that the microscope is on, and ready for operation

ZEN black -> **Locate Tab,**



View of **Tabs** on ZEN black, top left side

and set an imaging protocol for transmitted light so you can see the sample via the eye pieces.



Microscope Control example for transmission

Software


ZEN black

NAME

ZEISS

DEVELOPER

3 Change the objective to the lowest magnification, typically 10X



*You can change objectives in the **Microscope Control** panel*

4 Prepare a clean microscopy slide having a few marker points, *marker-slide*. You can keep track of the surface on which you made the marks by writing TOP and DOWN


Note

The slide can be reused




Example of marker-slide used to move to find the sample plane by focusing with the low magnification objective

- 5 Place the *marker-slide* on the sample holder so you can focus on the appropriate sample plane



Marker-slide on the slide holder


- 6 Focus on the sample via the eye pieces




Marker spot in focus as viewed via the eye pieces

Place glasses for reflection mode

- 7 Remove the sample holder carefully and place it a side



- 8 Remove the objective lens and carefully place it a side



- 9 Prepare a clean microscopy slide of size approximately 2x2cm



Small slide can be made by cutting a larger slide into the appropriate size

10


Place the 2x2cm clean glass plate on the objective turret.



Safety information

To avoid adding fingerprints to the slide you can wear protective gloves.

Make sure not to touch the copper pin with with glass. This pin is used to read the objectives collar, and it is delicate.




Place back the sample holder (which still has on it the marker sample)

11




12 Close the imaging chamber



Acquisition Settings

13 Go to ZEN black -> **Acquisition Tab**



14 Set the acquisition settings

Step 14 includes a Step case.

You have settings

You need settings

step case


You have settings

15 Load the acquisition settings by using the **Experiment Manager**

In our case these settings are stored as "insert alignment"

Image Acquisition

16 Go to **Live mode** and make sure image is not saturated




- 17 Find where the fibre of the HXP is connected to the microscope body. This will be on the left hand side of the microscope's body, after the on/off button




18

While in **Live mode**, carefully take out the fibre from the HXP that is connected to microscope body until you see the sharp edges of the fibre in the field of view



19

Use the levelling screws of the holder to match the 2 sets of reflections to one another



- 20 Take an image **snap**, and save the image for your own record
- 21 Place back the HXP fibre into position so you do not see the edges of the fibre while in **live** view

Place objective back into the turret

- 22 Remove the sample holder carefully so you don't touch/move the levelling screws



Note

From now on be very careful of not moving the levelling screws

23 Remove the glass slide placed at the base of the objective port

24 Place back the objective

25 Place back the sample holder

Safety information

Be careful not to touch the levelling screws and not to hit the objective lens with the holder.