

May 24, 2024 Version 1

Sucrose lysis buffer V.1

Forked from DNA EXTRACTION Protocol Template

This protocol is a draft, published without a DOI.



Colleen Kellogg¹

¹Hakai Institute

Better Biomolecular Ocea...

Jaanasheen Bala
UBC

OPEN ACCESS



Protocol Citation: Colleen Kellogg 2024. Sucrose lysis buffer. protocols.io <https://protocols.io/view/sucrose-lysis-buffer-dc5u2y6w>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol

Created: November 07, 2023

Last Modified: May 24, 2024

Protocol Integer ID: 99220

Disclaimer

Draft!

Abstract

This protocol describes the preparation of sucrose lysis buffer to preserve DNA on sterivex filters. The protocol is part of the Hakai Institute's pipeline to analyze microbial and environmental DNA from seawater samples and is implemented as a standard procedure for ongoing sampling programs.

Materials

DESCRIPTION e.g. filter	PRODUCT NAME AND MODEL Provide the official name of the product	MANUFACTURER Provide the name of the manufacturer of the product.	QUANTITY Provide qua
Durable equipment			
Content Cell	Content Cell	Content Cell	Content Cell
Content Cell	Content Cell	Content Cell	Content Cell
Consumable equipment			
Content Cell	Content Cell	Content Cell	Content Cell
Content Cell	Content Cell	Content Cell	Content Cell
Chemicals			
Content Cell	Content Cell	Content Cell	Content Cell
Content Cell	Content Cell	Content Cell	Content Cell

Before start

Read background information, MIOP and BePOP-OBON information under the "Guidelines" tab.



METHODS

1

Note

Wear gloves.

Equipment:

1. 1 Molarity (M) Tris pH 8 (made previously or purchased as solution)
2. 0.5 Molarity (M) EDTA pH 8 (made previously or purchased as solution)
3. Sucrose
4. MilliQ
5. 500 mL bottle top filtration unit

Final concentrations of chemicals in SLB:

EDTA: 40 Molarity (M)

Tris: 50 Molarity (M)

Sucrose: 0.75 Molarity (M)

- 2 Calculate how much Sucrose powder you need (using the molecular weight on the bottle, MW or FW) for a 0.75 Molarity (M) solution of 500 mL .

Calculation of Sucrose:

$$\frac{Y \text{ g}}{\text{mol}} \times \frac{0.75 \text{ mol}}{\text{L}} \times 0.5 \text{ L} = Z \text{ grams of Sucrose to add in step 2}$$

Where Y is the molecular weight of the sucrose (MW or FW) from the bottle.

- 3 Weigh out and add appropriate amount of Sucrose calculated in step 1 and place in a bottle or beaker.

- 4 Add 40 mL of 0.5 Molarity (M) EDTA to the beaker.

Calculation of Tris or EDTA:

Use the equation $C1V1 = C2V2$

Eg for EDTA:

$$(0.5M)(X \text{ mL}) = (0.04M)(500 \text{ mL})$$

Solve for X

$$((0.04M)(500 \text{ mL})/(0.5M)) = 40 \text{ mL of } 0.5M \text{ EDTA}$$

- 5 Add 25 mL of 1 Molarity (M) Tris to the beaker.

Calculation of Tris or EDTA:

Use the equation $C1V1 = C2V2$

Eg for EDTA:

$$(0.5M)(X \text{ mL}) = (0.04M)(500 \text{ mL})$$

Solve for X

$$((0.04M)(500 \text{ mL})/(0.5M)) = 40 \text{ mL of } 0.5M \text{ EDTA}$$

- 6 Add milliQ water to about the 400 mL line.

- 7 Add a stir bar and dissolve all powder.



- 8 Top up water to 500 mL. (no need to pH this one!)
- 9 Filter-sterilize and label bottle.