



VERSION 1

MAR 26, 2024

🌐 GeoMx Digital Spatial Profiler (DSP) Protocol - University of Minnesota TMCs V.1

Laura J Niedernhofer¹, David A Bernlohr¹

¹University of Minnesota, Minneapolis, MN

Cellular Senescence Network (SenNet) Method Development Community

UMN SenNet



Allie Pybas
UMN

OPEN  ACCESS



ABSTRACT

The nanoString GeoMx[®] Digital Spatial Profiler (DSP) is a platform that allows high-plex profiling at the protein and RNA level, providing spatial and temporal assessment of tumors in frozen or formalin-fixed paraffin-embedded limited tissue sample.

The following protocols/manuals were used at the University of Minnesota TMCs in collaboration with the University of Minnesota Genomics Center.

DOI:

dx.doi.org/10.17504/protocols.io.8epv5r7k6g1b/v1

Protocol Citation: Laura J Niedernhofer, David A Bernlohr 2024. GeoMx Digital Spatial Profiler (DSP) Protocol - University of Minnesota TMCs . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8epv5r7k6g1b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 22, 2024


Last Modified: Mar 26, 2024

PROTOCOL integer ID: 97153

Funders Acknowledgement:
NIH
Grant ID: 5U54AG076041-03
NIH
Grant ID: 5U54AG079754-02

Slide Preparation

1 GeoMx DSP Manual Slide Preparation User Manual

 MAN-10150-01.pdf 2.6MB

ROI Acquisition

2 GeoMx DSP Instrument User Manual

 MAN-10152-01.pdf 6.1MB

Library Preparation

3 GeoMx DSP NGS Readout User Manual

 MAN-10153-01.pdf 5.9MB

Note

Sequence with the read format 27,8,8,27

FASTQ Generation

- 4 BCL data from Illumina sequencer is demultiplexed and converted into FASTQ format using bcl2fastq version 2.20.0