



MAR 15, 2024

Bulk RNA sequencing

anita.adami¹

¹Laboratory of Molecular Neurogenetics, Department of Experimental Medical Science, Wallenberg Neuroscience Center and Lund Stem Cell Center, BMC A11, Lund University, 221 84 Lund, Sweden.

ASAP Collaborative Research Network

Jakobsson



Raquel Garza
Lund University

ABSTRACT

This protocol describes how to perform bulk RNA sequencing

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.36wgqjqbkvk5/v1

Protocol Citation: anita.adami
2024. Bulk RNA sequencing.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.36wgqjqbkvk5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 03, 2023

Last Modified: Mar 15, 2024

PROTOCOL integer ID: 78066

Keywords: ASAPCRN

Funders Acknowledgement:

Aligning Science Across
Parkinson's through the Michael
J. Fox Foundation for Parkinson's
Research

Grant ID: ASAP-000520

Swedish Research Council

Grant ID: 2018-02694

Swedish Brain Foundation

Grant ID: FO2019-0098

Cancerfonden

Grant ID: 190326

Barncancerfonden

Grant ID: PR2017-0053

NIHR Cambridge Biomedical
Research Centre

Grant ID: NIHR203312

Swedish Society for Medical
Research

Grant ID: S19-0100

National Institutes of Health

Grant ID: HG002385

Swedish Research Council

Grant ID: 2021-03494

Swedish Research Council

Grant ID: 2020-01660

Bulk RNA sequencing

- 1 Total RNA was isolated from nuclei ([Nuclear Extraction from Tissue for FACS](#)), cell culture samples, or tissue using the RNeasy Mini Kit (Qiagen).
- 2 Libraries were generated using Illumina TruSeq Stranded mRNA library prep kit (poly-A selection) and sequenced on a NextSeq500 (PE 2×150bp).