




MAR 27, 2023

Mitochondrial complex activity assays

michela.deleidi¹, María José Pérez J.¹, Hariam Raji¹,
Pascale Baden¹, Federico Bertoli¹

¹German Center for Neurodegenerative Diseases (DZNE), Tübingen,
72076 Germany



Federico Bertoli

ABSTRACT

Mitochondria complex activity assays measure the activity levels of the different complexes of the mitochondrial electron transport chain (ETC).

ATTACHMENTS

[676- 1425.docx](#)

OPEN ACCESS

DOI:

dx.doi.org/10.17504/protocols.io.eq2ly753elx9/v1

Protocol Citation: michela.deleidi, María José Pérez J., Hariam Raji, Pascale Baden, Federico Bertoli 2023. Mitochondrial complex activity assays. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.eq2ly753elx9/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 17, 2023

Last Modified: Mar 27, 2023

PROTOCOL integer ID:
78971

Keywords: Mitochondrial complex activity assays, mitochondrial electron transport chain

MATERIALS

Materials

- pyruvate
- malate
- ADP
- Succinate
- rotenone
- antimycin A
- TMPD (N,N,N',N'-tetramethyl-p-phenylenediamine dihydrochloride, Santa Cruz Biotechnology)
- ascorbic acid
- azide
- Qproteome Mitochondrial isolation kit



MitoCheck® Complex I Activity Assay Kit **Cayman Chemical**
Company Catalog #700930



Qproteome Mitochondria Isolation Kit **Qiagen Catalog #37612**


MAS buffer

A	B
Sucrose	70 mM
Mannitol	220 mM
KH ₂ PO ₄	5 mM
MgCl ₂	5 mM
EGTA	1 mM
HEPES pH 7.4	2 mM

Mitochondrial complex activity assays










35m

- 1 Isolate mitochondria from HEK cells, iPSC-derived neurons, or midbrain organoids using the Qproteome Mitochondrial isolation kit (QIAGEN, Cat. No. / ID: 37612) according to manufacturer's instructions.
- 2 Measure Complex I (NADH oxidase/coenzyme Q reductase) using the MitoCheck Complex I Activity Assay kit (Cayman Chemical, cat# 700930).

- 3 Determine the rate of NADH oxidation, which is proportional to CI activity, by a decrease in absorbance at 340 nm over  00:15:00 in the presence of ubiquinone and potassium cyanide to inhibit complex IV and prevent oxidation of ubiquinone. 15m







Note


To assess CI, CII, and CIV function, we used a respirometry approach based on XFp Extracellular Flux Analysis and then proceed with steps 4-9.

- 4 To this end, resuspend  3 mg of purified fresh mitochondria in  200 μ L of MAS buffer ( 70 millimolar (mM) sucrose,  220 millimolar (mM) mannitol,  5 millimolar (mM) KH_2PO_4 ,  5 millimolar (mM) MgCl_2 ,  1 millimolar (mM) EGTA,  2 millimolar (mM) HEPES  7.4) and seed in XFpSeahorse microplates. 5m

- 5 Centrifuge the plate at  2000 x g, 4°C, 00:05:00 . 5m



- 6 Measure the OCR before and after the serial addition of pyruvate + malate ( 5 millimolar (mM) each) + ADP 3,5 mM or 1 mM Succinate +  4 micromolar (μ M) rotenone,  4 micromolar (μ M) rotenone +  8 micromolar (μ M) antimycin A, 0,5 mM TMPD (N,N,N',N'-tetramethyl-p-phenylenediamine dihydrochloride, Santa Cruz Biotechnology) +  1 millimolar (mM) ascorbic acid, and  50 millimolar (mM) azide. 15m

- 7 Following each injection, record three measurements for a total period of  00:15:00 . 15m

- 8 Calculate Complex I-, II-, and IV-dependent respiration by subtracting OCR values from the substrates (Pyruvate + malate + ADP for CI, Succinate + rotenone for CII and TMPD + ascorbic acid for CIV) subtracted from the ones from the inhibitors (rotenone for CI, antimycin A + rotenone for CII and azide for CIV).

- 9 Normalize the experimental values to the protein content per well via a BCA assay.

